Precision Xtra®

Blood Glucose Test Strips

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Read This First


IMPORTANT: Read these instructions for use and the User's Manual supplied with your meter before you monitor your blood glucose. Failure to follow instructions will cause incorrect results.

Blood Glucose Test Strips with TrueMeasure Technology

English

What are my test strips for?

Precision Xtra blood glucose test strips are for use with the Precision Xtra Meter. The Precision Xtra blood glucose test strips are intended for *in-vitro* diagnostic (outside the body) use in the quantitative measurement of glucose in fresh capillary whole blood for self-testing by lay users (e.g., from the (a) finger, (b) forearm, (c) upper arm or (d) base of thumb), or by healthcare professionals. It is intended for use by healthcare professionals and people with diabetes mellitus as an aid in monitoring the effectiveness of a diabetes control program. It is not intended for the diagnosis of or screening for diabetes mellitus, and it is not intended for use on neonates , arterial or venous blood samples.

What's in my test strip box?

- Test strips individually wrapped in foil packets
- Instructions for use

What else do I need that is not in my test strip box?

- Precision Xtra Meter
- User's Manual
- MediSense Glucose and Ketone Control Solutions
- Lancing device and disposable lancets

How do I check my meter?

Control Solutions are used to check the performance of meters, test strips, and your testing technique. Perform a control solution test when using the meter for the first time, when opening a new box of strips, or when you want to confirm that your meter and test strips are working properly.

When a control solution test is done, you should get a result within the expected range printed on these instructions for use.

If a control solution test result falls outside this range, repeat the test.

Results that fall outside the range may be caused by:

- error in performing the test
- expired test strip
- expired or contaminated control solution
- test strip deterioration
- meter malfunction

WARNING: If you continue to get MediSense control solution test results that fall outside the range printed on these test strip instructions for use, the system may not be functioning properly. DO NOT use the system to test your blood until you get a control solution test result within the range printed on these test strip instructions for use. Contact Customer Service at 1-800-527-3339. For information on how to obtain control solutions, please contact Customer Service at 1-800-527-3339.

How do I obtain a blood drop?

- Before you obtain a blood drop, make sure the sample site is clean, dry, and warm. To warm the sample site, wash it in warm water, rub the skin vigorously for a few seconds, or apply a warm pad to it.
- Hang your arm down before lancing your finger or base of the thumb, to help blood flow.
- To obtain a blood drop from the arm, use a fleshy area away from bone. Avoid areas where there is a lot of hair.
- Apply the blood drop to the test strip immediately.

How do I monitor my blood glucose?

1. Remove the test strip from its foil packet. Open the test strip packet by tearing at the notch.
2. Insert the contact bars at the end of the test strip into the test port of the meter.

Gently push the test strip in until it stops. The meter turns on automatically.

Note: Check that LOT 45001 appears on the meter display window. If LOT 45001 is not displayed, contact Customer Service.

3. Obtain a blood drop. Follow the instructions for use packaged with the lancing device.
4. Touch the blood drop to the white target area at the end of the test strip.

The blood is drawn into the test strip.

IMPORTANT: After performing a blood glucose test, wash your hands and test site (for example, your forearm) with soap and water. Thoroughly dry your hands and test site.

What if the countdown does not start? If the countdown does not start, you may not have applied enough blood to the test strip. Apply a second drop of blood to the test strip within 5 seconds of the first drop. If the countdown still does not start or if more than 5 seconds have passed, discard the test strip, turn off your meter, and repeat steps 1-4.

You can use the opened foil packet to remove and discard your used blood glucose test strip. Discard the test strip properly.

What does my result mean?

The expected glucose range for a non-diabetic, non-pregnant fasting adult is less than 100 mg/dL.¹ Two hours after meals, levels should be less than 140 mg/dL.¹

Consult your healthcare professional to determine the range that is appropriate for you.

IMPORTANT - What if my result is unusually high or low?

Low or high blood glucose results can indicate a potentially serious medical condition. If your blood glucose result is unusually low or high, or if you do not feel the way your result indicates, repeat the test with a new test strip. If your blood glucose result is less than 50 mg/dL or higher than 300 mg/dL, you should contact your healthcare professional and follow their treatment advice. You may also use a control solution to check the performance of your system.

Follow your healthcare professional's advice before you make any changes to your diabetes medication program.

IMPORTANT - How do I take care of my test strips?

- Use the test strip immediately after opening its foil packet.
- Your test strips should be stored at a temperature between 39° – 86° F (4° – 30° C). Storage outside this range may cause incorrect results. Keep away from direct sunlight and heat.
- Use each test strip once and then discard it. It must not be used on more than one person including other family members due to the risk of spreading infection. All parts of the strip are considered biohazardous and can potentially transmit infectious diseases.
- Do not use out-of-date test strips. Check the expiration date printed on the test strip box and on every test strip foil packet. If only the year and month are printed on the test strip, then the expiration date is the last day of the month. For example, "EXP 2014/03" means the test strip expires on March 31, 2014.
- Do not use a test strip that is wet, bent, scratched or damaged.
- Do not use the test strip if its foil packet has a puncture or tear in it.
- Observe caution when using around children. Small parts may constitute a choking hazard.

What else do I need to know?

- The Precision Xtra meter can read blood glucose levels between 20 and 500 mg/dL.
- Use Precision Xtra Blood Glucose Test Strips at temperatures between 59° – 104° F (15° – 40° C) and 10% and 90% relative humidity (the amount of moisture in the air) for best results.
- Clinical testing demonstrates that altitudes up to 7200 feet (2195 meters) above sea level do not affect results.

IMPORTANT - If you choose to use a sample from the forearm, upper arm, or base of the thumb:

- Consult your healthcare professional before you use any one of these sites to monitor your blood glucose.
- Contact Customer Service at 1-800-527-3339 for further information. Outside the United States, contact your local Abbott Diabetes Care office or distributor for further information.
- Sampling from any one of these alternate sites may cause minor bruising and may leave marks that go away in a short time.
- There may be times when alternate site results are different from fingertip results. This happens when glucose levels change rapidly (for example, after you eat a meal, take insulin, or during or after exercise).
- Use alternate sites to monitor before, or more than two hours after, you eat a meal, take insulin, or exercise.
- Do **not** use blood samples from alternate sites when:
 1. You think your blood sugar is low,
 2. You have been diagnosed with hypoglycemic unawareness,
 3. The results from alternate sites do not match the way you feel, or
 4. It is within two hours of eating a meal, taking insulin, or exercising.

Are there important messages that I need to know about?

The following messages may mean you have obtained a blood glucose result that requires immediate attention or there may be a problem with the test strip:

- **LO** means your blood glucose may be lower than 20 mg/dL.
- **HI** means your blood glucose may be higher than 500 mg/dL.
- **Test Error 2** or **E-3** means there may be a test error or means your blood glucose may be too low to be read by the system.
- **Test Error 4** or **E-4** means your blood glucose may be too high to be read by the system.

If any of these messages show, repeat the test with a new test strip. If the same message shows again or the result does not reflect how you feel, contact your healthcare professional **immediately**. You may also use a Control Solution to check the performance of your system. Follow your healthcare professional's advice before you make any changes to your diabetes medication program.

IMPORTANT: The correct measurement units must show on your meter with every result. See below for the measurement units used in your country.

US	mg/dL
----	-------

If you have any questions and/or need assistance, please contact Customer Service.

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1-800-527-3339

Referencias

1. American Diabetes Association Position Statement: Diagnosis and classification of diabetes mellitus. Diabetes Care 2012;35(Suppl.1):S64-S71.

Referencias

1. American Diabetes Association Position Statement: Diagnosis and classification of diabetes mellitus. Diabetes Care 2012;35(Suppl.1):S64-S71.

ENGLISH

Important Information for Healthcare Professionals

Note: Capillary blood may be collected into heparin-containing or potassium EDTA-containing capillary tubes and used within 30 minutes. Do not use tubes containing fluoride or oxalate.

Limitations of Procedure

- This test strip is not designed for use with arterial, venous, neonatal, serum or plasma samples.
- Hematocrit range is 30%-60%.
- High levels of acetaminophen, up to 25 mg/dL, will not affect results.
- Test results may be erroneously low if the patient is severely dehydrated, or severely hypotensive, in shock or in a hyperglycemic-hypermolar state (with or without ketosis). Similar observations have been reported in the literature for other blood glucose monitoring systems.
- Extremely high levels of the following substances at the following concentrations do not affect results: uric acid, 24 mg/dL; ascorbic acid, 4 mg/dL; unconjugated bilirubin, 40 mg/dL; cholesterol, 500 mg/dL; and triglycerides, 1500 mg/dL.
- Do not use during xylose absorption testing.

Checking the System

Use appropriate levels of MediSense Control Solutions to do control solutions tests with each new box of test strips.

Test Principle

When the blood sample is applied to the test strip, the glucose in the blood reacts with the chemicals on the test strip, producing a small electrical current. This current is measured and a result is then displayed by the meter. The size of the current depends on the amount of glucose in the blood sample.

Composition

Glucose Dehydrogenase (<i>GDH-NAD Pseudomonas sp.</i>)	≥ 0.03 U
NAD+ (as sodium salt)	≥ 1.0 µg
Phenanthroline quinone	≥ 0.02 µg
Non-reactive ingredients	≥ 16.3 µg

Performance Characteristics

The performance of Precision Xtra Blood Glucose Test Strips has been evaluated in both laboratory and clinical studies.

Assay Range: 20–500 mg/dL **Test Time:** 5 seconds

Calibration Reference

The Precision Xtra Blood Glucose Test Strip is calibrated against the YSI Glucose Analyser (YSI Inc.). The YSI whole blood glucose results are multiplied by 1.12 to provide plasma-equivalent glucose values for the calibration of Precision Xtra Blood Glucose Test Strips.

Precision

Precision testing shows that results typically vary by no more than 2.7% to 4.0%. Please see Table 1.

Accuracy

Capillary blood glucose results obtained by lay users at 3 clinical centers were compared with those obtained using the YSI Analyzer, a laboratory instrument. Please see Table 2.

These studies show that the Precision Xtra System compares well with the laboratory reference method.

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© 2011 Abbott
Abbott Diabetes Care Inc.