
### Read This First

**IMPORTANT: Read these instructions for use and the User's Manual supplied with your meter before you monitor your blood ketone. Failure to follow instructions will cause incorrect results.**

### English

#### What are my test strips for?

Precision Xtra Blood β-Ketone Test Strips are for use with the Precision Xtra Meter. The Precision Xtra Blood β-Ketone Test Strips are intended for *in-vitro* diagnostic (outside the body) use and are for self-testing or healthcare professional use. The test strips are intended to quantitatively measure blood Ketone (Beta-Hydroxybutyrate) in fresh capillary whole blood from the fingertip. Healthcare professionals may also use venous whole blood samples, provided the samples are used within 30 minutes of collection. The test strips are not intended for use in diagnosis or screening of diabetes mellitus, but are to be used as an aid in monitoring the effectiveness of diabetes control programs.

#### What's in my test strip box?

- Test strips individually wrapped in foil packets
- Instructions for use

#### What else do I need that is not in my test strip box?

- Precision Xtra Meter
- MedSense Glucose and Ketone Control Solutions
- User's Manual
- Lancing device and disposable lancets

#### How do I check my meter?

Control Solutions are used to check the performance of meters, test strips, and your testing technique. Perform a control solution test when using the meter for the first time, when opening a new box of test strips, or when you want to confirm that your meter and test strips are working properly. When a control solution test is done, you should get a result within the expected range printed on these instructions for use. If a control solution test result falls outside this range, repeat the test.

Results that fall outside the range may be caused by:

- error in performing the test
- expired test strip
- expired or contaminated control solution
- test strip deterioration
- meter malfunction

**WARNING:** If you continue to get MedSense control solution test results that fall outside the range printed on these instructions for use, the system may not be functioning properly. DO NOT use the system to test your blood until you get a control solution test result within the range printed on these instructions for use. Contact Customer Service at 1-800-527-3339.


For information on how to obtain Control Solutions, please contact Customer Service at 1-800-527-3339.

#### How do I obtain a blood drop?

- Before you obtain a blood drop, make sure your fingertips are clean, dry, and warm. To warm your fingertips, wash your hands in warm water.
- Hang your arm down before lancing your fingertip, to help blood flow.
- Apply the blood drop to the test strip immediately.

#### How do I monitor my blood ketone?

1. Remove the test strip from its foil packet. Open the test strip packet by tearing at the notch.
2. Insert the contact bars at the end of the test strip into the test port of the meter. Gently push the test strip in until it stops. The meter turns on automatically. **Note:** Check that LOT 75001 or CODE 75001 appears on the meter display window. If LOT 75001 or CODE 75001 is not displayed, contact Customer Service.
3. Obtain a blood drop. Follow the instructions for use packaged with the lancing device.
4. Touch the blood drop to the white target area at the end of the test strip. The blood is drawn into the test strip.

**What if the countdown does not start?** The meter displays the blood ketone result in 10 seconds. If the countdown does not start, you might not have applied enough blood to the test strip. Apply a second drop of blood to the test strip within 30 seconds of the first drop. If the countdown still does not start, or if more than 30 seconds have passed, discard the test strip, turn off your meter, and repeat steps 1-4.

You can use the opened foil packet to remove and discard your used blood β-Ketone test strip. Discard the test strip properly.

#### What does my result mean?

The blood ketone test measures Beta-Hydroxybutyrate (β-OHB), the most important of the three ketone bodies in the blood<sup>1</sup>. Normally, levels of β-OHB are expected to be less than 0.6 mmol/L<sup>2</sup>. β-OHB levels may increase if a person fasts, exercises vigorously or has diabetes and becomes ill<sup>1,3</sup>.

If your blood ketone result is 0.0 mmol/L and your blood glucose result is higher than 300 mg/dL, repeat both the ketone and the glucose tests with new test strips. If you get similar results again, contact your healthcare professional. Follow your healthcare professional's advice before you make any changes to your diabetes medication program.

If your blood ketone result is between 0.6 and 1.5 mmol/L and your blood glucose result is higher than 300 mg/dL, this may indicate development of a problem that could require medical assistance. Follow your healthcare professional's instructions.

If your blood ketone result is higher than 1.5 mmol/L and your blood glucose result is higher than 300 mg/dL, contact your healthcare professional **immediately** for advice and assistance. You may be at risk of developing diabetic ketoacidosis (DKA)<sup>4,6</sup>.

#### IMPORTANT - How do I take care of my test strips?

- Use the test strip immediately after opening its foil packet.
- Your test strips should be stored at a temperature between 39° – 86°F (4° – 30°C). Storage outside this range may cause incorrect results. Keep away from direct sunlight and heat.
- Use each test strip once and then discard it. It must not be used on more than one person including other family members due to the risk of spreading infection. All parts of the strip are considered biohazardous and can potentially transmit infectious diseases.
- Do not use out-of-date test strips. Check the expiration date printed on the test strip box and on every test strip foil packet. If only the year and month are printed on the test strip, then the expiration date is the last day of the month. For example, "EXP 2014/03" means the test strip expires on March 31, 2014.
- Do not use a test strip that is wet, bent, scratched or damaged.
- Do not use the test strip if its foil packet has a puncture or tear in it.
- Observe caution when using around children. Small parts may constitute a choking hazard.

#### What else do I need to know?

- The Precision Xtra meter can read blood β-Ketone levels between 0.0 and 8.0 mmol/L.
- Use Precision Xtra Blood β-Ketone Test Strips at temperatures between 64° – 86°F (18° – 30°C) and 10% and 90% relative humidity (the amount of moisture in the air) for best results.
- Clinical testing demonstrates that altitudes up to 7200 feet (2195 meters) above sea level do not affect results.

#### Are there important messages that I need to know about?

The following messages may indicate you have obtained a blood ketone result that requires immediate action or there may be a problem with the test strip:

- **HI** means your meter has determined that your blood ketone result is higher than 8.0 mmol/L.
- **Test Error 2** or **Test Error 4** or **E-3** or **E-4** means there may be a test error.

If any of these messages show, repeat the test with a new test strip. If the same message shows again, contact your healthcare professional **immediately**. You may also use a MedSense Control Solution to check the performance of your system. Follow your healthcare professional's advice before you make any changes to your diabetes medication program.

**If you have any questions and/or need assistance, please contact Customer Service.**

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**1-800-527-3339**

## ENGLISH

### Important Information for Healthcare Professionals

**Note: Venous whole blood samples may be collected into sodium or lithium heparin tubes or EDTA tubes and used within 30 minutes. Do not use tubes containing fluoride or oxalate.**

#### Limitations of Procedure

- This test strip has not been evaluated for alternative site testing.
- This test strip is not designed for use with arterial, neonatal, serum or plasma samples.
- Hematocrit range is 30%-60%.
- Test Results may be erroneously low if the patient is severely dehydrated, or severely hypotensive, in shock or in a hyperglycemic-hypermolar state.
- The following substances have no significant effect on blood ketone results:
  - Captopril up to 500 µg/dL
  - Dopamine up to 90 µg/dL
  - Acetaminophen up to 25 mg/dL
  - Ascorbic acid up to 4 mg/dL
  - Cholesterol up to 500 mg/dL
  - L-DOPA up to 600 µg/dL
  - Gentisic acid up to 1.8 mg/dL
  - Uric acid up to 24 mg/dL
  - Unconjugated bilirubin up to 20 mg/dL
  - Triglycerides up to 1875 mg/dL

#### Checking the System

Use appropriate levels of MediSense Control Solutions to do control solutions tests with each new box of test strips.

#### Test Principle

When the blood sample is applied to the test strip, the β-OHB in the blood reacts with the chemicals on the test strip, producing a small electrical current. This current is measured and a result is then displayed by the meter. The size of the current depends on the amount of β-OHB in the blood sample.

#### Composición

β-Hydroxybutyrate Dehydrogenase ( <i>Pseudomonas sp</i> )	≥	0.03 U
NAD (Free acid form)	≥	1.67 µg
Phenanthroline quinone	≥	0.29 µg
Non-reactive ingredients	≥	19.51 µg

**Assay Range:** 0.0 - 8.0 mmol/L

**Test Time:** 10 seconds

**Sample Volume:** 1.5 µL

#### Calibration Reference

The Precision Xtra Blood β-Ketone Test Strip is calibrated to reflect plasma β-hydroxybutyrate using the Randox assay kit (RB1007).

#### Precision

Precision testing shows that results typically vary by no more than 3.1% to 3.8%. Results were obtained in a laboratory study using venous whole blood samples (n = 20 per level). Please see Table 1.

#### Accuracy

Accuracy testing shows that results are comparable between trained operators and lay users. Accuracy was assessed at three clinics and four hospitals by comparing whole blood β-OHB results with plasma results obtained using a reference laboratory instrument. Please see Table 2. These studies show that the Precision Xtra System compares well with the laboratory reference method.

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#### References / Referencias

- Schade DS, Eaton RP. Metabolic and clinical significance of ketosis. Special Topics in Endocrinology and Metabolism 1982;4:1-27.
- Wiggam MI, O’Kane MJ, Harper R, Atkinson AB, Hadden DR, Trimble ER, Bell PM. Treatment of diabetic ketoacidosis using normalization of blood 3-hydroxybutyrate concentration as the endpoint of emergency management. Diabetes Care 1997;20:1347-52.
- Harano Y, Kosugi K, Hyosu T, Suzuki M, Hidaka H, Kashiwagi A, Uno S, Shigeta Y. Ketone bodies as markers for Type 1 (insulin-dependent) diabetes and their value in the monitoring of diabetes control. Diabetologia 1984;26:343-8.
- Ubukata E. Diurnal variation of blood ketone bodies in insulin-dependent diabetes mellitus and noninsulindependent diabetes mellitus patients: The relationship to serum C-peptide immunoreactivity and free insulin. Ann Nutr Metab 1990;34:333-42.
- Luzi L, Barrett EJ, Groop LC, Ferrannini E, DeFronzo RA. Metabolic effects of low-dose insulin therapy on glucose metabolism in diabetic ketoacidosis. Diabetes 1988;37:1470-77.
- Hale PJ, Crase J, Nattrass M. Metabolic effects of bicarbonate in the treatment of diabetic ketoacidosis. Br Med J 1984;289:1035-8.
- Tietz Textbook of Clinical Chemistry. Burtis CA, Ashwood ER, eds. 3rd Edition, Philadelphia, PA: W.B. Saunders Co. (1999) p. 1790-1845.
- National Cholesterol Education Program. ATP III Guidelines At-A-Glance Quick Desk Reference. NIH Publication No. 01-3305. May 2001.

This product(s) and/or its use are protected by one or more of the following patents: US5,509,410; US5,628,890; US5,727,548; US6,129,823; US6,736,957; US6,764,581; US6,939,450; US7,905,999; US6,377,894; US6,600,997; US6,773,671; US5,682,884; US6,591,125; US7,058,437; US7,504,019; US7,740,581; US7,922,883; US8,118,993; EP1,009,850B1; EP1,119,637B8; EP1,023,455B1; EP1135679B1; EP1801229B1; CA2302448C; CA2346415C; CA2351796C; CA2353670C and CA2305800C. Additional patents may be issued and/or pending.

Precision Xtra is a trademark of the Abbott Group of Companies in various jurisdictions.

© 2011 Abbott

ART23419 Rev. B 06/12

Manufacturer / Fabricante:

Abbott Diabetes Care Inc.  
1360 South Loop Road  
Alameda, CA 94502

### LOT

Expected Results for Use with Medisense Control Solutions / Resultados esperados para uso con las soluciones de control MediSense

**LO** Low / Bajo:

**MID** Mid / Medio:

**HI** High / Alto: