



# FreeStyle *Libre*

FLASH GLUCOSE MONITORING SYSTEM



FreeStyle LibreLink app  
A FreeStyle Libre product



## Interactive Tutorial



Next >

# Important information about the FreeStyle Libre System

## Important Safety Information

### Indications for Use

The FreeStyle Libre Flash Glucose Monitoring System is a continuous glucose monitoring (CGM) device indicated for the management of diabetes in persons age 18 and older. It is designed to replace blood glucose testing for diabetes treatment decisions.

The System detects trends and tracks patterns aiding in the detection of episodes of hyperglycemia and hypoglycemia, facilitating both acute and long-term therapy adjustments. Interpretation of the System readings should be based on the glucose trends and several sequential readings over time. The System is intended for single patient use and requires a prescription.

### Contraindications



**MRI/CT/Diathermy:** The System must be removed prior to Magnetic Resonance Imaging (MRI), Computed Tomography (CT) scan, or high-frequency electrical heat (diathermy) treatment. The effect of MRI, CT scans, or diathermy on the performance of the System has not been evaluated. The exposure may damage the Sensor and may impact proper function of the device which could cause incorrect readings.

### WARNINGS:

- **Do not ignore symptoms that may be due to low or high blood glucose:** if you are experiencing symptoms that are not consistent with your glucose readings, consult your health care professional.
- **Checking Sensor glucose readings with a blood glucose meter:** Under the following conditions, Sensor glucose readings may not be accurate and you should conduct a fingerstick test using a blood glucose meter. You should not use Sensor glucose readings to make a diabetes treatment decision:
  - If you suspect that your reading may be inaccurate for any reason
  - When you are experiencing symptoms that may be due to low or high blood glucose
  - When you are experiencing symptoms that do not match the Sensor glucose readings
  - During times of rapidly changing glucose (more than 2 mg/dL per minute), when interstitial fluid glucose levels as measured by the Sensor may not accurately reflect blood glucose levels
  - When the Sensor glucose reading does not include a Current Glucose number or Glucose Trend Arrow
  - In order to confirm hypoglycemia or impending hypoglycemia as reported by the Sensor

- When you see the  symbol, you must check your blood glucose with a blood glucose meter before making any treatment decisions. Sensor readings may not accurately reflect blood glucose levels.
- If you are using the FreeStyle LibreLink app, you must also have access to a blood glucose monitoring system as the App does not provide one.
- **Hypoglycemic unawareness:** The System has not been evaluated for use in patients with hypoglycemic unawareness and will not automatically alert you of a hypoglycemic event without you scanning your Sensor.
- **No alarms without a Sensor scan:** The System does not have alarms that will automatically notify you when you are having a severe low (hypoglycemic) or high (hyperglycemic) glucose event unless you scan your Sensor. For example, the System does not have an alarm that can alert or wake you when you are sleeping in the case of low or high glucose.
- **Choking hazard:** The System contains small parts that may be dangerous if swallowed.

## Cautions and Limitations

Below are important cautions and limitations to keep in mind so you can use the System safely. They are grouped into categories for easy reference.



### What to know about Alarms/Alerts:

- There are NO alarms or alerts unless you scan the Sensor.



### What to know before using the System:

- Review all product information before use.
- Take standard precautions for transmission of blood borne pathogens to avoid contamination.



### Who should not use the System:

- **Do not use the System in people less than 18 years of age.** The System is not approved for use in people under 18 years of age and Sensor readings in this population may be inaccurate. In general, continuous glucose monitoring systems are recognized to be less accurate in children than in adults.
- **Do not use the System in critically ill patients.** The System is not approved for use in these patients. It is not known how different conditions or medications common to the critically ill population may affect performance of the System. Sensor glucose readings may be inaccurate in critically ill patients.
- **Do not use the System in pregnant women or persons on dialysis.** The System is not approved for use in pregnant women or persons on dialysis and has not been evaluated in these populations.
- Performance of the System when used with other implanted medical devices, such as pacemakers, has not been evaluated.



### What should you know about wearing a Sensor:

- After the start-up period, the Sensor can be worn for up to the wear duration specified in your Sensor Kit's product insert.
- Some individuals may be sensitive to the adhesive that keeps the Sensor attached to the skin. If you notice significant skin irritation around or under your Sensor, remove the Sensor and stop using the System. Contact your health care professional before continuing to use the System.
- Intense exercise may cause your Sensor to loosen due to sweat or movement of the Sensor. Remove and replace your Sensor if it starts to loosen and follow the instructions to select an appropriate application site.
- Do not reuse Sensors. The Sensor and Sensor Applicator are designed for single use. Reuse may result in no glucose readings and infection. Not suitable for re-sterilization. Further exposure to irradiation may cause inaccurate results.
- If a Sensor breaks inside your body, call your health care professional.



### **How to Store the Sensor Kit:**

- Store the Sensor Kit between 39°F and 77°F. Storage outside of this range may cause inaccurate Sensor glucose readings. While you don't need to keep your Sensor Kit in a refrigerator, you can as long as the refrigerator is between 39°F and 77°F. Do not freeze.
- Store the Sensor Kit between 10-90% non-condensing humidity.




### **When not to use the System:**

- Do NOT use if the Sensor Kit package, Sensor Pack or Sensor Applicator appear to be damaged or already opened due to risk of no results and/or infection.
- Do NOT use if Sensor Kit contents are past expiration date.
- Do NOT use if the Reader appears to be damaged due to risk of electric shock and/or no results.



### **What to know before you Apply the Sensor:**

- The Sensor Pack and Sensor Applicator are packaged as a set (separately from the Reader) and have the same Sensor code. Check that the Sensor codes match before using your Sensor Pack and Sensor Applicator. Do not use Sensor Packs and Sensor Applicators with different Sensor codes together as this will result in incorrect glucose readings.
- Clean the application site and ensure that it is dry prior to Sensor insertion. This helps the Sensor stay attached to your body.
- Clean hands prior to Sensor handling/insertion to help prevent infection.
- Change the application site for the next Sensor application to prevent discomfort or skin irritation.
- Sensor placement is not approved for sites other than the back of the arm. If placed in other areas, the Sensor may not function properly.



- Select an appropriate Sensor site to help the Sensor stay attached to the body and prevent discomfort or skin irritation. Avoid areas with scars, moles, stretch marks, or lumps. Select an area of skin that generally stays flat during normal daily activities (no bending or folding). Choose a site that is at least 1 inch away from an insulin injection site.



#### **When is Sensor Glucose different from Blood Glucose:**

- Physiological differences between the interstitial fluid and capillary blood may result in differences in glucose readings between the System and results from a fingerstick test using a blood glucose meter. Differences in glucose readings between interstitial fluid and capillary blood may be observed during times of rapid change in blood glucose, such as after eating, dosing insulin, or exercising.



#### **What to know about interfering substances such as Vitamin C and Aspirin:**

- Taking ascorbic acid (vitamin C) while wearing the Sensor may falsely raise Sensor glucose readings. Taking salicylic acid (used in some pain relievers such as aspirin and some skin care products) may slightly lower Sensor glucose readings. The level of inaccuracy depends on the amount of the interfering substance active in the body.
- Test results did not indicate interference for methyldopa (used in some drugs to treat high blood pressure) or tolbutamide (infrequently used in some drugs to treat diabetes in the US) at maximum circulating levels. However, concentrations of potential interferents in interstitial fluid are unknown compared to circulating blood.



#### **What to know about X-Rays:**

- The Sensor should be removed prior to exposing it to an X-ray machine. The effect of X-rays on the performance of the System has not been evaluated. The exposure may damage the Sensor and may impact proper function of the device to detect trends and track patterns in glucose values during the wear period.



#### **When to remove the Sensor:**

- If the Sensor is becoming loose or if the Sensor tip is coming out of your skin, you may get no readings or unreliable readings, which may not match how you feel. Check to make sure your Sensor has not come loose. If it has come loose, remove it and apply a new one.
- If you believe your glucose readings are not correct or are inconsistent with how you feel, perform a blood glucose test on your finger to confirm your glucose. If the problem continues, remove the current Sensor and apply a new one.



#### **What to do if you are dehydrated:**

- Severe dehydration and excessive water loss may cause inaccurate Sensor glucose readings. If you believe you are suffering from dehydration, consult your health care professional immediately.



#### **What to know about the Reader's Built-in Meter:**

- The FreeStyle Libre Reader has a built-in blood glucose meter that is designed to be used only with FreeStyle Precision Neo blood glucose test strips and MediSense Glucose and Ketone Control Solution. Using other test strips with the Reader's built-in meter will produce an error or cause the Reader's built-in meter to not turn on or start a test. The Reader's built-in meter does not have ketone testing functionality.
- The Reader's built-in meter is not for use on people who are dehydrated, hypotensive, in shock, or for individuals in hyperglycemic-hyperosmolar state, with or without ketosis.
- The Reader's built-in meter is not for use on neonates, in critically-ill patients, or for diagnosis or screening of diabetes.
- See Using the Reader's Built-in meter section of the Reader Kit User's Manual for additional important information on the use of the Reader's built-in meter.



#### **Where to charge your Reader:**

- Be sure to select a location for charging that allows the power adapter to be easily unplugged. Do NOT block access to the charger due to the potential risk of electrical shock.



#### **What to know about FreeStyle LibreLink:**

- FreeStyle LibreLink installed on a smartphone is intended for use by a single person. It must not be used by more than one person due to the risk of misinterpreting glucose information.
- FreeStyle LibreLink and FreeStyle Libre Readers do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.



◀ Previous

Acknowledge





# Welcome to Your System!

The FreeStyle Libre System is a continuous glucose monitoring (CGM) System designed to replace blood glucose testing for diabetes treatment decisions. This Interactive Tutorial will help you learn how to set up and use your new System. Topics include:

- Product overview
- Reader setup & use
- App setup & use
- Sensor application & start up
- Product use & treatment decisions guide

Pay special attention to and throughout this tutorial. Click the icons to view important considerations about using the System.

For more details, refer to the User's Manual and the Quick Reference Guide.

[◀ Previous](#)[Next ▶](#)



## Product Overview

[System Overview](#) >

[Getting to Know the Reader](#) >

[Getting to Know FreeStyle LibreLink App](#) >

[Sensor Kit](#) >

[Sensor Glucose Readings](#) >

IMPORTANT: Safety information about the System is in the User's Manual and the Quick Reference Guide. Read all of the information in the User's Manual, the Quick Reference Guide and the FreeStyle Precision Neo blood glucose test strip instructions for use before using your System.



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

[◀ Previous](#)

[Next ▶](#)

## System Overview

The System has two main parts: a disposable Sensor and either a handheld Reader or mobile app. The Sensor is worn on the back of the upper arm and does not need to be calibrated with blood glucose values. The Reader or mobile app is used to wirelessly scan the Sensor and get glucose readings. The Reader also has a built-in blood glucose meter, which works with FreeStyle Precision Neo blood glucose test strips.

**IMPORTANT:** Safety information about the System is in the User's Manual and the Quick Reference Guide. Read all of the information in the User's Manual, the Quick Reference Guide and the FreeStyle Precision Neo blood glucose test strip instructions for use before using your System.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

Previous

Next



## Getting to Know the Reader

The Reader scans the Sensor and provides glucose readings. It can store approximately 90-days of glucose history as well as notes about activities, insulin, meals, and exercising. When scanning, the Sensor automatically transmits data to the Reader.

The Reader is compact, lightweight, and easy to hold. It has a backlit color touchscreen and uses a rechargeable battery. The Reader is not waterproof.



**USB Port**

Used to charge the Reader and connect it to a computer.

**Test Strip Port**

Insert a test strip here to use the built-in meter.

**Touch Screen**

**Home Button**

Turns the Reader on/off and takes you to the Home screen from any other screen.



◀ Previous

Next ▶

## Getting to Know FreeStyle LibreLink App

The App performs some similar functions to the Reader. You can use it to start a Sensor, get glucose readings from the Sensor, and store your glucose history and notes you enter.

FreeStyle LibreLink app is available for download from the App Store.

The App is not compatible with all smartphones. Before upgrading your smartphone or its operating system, check [www.FreeStyleLibre.com](http://www.FreeStyleLibre.com).



# Sensor Kit

The Sensor automatically measures and continuously stores glucose readings for 8 hours. The Sensor Kit has two parts: a Sensor Pack and a Sensor Applicator. Once you have assembled the Sensor, you will apply it to the back of your upper arm. It has a small, flexible, 5mm-long filament that is inserted just under the skin. After the start-up period, the Sensor can be worn for up to the wear duration specified in your Sensor Kit's product insert.

**IMPORTANT:** The Sensor is water-resistant in up to 3 feet (1 meter) of water. Do not immerse longer than 30 minutes.

**Sensor Applicator**  
Applies the Sensor to the body.



**Sensor Pack**  
Used with the Sensor Applicator to prepare the Sensor for use.



**Sensor**  
Measures your glucose while on your body (only visible after applied).




[◀ Previous](#) | [Next ▶](#)



# Sensor Glucose Readings

The Sensor Glucose Reading screen appears after you use your device to scan your Sensor. Your Reading includes your Current Glucose, a Glucose Trend Arrow indicating which way your glucose is going, and a graph of your current and stored glucose readings.

## What you see on the Reader



◀ Previous | Next ▶



# Sensor Glucose Readings

## What you see on the App

### Message

Tap for more information.

### Current Glucose

Glucose from your latest scan.

### Check Blood Glucose symbol

Tap for more information.

### Food Note



### Glucose Trend Arrow

Direction your glucose is going.

### Rapid-Acting Insulin Note

### Glucose Graph

Graph of your current and stored glucose readings.

### Add Notes

Tap to add notes to your glucose reading.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



## Reader

First Time Reader Setup >

Reader Home Screen >

Setting Reminders >

Changing the Reader Settings >

Set the **Current Date** using the arrows on the touchscreen.

Touch **next** to continue.



Set the **Current Time**. Touch **next** to continue.

Touch + or - to set your **Target Glucose Range**. Work with  your health care professional to determine your Target Glucose Range. Touch **next** to continue.

The Reader now displays important information about how to understand the **Glucose Trend Arrow**. Touch **next** to move through the next topics.

When the setup is complete, touch **done** to return to the Home screen.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶

# First Time Reader Setup

Complete the setup to use the Reader to check your Sensor glucose readings or use the Reader's built-in meter.

## How to do it:

Press the Home Button to turn on the Reader.

If prompted, use the touchscreen to select your preferred language for the Reader. Touch **OK** to continue. 

Set the **Current Date** using the arrows on the touchscreen. Touch **next** to continue.

Set the **Current Time**. Touch **next** to continue. 

Touch + or - to set your **Target Glucose Range**. Work with your health care professional to determine your Target Glucose Range. Touch **next** to continue. 

The Reader now displays important information about how to understand the **Glucose Trend Arrow**. Touch **next** to move through the next topics.

When the setup is complete, touch **done** to return to the Home screen.



 Previous | Next 



## First Time Reader

Complete the setup to use the Reader to take glucose readings or use the Reader to view glucose readings.

### How to do it:

Press the Home Button to turn on the Reader.

If prompted, use the touchscreen to select your preferred language for the Reader. Touch **OK** to continue.



Set the **Current Date** using the arrows on the touchscreen. Touch **next** to continue.

Set the **Current Time**. Touch **next** to continue.



Touch + or - to set your **Target Glucose Range**. Work with your health care professional to determine your Target Glucose Range. Touch **next** to continue.



The Reader now displays important information about how to understand the **Glucose Trend Arrow**. Touch **next** to move through the next topics.

When the setup is complete, touch **done** to return to the Home screen.

### Note

Use the pad of your finger. Do NOT use your fingernail or any other object on the screen.





# First Time Reader

Complete the setup to use the Reader to view your glucose readings or use the Reader to set the time and date correctly. These values affect the Reader data and settings.

## How to do it:

Press the Home Button to turn on the Reader.

If prompted, use the touchscreen to select your preferred language for the Reader. Touch **OK** to continue.

Set the **Current Date** using the arrows on the touchscreen. Touch **next** to continue.

Set the **Current Time**. Touch **next** to continue.

Touch + or - to set your **Target Glucose Range**. Work with your health care professional to determine your Target Glucose Range. Touch **next** to continue.

The Reader now displays important information about how to understand the **Glucose Trend Arrow**. Touch **next** to move through the next topics.

When the setup is complete, touch **done** to return to the Home screen.

## Caution



It is very important to set the time and date correctly. These values affect the Reader data and settings.





## First Time Reader

Complete the setup to use the Reader. Your Target Glucose Range is displayed on glucose graphs on the Reader and used to calculate your Time In Target.

### How to do it:

Press the Home Button to turn on the Reader.

If prompted, use the touchscreen to select your preferred language for the Reader. Touch **OK** to continue.

Set the **Current Date** using the arrows on the touchscreen. Touch **next** to continue.


Set the **Current Time**. Touch **next** to continue.

Touch + or - to set your **Target Glucose Range**. Work with your health care professional to determine your Target Glucose Range. Touch **next** to continue.

The Reader now displays important information about how to understand the **Glucose Trend Arrow**. Touch **next** to move through the next topics.

When the setup is complete, touch **done** to return to the Home screen.

### Note




◀ Previous

Next ▶



# Reader Home Screen

The Reader Home Screen provides access to information about glucose readings and the System. Press the Home Button to go to the Home Screen from any other screen.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶

# Setting Reminders

Use Reminders to remember when to check glucose, take insulin, or as a general alarm.

## How to do it:

Touch the symbol on the Home Screen.

Touch to select which **Type** of reminder to set: Check Glucose, Take Insulin, or Alarm.

Touch to select how often the Reminder needs to **Repeat**: Daily, Once, or Timer.

Select the Reminder **Time**, using the arrows on the touchscreen.

Touch **save** to save all Reminder settings.

From the Reminders screen, you can turn the reminder **On/Off** or **add new** reminders.


Touch **done** to return to the Home Screen.

When reminders are On, the next reminder time appears next to the symbol on the Home Screen. For example,

8:30am



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM



◀ Previous

Next ▶



# Setting Reminders

Use Reminders to remember when to check glucose, take insulin, or as a general alarm.

## Note

You can set reminders for a specific time (e.g. 8:30 am) or as a timer (e.g. 3 hours from the current time).



### How to do it:

Touch the symbol on the Home Screen.

Touch to select which **Type** of reminder to set: Check Glucose, Take Insulin, or Alarm.



Touch to select how often the Reminder needs to **Repeat**: Daily, Once, or Timer.


Select the Reminder **Time**, using the arrows on the touchscreen.

Touch **save** to save all Reminder settings.

From the Reminders screen, you can turn the reminder **On/Off** or **add new** reminders.

Touch **done** to return to the Home Screen.

When reminders are On, the next reminder time appears next to the symbol on the Home Screen. For example, 8:30am



◀ Previous

Next ▶

# Changing the Reader Settings

Many System features can be customized from the Settings menu.

## How to do it:

Touch the Settings Symbol  on the Home Screen.

Touch the arrows to scroll up or down. Touch the setting you want to change: Sounds, Target Range, Control Solution Test, Time & Date, Language, System Status, Reader Basics and Dose Increment. Touch **OK** to save.

Touch **System Status** to view System information including current Sensor end date, serial numbers for Sensor and Reader, number of Sensors used, number of tests performed using test strips, Event Logs which may be used by Customer Service to help troubleshoot your System, and perform a Reader Test.

Touch **Reader Basics** to access important information about the Reader.



 Previous | Next 



## App

First Time App Setup



App Home Screen



Setting Reminders



Changing App Settings



Confirm your country and tap **NEXT**.

You need a LibreView account to use the App. Follow onscreen instructions to review legal information and create a new account or login to your existing account.

LibreView Data Management Software is developed and distributed by Newyu, Inc. Use of FreeStyle LibreLink requires registration with LibreView, a service provided by Abbott and Newyu, Inc.

Confirm your glucose unit of measure and tap **NEXT**.

Set your Target Glucose Range and tap **NEXT**. Work with  your health care professional to determine your Target Glucose Range.

Select how you count carbohydrates and tap **NEXT**. 

The App now displays some important information. Tap **NEXT** to move through the screens.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

# LibreLink



 Abbott



Already have a LibreView account? [Sign In](#)

**GET STARTED NOW**

< Previous

Next >

# First Time App Setup

## How to do it:

Check that your smartphone is connected to a network (WiFi or cellular). Install FreeStyle LibreLink from the App store and open the App.

Swipe left to view some helpful tips or tap **GET STARTED NOW**.

Confirm your country and tap **NEXT**.

You need a LibreView account to use the App. Follow onscreen instructions to review legal information and create a new account or login to your existing account.


LibreView Data Management Software is developed and distributed by Newyu, Inc. Use of FreeStyle LibreLink requires registration with LibreView, a service provided by Abbott and Newyu, Inc.

Confirm your glucose unit of measure and tap **NEXT**.

Set your Target Glucose Range and tap **NEXT**. Work with your health care professional to determine your Target Glucose Range.

Select how you count carbohydrates and tap **NEXT**.

The App now displays some important information. Tap **NEXT** to move through the screens.



< Previous

Next >



# First Time App Setup

## How to do it:

Check that your smartphone is connected to Wi-Fi (or cellular). Install FreeStyle LibreLink from the App store and open the App.

Swipe left to view some helpful tips or tap **GET STARTED NOW**.

Confirm your country and tap **NEXT**.

You need a LibreView account to use the App. Follow onscreen instructions to review legal information and create a new account or login to your existing account.

LibreView Data Management Software is developed and distributed by Newyu, Inc. Use of FreeStyle LibreLink requires registration with LibreView, a service provided by Abbott and Newyu, Inc.

Confirm your glucose unit of measure and tap **NEXT**.

Set your Target Glucose Range and tap **NEXT**. Work with your health care professional to determine your Target Glucose Range.

Select how you count carbohydrates and tap **NEXT**.


The App now displays some important information. Tap **NEXT** to move through the screens.

## Note

Your Target Glucose Range is displayed on glucose graphs in the App and used to calculate your Time in Target.



  
**FreeStyle  
LibreLink**



 **Abbott**



Already have a LibreView account? [Sign In](#)

**GET STARTED NOW**

  
**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

 Previous

Next 



# First Time App Setup

## How to do it:

Check that your smartphone is connected to Wi-Fi (or cellular). Install FreeStyle LibreLink from the App store and open the App.

Swipe left to view some helpful tips or tap **GET STARTED NOW**.

Confirm your country and tap **NEXT**.

You need a LibreView account to use the App. Follow onscreen instructions to review legal information and create a new account or login to your existing account.

LibreView Data Management Software is developed and distributed by Newyu, Inc. Use of FreeStyle LibreLink requires registration with LibreView, a service provided by Abbott and Newyu, Inc.

Confirm your glucose unit of measure and tap **NEXT**.

Set your Target Glucose Range and tap **NEXT**. Work with your health care professional to determine your Target Glucose Range.

Select how you count carbohydrates and tap **NEXT**.


The App now displays some important information. Tap **NEXT** to move through the screens.

## Note

The carbohydrate unit will be used in any food notes you enter in the App.



FreeStyle  
LibreLink



Abbott



Already have a LibreView account? [Sign In](#)

**GET STARTED NOW**



FreeStyle  
Libre

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# App Home Screen

The App Home Screen provides access to information about glucose readings and the App. To return to the Home Screen from another screen, go to the Main Menu and tap **Home**.

## Main Menu

Tap to access the Home Screen, Logbook, other history options, and the Share option. You can also access Settings, Help, and other information.

## Scan Button

Tap this symbol or the button at the bottom of the screen when you're ready to scan your Sensor.

## Glucose Graph


Graph of your stored Sensor glucose readings.

## Glucose Information

Your Time In Target, information about your last scan, and average glucose for the last 24 hours.

## Scan Button

Tap this button or the symbol at the top of the screen when you're ready to scan your Sensor.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Setting Reminders

Use Reminders to remember when to check glucose, take insulin, or as a general alarm.

## How to do it:

Go to the Main Menu and tap Reminders. Tap **ADD REMINDER**.

Name your reminder.

Tap the time fields to set the time for the reminder.



Tap **DONE**. You will now see your reminder on the list along with the time you will receive it.

### Add Reminder

Reminder Name  
Exercise

2	57
3	58
4	59 AM
5	00 PM
6	01
7	02
8	03

Repeating

<input type="checkbox"/> All	<input type="checkbox"/> Sunday
<input checked="" type="checkbox"/> Monday	<input type="checkbox"/> Tuesday
<input checked="" type="checkbox"/> Wednesday	<input type="checkbox"/> Thursday
<input checked="" type="checkbox"/> Friday	<input type="checkbox"/> Saturday

**CANCEL** **DONE**



**FreeStyle**  
**Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Setting Reminders

Use Reminders to remember when to take insulin, or as a general alarm.

## How to do it:

Go to the Main Menu and tap Reminders. Tap ADD REMINDER.

Name your reminder.

Tap the time fields to set the time for the reminder.

Tap **DONE**. You will now see your reminder on the list along with the time you will receive it.

### Note



If you would like the reminder to repeat, tap the slider to the right.

You can also select which days you would like to receive the reminder.

Reminder

1	57
2	58
3	59 AM
4	00 PM
5	01
6	02
7	03
8	04

Repeating

<input type="checkbox"/> All	<input type="checkbox"/> Sunday
<input checked="" type="checkbox"/> Monday	<input type="checkbox"/> Tuesday
<input checked="" type="checkbox"/> Wednesday	<input type="checkbox"/> Thursday
<input checked="" type="checkbox"/> Friday	<input type="checkbox"/> Saturday

**CANCEL** **DONE**



FreeStyle  
**Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Setting Reminders

Use Reminders to remember when to scan your Sensor, to remind you to take insulin, or as a general alarm.

## How to do it:

Go to the Main Menu and tap **REMINDER**.

Name your reminder.

Tap the time fields to set the time.

Tap **DONE**. You will now see your reminder on the list along with the time you will receive it.

### Note



- There is one default reminder to help you remember to scan your Sensor. This Scan Sensor reminder can be changed or disabled but cannot be deleted.
- To turn off a reminder, tap the slider to the left.
- To delete a reminder, swipe the reminder and tap the symbol. The Scan Sensor reminder cannot be deleted.
- Your reminders will be received as notifications that you can swipe or tap to dismiss.

Reminder

AM  
PM

6	01
7	02
8	03

Repeating

<input type="checkbox"/> All	<input type="checkbox"/> Sunday
<input checked="" type="checkbox"/> Monday	<input type="checkbox"/> Tuesday
<input checked="" type="checkbox"/> Wednesday	<input type="checkbox"/> Thursday
<input checked="" type="checkbox"/> Friday	<input type="checkbox"/> Saturday

**CANCEL** **DONE**



FreeStyle  
**Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Changing App Settings

## How to do it:

---

Go to the Main Menu to access the App settings.

---

**Unit of Measurement** - View the glucose unit of measure used in the App.

**Target Glucose Range** - Set the target range you want to be displayed on the App glucose graph. It is also used to calculate your Time In Target. Tap **SAVE** when you are done.

**Carbohydrate Units** - Choose grams or servings for food notes that you enter. Tap **SAVE** when you are done.

**Text to Speech** - Turn on Text to Speech to have the glucose reading read aloud when you scan the Sensor. You will hear only your current glucose value and trend arrow direction along with whether you need to do a blood glucose test before making treatment decisions. Additional information, such as the glucose graph and any message, is available on your My Glucose screen. Always review your My Glucose screen to get complete information. Tap **SAVE** when you are done.

---



◀ Previous

Next ▶



## Sensor

- [Application Site Selection >](#)
- [Preparing the Application Site >](#)
- [Preparing the Sensor Applicator >](#)
- [Applying Your Sensor >](#)
- [Starting Your Sensor >](#)
- [Removing Your Sensor >](#)

cm) away from an insulin injection site. To prevent discomfort or skin irritation, you should select a different site other than the one most recently used.



Only apply Sensor to the back of your upper arm.



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶




## Application Site Selection

Apply Sensors only on the back of your upper arm. If placed in other areas, the Sensor may not function properly and could give inaccurate readings. The application of the Sensor is not approved for other sites. Avoid areas with scars, moles, stretch marks, or lumps.



Select an area of skin that generally stays flat during your normal daily activities (no bending or folding). Choose a site that is at least 1 inch (2.5 cm) away from an insulin injection site. To prevent discomfort or skin irritation, you should select a different site other than the one most recently used.



Only apply Sensor to the back of your upper arm.



◀ Previous | Next ▶



## Application Site


Apply Sensors only on the back of your upper arm. If placed in other areas, the Sensor may not function properly and could affect readings. The application of the Sensor is not approved for other sites. Avoid areas with hair, moles, stretch marks, or lumps.

Select an area of skin that generally stays flat during your normal daily activities (no bending or folding). Choose a site that is at least 1 inch (2.5 cm) away from an insulin injection site. To prevent discomfort or skin irritation, you should select a different site other than the one most recently used.

### Caution



Intense exercise may cause your Sensor to loosen due to sweat or movement of the Sensor. Remove and replace your Sensor if it starts to loosen and follow the instructions to select an appropriate application site.



Only apply Sensor to the back of your upper arm.




◀ Previous | Next ▶



## Preparing the Application Site

Clean application site with an alcohol wipe and allow site to dry before proceeding. This helps the Sensor stay attached to your body.

[Previous](#)[Next](#)




## Preparing the Application Site

Clean application site with an alcohol swab. Allow site to dry before proceeding. Make sure the Sensor stay attached to your skin.

### Note



The area **MUST** be clean and dry, or the Sensor may not stick to the site.



◀ Previous

Next ▶



# Preparing the Sensor Applicator

To ensure glucose readings are accurate, make certain the Sensor Pack and Sensor Applicator codes match.

## How to do it:



Open the Sensor Pack by peeling the lid off completely.



Unscrew the cap from the Sensor Applicator and set the cap aside.

Place the Sensor Pack on a flat hard surface and line up the dark mark on the Sensor Applicator with the dark mark on the Sensor Pack. Press firmly down on the Sensor Applicator until it comes to a stop.

Lift the Sensor Applicator out of the Sensor Pack.




The Sensor Applicator is prepared and ready to apply the Sensor.



**FreeStyle**  
**Libre**

FLASH GLUCOSE MONITORING SYSTEM



◀ Previous

Next ▶



## Preparing the Sensor

To ensure glucose readings are accurate, make sure the Sensor Pack and Sensor Applicator codes match.

### How to do it:

Open the Sensor Pack by peeling back the seal completely.

Unscrew the cap from the Sensor Applicator and set the cap aside.

Place the Sensor Pack on a flat hard surface and line up the dark mark on the Sensor Applicator with the dark mark on the Sensor Pack. Press firmly down on the Sensor Applicator until it comes to a stop.


Lift the Sensor Applicator out of the Sensor Pack.

The Sensor Applicator is prepared and ready to apply the Sensor.

### Caution



The Sensor Pack and Sensor Applicator are packaged as a set (separately from the Reader) and have the same Sensor code. Check that the Sensor codes match before using the Sensor Pack and Sensor Applicator. Do not use Sensor Packs and Sensor Applicators with different Sensor codes together as this will result in incorrect glucose readings.



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



## Preparing the Sensor

To ensure glucose readings are accurate, make sure to correctly assemble the Sensor Pack and Sensor Applicator. Make sure the Sensor Pack and Sensor Applicator lot numbers match.

### How to do it:

Open the Sensor Pack by peeling the lid off completely.

Unscrew the cap from the Sensor Applicator and set the cap aside.

Place the Sensor Pack on a flat hard surface and line up the dark mark on the Sensor Applicator with the dark mark on the Sensor Pack. Press firmly down on the Sensor Applicator until it comes to a stop.


Lift the Sensor Applicator out of the Sensor Pack.

The Sensor Applicator is prepared and ready to apply the Sensor.

### Caution



Do NOT use if the Sensor Pack or the Sensor Applicator seem to be damaged or already opened. Do NOT use if past expiration date.



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



## Preparing the Sensor

To ensure glucose readings are accurate, make sure the Sensor Pack and Sensor Applicator are matched.

### Caution



The Sensor Applicator now contains a needle. Do NOT touch inside the Sensor Applicator or put it back into the Sensor Pack.

#### How to do it:

Open the Sensor Pack by peeling the lid off completely.



Unscrew the cap from the Sensor Applicator and set the cap aside.




Place the Sensor Pack on a flat hard surface and line up the dark mark on the Sensor Applicator with the dark mark on the Sensor Pack. Press firmly down on the Sensor Applicator until it comes to a stop.

Lift the Sensor Applicator out of the Sensor Pack.



The Sensor Applicator is prepared and ready to apply the Sensor.



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶

# Applying Your Sensor

## How to do it:

Place the Sensor Applicator over the prepared site and push down firmly to apply the Sensor to your body.



Gently pull the Sensor Applicator away from your body. The Sensor should now be attached to your skin.



Make sure the Sensor is secure after application. Put the cap back on the Sensor Applicator. Discard the used Sensor Pack and Sensor Applicator according to local regulations.



Previous | Next



# Applying Your Sensor

## How to do it:

Place the Sensor Applicator over the prepared site and push down firmly to apply the Sensor to your body.


Gently pull the Sensor Applicator away from your body. The Sensor should now be attached to your skin.

Make sure the Sensor is secure after application. Put the cap back on the Sensor Applicator. Discard the used Sensor Pack and Sensor Applicator according to local regulations.

### Caution



Do NOT push down on the Sensor Applicator until placed over prepared site to prevent unintended results or injury.



◀ Previous | Next ▶



# Applying Your Sensor

## How to do it:


Place the Sensor Applicator over the site and push down firmly to attach to your body.

Gently pull the Sensor Applicator away from your body. The Sensor should now be attached to your skin.

Make sure the Sensor is secure after application. Put the cap back on the Sensor Applicator. Discard the used Sensor Pack and Sensor Applicator according to local regulations.

### Note

Applying the Sensor may cause bruising or bleeding. If there is bleeding that does not stop, remove the Sensor, and contact your health care professional.



◀ Previous

Next ▶

# Starting Your Sensor


Start your Sensor with the Reader if you would like to use the Reader to check your glucose. If you want to also use the App with the Sensor, scan the Sensor with the App after starting it with the Reader.

## How to do it with the Reader:

Press the Home Button to turn on the Reader.

Touch Start New Sensor.

Hold the Reader within 1.5 inches (4 cm) of the Sensor to scan it. This starts your Sensor. If sounds are turned on, the Reader beeps when the Sensor has been successfully activated. The Sensor can be used to check your glucose after the start-up period. During the start-up period, you can scan the Sensor to check its status.



◀ Previous

Next ▶



# Starting Your Sensor

Start your Sensor with the Reader. You can also use the Reader to check your Sensor's status. You can also use the App with the Sensor. You can also use the App with the Sensor after starting it with the Reader.

## How to do it with the Reader:

Press the Home Button to turn on the Reader.

Touch Start New Sensor.

Hold the Reader within 1.5 inches (4 cm) of the Sensor to scan it. This starts your Sensor. If sounds are turned on, the Reader beeps when the Sensor has been successfully activated. The Sensor can be used to check your glucose after the start-up period. During the start-up period, you can scan the Sensor to check its status.




**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

## Note

- If the Sensor is not successfully scanned within 15 seconds, the Reader displays a prompt to scan the Sensor again. Touch **OK** to return to the Home Screen and touch **Start New Sensor** to scan your Sensor.
- Remember that FreeStyle LibreLink and Readers do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.



0:00 / 0:10



◀ Previous

Next ▶



# Starting Your Sensor

Start the Sensor with the App if you would like to only use the App to check your glucose. If you want to use both the Reader and the App to check your glucose, you must start the Sensor with the Reader first.

## How to do it with the App:

Tap the scan button .



Hold the top of your iPhone near the Sensor. Hold still until you hear a tone and/or feel a vibration. This completes the scan.



The Sensor can be used to check your glucose after the start-up period. While the Sensor is starting up, you can navigate away from the App.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶

# Starting Your Sensor

Start the Sensor with the App or Reader. You can either tap the blue box on the Home Screen or the scan button at the top right.

You can either tap the blue box on the Home Screen or at the top right.

## How to do it with the App:

Tap the scan button .



Hold the top of your iPhone near the Sensor. Hold still until you hear a tone and/or feel a vibration. This completes the scan.



The Sensor can be used to check your glucose after the start-up period. While the Sensor is starting up, you can navigate away from the App.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

Previous

Next



# Starting Your Sensor

Start the Sensor with the App. Open the App to check your glucose. If you have a Reader, connect the Reader and the App to check your glucose. Connect the Reader with the App to check your glucose. Connect the Reader with the App to check your glucose. Connect the Reader with the Reader first.

## How to do it with the App:

Tap the scan button .

Hold the top of your iPhone near the Sensor. Hold still until you hear a tone and/or feel a vibration. This completes the scan.

The Sensor can be used to check your glucose after the start-up period. While the Sensor is starting up, you can navigate away from the App.

### Note



- If your Sensor is not successfully scanned, you may receive a Scan Error. Tap the scan button and scan again.
- Remember that FreeStyle LibreLink and Readers do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶


# Removing Your Sensor

The Sensor automatically stops working after the wear duration specified in the Sensor Kit's product insert and must be replaced. Replace the Sensor if you notice any irritation or discomfort at the application site or if your device reports a problem with the Sensor currently in use. Taking action early can keep small problems from turning into larger ones.

## How to do it:

Pull up the edge of the adhesive that keeps your Sensor attached to your skin. Slowly peel away from your skin in one motion. 

Discard the used Sensor following directions from your health care professional. See the Maintenance and Disposal section of the Reader Kit User's Manual.



Previous

Next



## Removing Your Sensor

The Sensor automatically stops working after the wear duration specified in the User's Manual. The Sensor is a single-use device and must be replaced. If you notice any irritation or discomfort at the application site or if your device reports a problem with the Sensor currently in use. Taking action early can keep small problems from turning into larger ones.


### How to do it:

Pull up the edge of the adhesive that keeps your Sensor attached to your skin. Slowly peel away from your skin in one motion. 

Discard the used Sensor following directions from your health care professional. See the Maintenance and Disposal section of the Reader Kit User's Manual.

### Note

Any remaining adhesive residue on the skin can be removed with warm soapy water or isopropyl alcohol.



◀ Previous

Next ▶



## Product Use

- [Checking Glucose with a Sensor >](#)
- [Understanding Sensor Readings - Reader >](#)
- [Understanding Sensor Readings - App >](#)
- [Treatment Decisions Guide >](#)
- [Treatment Decisions Guide - Knowledge Check >](#)
- [Checking Glucose with a Test Strip >](#)
- [Adding Notes >](#)
- [Reviewing History >](#)

Play ➤



◀ Previous | Next ➤



# Checking Glucose with a Sensor

## How to do it with the Reader:

Press the Home Button to turn on the Reader or touch **Check Glucose** on the Home Screen.

Hold the Reader within 1.5 inches (4 cm) of the Sensor to scan it. The Sensor wirelessly sends glucose readings to the Reader. If sounds are turned on, the Reader beeps when the Sensor has been successfully scanned.

The Reader displays the current glucose reading along with a glucose graph and an arrow indicating the direction the glucose is going.

Play ➤



◀ Previous | Next ➤



# Checking Glucose

## How to do it with the Reader:

Press the Home Button to turn **Glucose** on the Home Screen.


Hold the Reader within 1.5 inches (4 cm) of the Sensor to scan it. The Sensor wirelessly sends glucose readings to the Reader. If sounds are turned on, the Reader beeps when the Sensor has been successfully scanned.

The Reader displays the current glucose reading along with a glucose graph and an arrow indicating the direction the glucose is going.

### Note



If the Sensor is not successfully scanned within 15 seconds, the Reader displays a prompt to scan the Sensor again. Touch **OK** to return to the Home Screen and touch **Check Glucose** to scan your Sensor.



Play ➤



◀ Previous

Next ➤



# Checking Glucose with a Sensor

## How to do it with the Reader:

Press the Home Button to turn on the Reader or touch **Check Glucose** on the Home Screen.

Hold the Reader within 1.5 inches (4 cm) of the Sensor to scan it. The Sensor wirelessly sends glucose readings to the Reader. If sounds are turned on, the Reader beeps when the Sensor has been successfully scanned.

The Reader displays the current glucose reading along with a glucose graph and an arrow indicating the direction the glucose is going.

Play ➤



FreeStyle Libre



◀ Previous | Next ➤



# Checking Glucose with a Sensor

## How to do it with the App:

Tap the scan button  .

Hold the top of your iPhone near the Sensor. Hold still  until you hear a tone and/or feel a vibration. This completes the scan.

The App displays the current glucose reading along with a glucose graph and an arrow indicating the direction the glucose is going.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Checking Glucose

## How to do it with the App:

Tap the scan button .

Hold the top of your iPhone near the Sensor. Hold still until you hear a tone and/or feel a vibration. This completes the scan.

The App displays the current glucose reading along with a glucose graph and an arrow indicating the direction the glucose is going.

### Note

If your Sensor is not successfully scanned, you may receive a Scan Error. Tap the scan button and scan again.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous


Next ▶

# Understanding Sensor Readings - Reader

[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

**Note:** While Sensor glucose readings are gathered in the System range of 40-500 mg/dL, the graph display range is 0-350 mg/dL for ease of review on screen. Glucose readings above 350 mg/dL are displayed at 350 mg/dL. For sequential readings above 350 mg/dL, a line is displayed at 350 mg/dL.

**Note:** The symbol may appear, indicating the Reader time was changed. Gaps in the graph may result or glucose readings may be hidden.



**FreeStyle**  
**Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[Previous](#)[Next](#)



# Understanding Sensor Readings - Reader

[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

The **Glucose Trend Arrow** gives an indication of the direction your glucose is going.

**Note:** The **Glucose Trend Arrow** may not always appear with your reading. When there is no Glucose Trend Arrow, you will see a  symbol, telling you to do a blood glucose test before making treatment decisions.



**Glucose is rising quickly** (more than 2 mg/dL per minute)



**Glucose is rising** (between 1 and 2 mg/dL per minute)



**Glucose is changing slowly** (less than 1 mg/dL per minute)



**Glucose is falling** (between 1 and 2 mg/dL per minute)



**Glucose is falling quickly** (more than 2 mg/dL per minute)

[Previous](#)[Next](#)




# Understanding Sensor Readings - Reader

[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

If your glucose is projected to be higher than 240 mg/dL or lower than 70 mg/dL within 15 minutes, you will see a or a message on the screen. You can touch the message button for more information and set a reminder to check your glucose.

When there is a message, you will see a symbol, telling you to do a blood glucose test before making treatment decisions.

**Note:** If you are not sure about a message or reading, contact your health care professional before you do anything.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[Previous](#)[Next](#)





# Understanding Sensor Readings - Reader

[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

If your glucose is higher than 240 mg/dL or lower than 70 mg/dL, you will see a **High Glucose** Or **Low Glucose** message on the screen. You can touch the message button for more information and set a reminder to check your glucose.

When there is a **Low Glucose** message, you will see a symbol, telling you to do a blood glucose test before making treatment decisions.

**Note:** If you are not sure about a message or reading, contact your health care professional before you do anything.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[◀ Previous](#)[Next ▶](#)

# Understanding Sensor Readings - Reader





[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

If LO appears on the Reader, your reading is lower than 40 mg/dL.

If HI appears on the Reader, your reading is higher than 500 mg/dL.

You can touch the message button for more information. If you get a LO or HI result on the Reader you will see a  symbol, telling you to do a blood glucose test before making treatment decisions.

**Note:** If you are not sure about a message or reading, contact your health care professional before you do anything.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM


[Previous](#)[Next](#)

# Understanding Sensor Readings - App

[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

**Note:** The glucose graph in the App will scale to 500 mg/dL to accommodate glucose readings above 350 mg/dL.

**Note:** The symbol may appear, indicating the smartphone's time was changed. Gaps in the graph may result or glucose readings may be hidden.



**Note:** Your current glucose value determines the background color on the My Glucose screen:

Orange High glucose (above 240 mg/dL)

Yellow Between the Target Glucose Range and high or low glucose level

Green Within the Target Glucose Range

Red Low glucose (below 70 mg/dL)



FreeStyle  
*Libre*  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Understanding Sensor Readings - App

[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

The **Glucose Trend Arrow** gives an indication of the direction your glucose is going.

**Note:** The **Glucose Trend Arrow** may not always appear with your reading. When there is no Glucose Trend Arrow, you will see a  symbol, telling you to do a blood glucose test before making treatment decisions.



**Glucose is rising quickly** (more than 2 mg/dL per minute)



**Glucose is rising** (between 1 and 2 mg/dL per minute)



**Glucose is changing slowly** (less than 1 mg/dL per minute)



**Glucose is falling** (between 1 and 2 mg/dL per minute)



**Glucose is falling quickly** (more than 2 mg/dL per minute)

[Previous](#)[Next](#)



# Understanding Sensor Readings - App


[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

If your glucose is projected to be higher than 240 mg/dL or lower than 70 mg/dL within 15 minutes, you will see a

**⚠ GLUCOSE GOING HIGH** or a **⚠ GLUCOSE GOING LOW** message on the screen. You can touch the **!** symbol for more information and set a reminder to check your glucose.

When there is a **⚠ GLUCOSE GOING LOW** message, you will see a symbol, telling you to do a blood glucose test before making treatment decisions.

**Note:** If you are not sure about a message or reading, contact your health care professional before you do anything.



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

[Previous](#)[Next](#)




# Understanding Sensor Readings - App

[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

If your glucose is higher than 240 mg/dL or lower than 70 mg/dL, you will see a **HIGH GLUCOSE** or **LOW GLUCOSE** message on the screen. You can touch the  symbol for more information and set a reminder to check your glucose.

When there is a **LOW GLUCOSE** message, you will see a  symbol, telling you to do a blood glucose test before making treatment decisions.

**Note:** If you are not sure about a message or reading, contact your health care professional before you do anything.

[Previous](#)[Next](#)



# Understanding Sensor Readings - App

[Sensor Glucose Reading](#)[Trend Arrows](#)[Going High/Low Reading](#)[High/Low Reading](#)[HI/LO Reading](#)

If LO appears, your reading is lower than 40 mg/dL.

If HI appears, your reading is higher than 500 mg/dL.

You can touch the symbol for more information. If you get a LO or HI result you will see a symbol, telling you to do a blood glucose test before making treatment decisions.

**Note:** If you are not sure about a message or reading, contact your healthcare professional before you do anything.




[◀ Previous](#) | [Next ▶](#)




# Treatment Decisions Guide

Using Sensor Glucose Readings for treatment decisions

## Reader



## App



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶

## Treatment Decisions Guide

**WARNING:** The System can replace blood glucose testing except in a few situations. These are the times when you need to do a blood glucose test before deciding what to do or what treatment decision to make as Sensor readings may not accurately reflect blood glucose levels.



**Do a blood glucose test** if you see the Check Blood Glucose symbol. The symbol means your Sensor glucose reading may not be accurate. For example, there may be times when you get a low glucose reading but you do not actually have low glucose.



**Do a blood glucose test** if you think your glucose readings are not correct or do not match how you feel. Do not ignore symptoms that may be due to low or high glucose. **Note:** The symbol will **NOT** appear in this situation.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶





# Treatment Decisions Guide – Example Scenarios

[When You Wake-Up](#)[Before Breakfast](#)[After Breakfast](#)[Lunch](#)[After Exercising](#)[Before Dinner](#)[After Dinner](#)**What you see:**

When you wake-up, your current glucose is 65 mg/dL and the trend arrow shows it is changing slowly → . There is a Low Glucose (Reader) / LOW GLUCOSE (App) message at the top of the screen and the symbol.

**What it means:**

Anytime you see the symbol, you should do a blood glucose test before deciding what to do.

**Reader****App**

**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[Previous](#)[Next](#)



# Treatment Decisions Guide – Example Scenarios

[When You Wake-Up](#)[Before Breakfast](#)[After Breakfast](#)[Lunch](#)[After Exercising](#)[Before Dinner](#)[After Dinner](#)


**What you see:** Before breakfast, your current glucose is 115 mg/dL. The graph shows that your glucose is going up and so does the trend arrow ↗.

**What it means:**

Consider what might be causing your glucose to go up and what you might do to prevent a high glucose.

*For example:*

- How much insulin should you take before your meal?
- Since you see ↗, should you consider taking a little more insulin?



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[◀ Previous](#)[Next ▶](#)





# Treatment Decisions Guide – Example Scenarios

[When You Wake-Up](#)[Before Breakfast](#)[After Breakfast](#)[Lunch](#)[After Exercising](#)[Before Dinner](#)[After Dinner](#)

**What you see:** After breakfast, your current glucose is 108 mg/dL. The trend arrow shows it is going down quickly ↓. There is a Glucose Going Low (Reader) / GLUCOSE GOING LOW (App) message at the top of the screen and the symbol.

**What it means:**

Anytime you see the symbol, you should do a blood glucose test before deciding what to do.

**Reader****App**

**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM


[Previous](#)[Next](#)



# Treatment Decisions Guide – Example Scenarios

[When You Wake-Up](#)[Before Breakfast](#)[After Breakfast](#)[Lunch](#)[After Exercising](#)[Before Dinner](#)[After Dinner](#)**What you see:**

When you checked your glucose before lunch, it was 90 mg/dL and rising. Before eating lunch, you took enough insulin to cover the meal and a little more since your trend arrow was ↗.


*Before Lunch**Reader**App***What it means:**

**Don't take a correction dose within 2 hours of your meal dose.** This may result in "insulin stacking" and low glucose. Consider what might be causing your glucose to go up and what you might do to prevent a high glucose.

*For example:*

- Has the insulin you took for your meal reached its full effect?
- Scan your Sensor again later.

90 minutes after lunch, your current glucose is 225 mg/dL. The graph shows that your glucose is still going up, and so does the trend arrow ↗.

*After Lunch**Reader**App*

**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

[◀ Previous](#)[Next ▶](#)



# Treatment Decisions Guide – Example Scenarios


[When You Wake-Up](#)[Before Breakfast](#)[After Breakfast](#)[Lunch](#)[After Exercising](#)[Before Dinner](#)[After Dinner](#)

**What you see:** After exercising, you are feeling shaky, sweaty, and dizzy – symptoms you generally get when you have low glucose. But, your current glucose is 204 mg/dL.

**What it means:**

Anytime you get a reading that doesn't match how you feel, do a blood glucose test.

**Note:** The Check Blood Glucose  symbol will **NOT** appear in this situation.

**Reader****App**

**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[Previous](#)[Next](#)

# Treatment Decisions Guide – Example Scenarios

[When You Wake-Up](#)[Before Breakfast](#)[After Breakfast](#)[Lunch](#)[After Exercising](#)[Before Dinner](#)[After Dinner](#)

**What you see:** Before dinner, your current glucose is 134 mg/dL. The graph shows that your glucose is going down and so does the trend arrow ↘.

**What it means:**

Consider what might be causing your glucose to go down and what you might do to prevent a low glucose.

*For example:*

- How much insulin should you take before your meal?
- Since you see ↘, should you consider taking a little less insulin?



**FreeStyle**  
**Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[◀ Previous](#)[Next ▶](#)




# Treatment Decisions Guide – Example Scenarios

[When You Wake-Up](#)[Before Breakfast](#)[After Breakfast](#)[Lunch](#)[After Exercising](#)[Before Dinner](#)[After Dinner](#)

**What you see:** After dinner, your current glucose is 215 mg/dL but there is no trend arrow. There is also the  symbol on the screen.

**What it means:**

Anytime you see the  symbol, you should do a blood glucose test before deciding what to do.

**Reader****App**

**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[◀ Previous](#)[Next ▶](#)



## Treatment Decisions Guide – Knowledge Check

1. What would you do if you scanned your Sensor and saw this  symbol with your reading?

- Do not treat based on this reading – check my blood glucose with a test strip
- Make a treatment decision



◀ Previous

Next ▶



## Treatment Decisions Guide – Knowledge Check

1. What would you do if you scanned your Sensor and saw this symbol with your reading?

Do not treat based on this reading – check my blood glucose with a test strip

Make a treatment decision

Incorrect. Don't make treatment decisions based on a Sensor glucose reading that includes this symbol. The symbol means your Sensor glucose reading may not be accurate.

Please click [Previous](#) and try again.



[◀ Previous](#)

[Next ▶](#)



## Treatment Decisions Guide – Knowledge Check

1. What would you do if you scanned your Sensor and saw this  symbol with your reading?



Do not treat based on this reading – check my blood glucose with a test strip



Make a treatment decision

Correct! Anytime you see this  symbol, check your blood glucose with a test strip before making treatment decisions. The  symbol means your Sensor glucose reading may not be accurate.



**FreeStyle**  
**Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



## Treatment Decisions Guide – Knowledge Check

2. What would you do if you are getting ready to eat lunch, you scan your Sensor and you get a reading that your glucose level is above your target range and changing slowly → .

- Do not treat based on this reading – check my blood glucose with a test strip
  
- Make a treatment decision, such as take insulin

[◀ Previous](#)[Next ▶](#)



## Treatment Decisions Guide – Knowledge Check

2. What would you do if you are getting ready to eat lunch, you scan your Sensor and you get a reading that your glucose level is above your target range and changing slowly → .



Do not treat based on this reading – check my blood glucose with a test strip

Make a treatment decision, such as take insulin

Incorrect. You could check your blood glucose with a test strip, but you don't have to, unless you have symptoms that don't match your reading.

Please click [Previous](#) and try again.



[◀ Previous](#)

[Next ▶](#)



## Treatment Decisions Guide – Knowledge Check

2. What would you do if you are getting ready to eat lunch, you scan your Sensor and you get a reading that your glucose level is above your target range and changing slowly → .

Do not treat based on this reading – check my blood glucose with a test strip



Make a treatment decision, such as take insulin

Correct! Based on your reading and what you are about to eat, consider how much insulin you should take.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

Previous

Next



## Treatment Decisions Guide – Knowledge Check

3. What would you do if your body was telling you that your glucose was low, but when you scanned your Sensor your reading was high?

- Do not treat based on this reading – check my blood glucose with a test strip
  
- Do something to lower glucose, such as take insulin

[◀ Previous](#)[Next ▶](#)



## Treatment Decisions Guide – Knowledge Check

3. What would you do if your body was telling you that your glucose was low, but when you scanned your Sensor your reading was high?

Do not treat based on this reading – check my blood glucose with a test strip



Do something to lower glucose, such as take insulin

Incorrect. If you feel low, trust your symptoms. Check your blood glucose with a test strip and treat accordingly. The Sensor may be inaccurate sometimes. Remember the  symbol will NOT appear in this situation.

Please click [Previous](#) and try again.



[◀ Previous](#)

[Next ▶](#)



## Treatment Decisions Guide – Knowledge Check

3. What would you do if your body was telling you that your glucose was low, but when you scanned your Sensor your reading was high?



Do not treat based on this reading – check my blood glucose with a test strip



Do something to lower glucose, such as take insulin

You're right! Anytime you get a reading that doesn't match how you feel, check your blood glucose with a test strip before you make a treatment decision. Remember the symbol will NOT appear in this situation.



**FreeStyle**  
**Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶

## Treatment Decisions Guide – Knowledge Check

4. “Insulin stacking” is when you take two or more doses of rapid-acting insulin too close together. Which of the following scenarios would avoid “insulin stacking”?

- You took your full breakfast insulin dose and ate breakfast. 75 minutes later your glucose was high, so you took a correction insulin dose.
  
- You took your full breakfast insulin dose and ate breakfast. 75 minutes later your glucose was high. You know insulin takes some time to work, so you decided to do nothing and scan again later.

[◀ Previous](#)[Next ▶](#)

## Treatment Decisions Guide – Knowledge Check

4. "Insulin stacking" is when you take two or more doses of rapid-acting insulin too close together. Which of the following scenarios would avoid "insulin stacking"?



- ⚡ You took your full breakfast insulin dose and ate breakfast. 75 minutes later your glucose was high, so you took a correction insulin dose.
- ⚡ You took your full breakfast insulin dose and ate breakfast. 75 minutes later your glucose was high. You know insulin takes some time to work, so you decided to do nothing and scan again later.

Incorrect. Taking multiple doses of rapid-acting insulin too close together may lead to "insulin stacking" and low glucose. Don't take a correction dose within 2 hours of your meal dose.

Please click [Previous](#) and try again.



[◀ Previous](#)

[Next ▶](#)



## Treatment Decisions Guide – Knowledge Check

4. "Insulin stacking" is when you take two or more doses of rapid-acting insulin too close together. Which of the following scenarios would avoid "insulin stacking"?

You took your full breakfast insulin dose and ate breakfast. 75 minutes later your glucose was high, so you took a correction insulin dose.



You took your full breakfast insulin dose and ate breakfast. 75 minutes later your glucose was high. You know insulin takes some time to work, so you decided to do nothing and scan again later.

That's correct! Since you did not take multiple doses of rapid-acting insulin too close together, you avoided "insulin stacking" and low glucose. Don't take a correction dose within 2 hours of your meal dose.



**FreeStyle**  
**Libre**

FLASH GLUCOSE MONITORING SYSTEM


◀ Previous

Next ▶


## Treatment Decisions Guide – Knowledge Check

5. What would you do if you scanned your Sensor before exercise, and you got this Reading?

Reader



App



- Do not treat based on this result – check my blood glucose with a test strip.
- Do something to bring my glucose up, like eat or take glucose tablets.




◀ Previous

Next ▶


## Treatment Decisions Guide – Knowledge Check

5. What would you do if you scanned your Sensor before exercise, and you got this Reading?

Reader



App



- Do not treat based on this result – check my blood glucose with a test strip.



- Do something to bring my glucose up, like eat or take glucose tablets.

Incorrect! Although your glucose reading is low, anytime you see this  symbol, check your blood glucose with a test strip before making treatment decisions. The  symbol means your Sensor glucose reading may not be accurate.

Please click [Previous](#) and try again.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[◀ Previous](#)


[Next ▶](#)




## Treatment Decisions Guide – Knowledge Check

5. What would you do if you scanned your Sensor before exercise, and you got this Reading?

Reader



App



Do not treat based on this result – check my blood glucose with a test strip.

Do something to bring my glucose up, like eat or take glucose tablets.

That's right! Anytime you see this symbol, check your blood glucose with a test strip before making treatment decisions.



**FreeStyle**  
**Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Treatment Decisions Guide – Knowledge Check

Congratulations – you have finished the knowledge assessment. Click Next to learn more about the System.

[◀ Previous](#)[Next ▶](#)



# Checking Glucose with a Test Strip

You can use the Reader's built-in meter to check your blood glucose, whether you are wearing a Sensor or not. Be sure to read the test strip instructions for use prior to using the built-in meter.

## How to do it:



Wash your hands with warm soapy water for accurate results. Thoroughly dry your hands. To warm the site, apply a warm dry pad or rub vigorously for a few seconds.

Check the FreeStyle Precision Neo test strip expiration date. Do not use expired test strips as they may give inaccurate results.

Open the foil test strip packet at the notch and tear down to remove the test strip. Use the test strip immediately after removing from the foil packet.

Insert the test strip with the three black lines at the end facing up. Push the strip in until it stops.



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Checking Glucose

You can use the Reader's built-in meter to check glucose levels. Follow the test strip instructions for use.

## How to do it:

Wash your hands with warm soapy water before testing to get accurate results. Thoroughly dry your hands with a clean cloth. Warm the site, apply a warm dry compress to the site and rub vigorously for a few seconds.

Check the FreeStyle Precision Neo Blood Glucose Test Strip's expiration date. Do not use expired strips. Expired strips may give inaccurate results.

Open the foil test strip packet. Hold the packet upright and fold it down to remove the test strip. Do not touch the test strip. Immediately after removing from the packet, insert the test strip into the Reader.

Insert the test strip with the three colored circles facing up. Push the strip firmly into the Reader until it is seated.

## Caution



Test on your fingers in accordance with the Intended Use.

### Intended Use

The FreeStyle Libre Reader's built-in meter is for use outside the body only (in vitro diagnostic use) in the quantitative measurement of glucose in fresh whole blood for self testing by lay users from the fingers. It is not intended to be used for testing neonatal blood samples or for the diagnosis or screening of diabetes.

The FreeStyle Libre Reader's built-in meter is indicated for the home (lay) user in the management of patients with diabetes. It is intended to be used by a single person and should not be shared.

The FreeStyle Precision Neo Blood Glucose Test Strips are for use with the FreeStyle Libre Reader's built-in meter to quantitatively measure glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips.

Sensor or not. Be sure to read



◀ Previous

Next ▶



# Checking Glucose

You can use the Reader's built-in sensor to check your glucose level. Follow the test strip instructions for use.

## How to do it:

Wash your hands with warm soapy water for accurate results. Thoroughly dry your hands. To warm the site, apply a warm dry pad or rub vigorously for a few seconds.

Check the FreeStyle Precision Neo test strip expiration date. Do not use expired test strips as they may give inaccurate results.

Open the foil test strip packet at the notch and tear down to remove the test strip. Use the test strip immediately after removing from the foil packet.

Insert the test strip with the three black lines at the end facing up. Push the strip in until it stops.

## Note



Do not use lotion or cream on the test site. Avoid moles, veins, bones, and tendons. Bruising may occur at the test site. If you get a bruise, consider selecting another site.



◀ Previous

Next ▶



# Checking Glucose

You can use the Reader's built-in meter to check your glucose level. Follow the test strip instructions for use.

## How to do it:

Wash your hands with warm soapy water for accurate results. Thoroughly dry your hands. To warm the site, apply a warm dry pad or rub vigorously for a few seconds.

Check the FreeStyle Precision Neo test strip expiration date. Do not use expired test strips as they may give inaccurate results.

Open the foil test strip packet at the notch and tear down to remove the test strip. Use the test strip immediately after removing from the foil packet.

Insert the test strip with the three black lines at the end facing up. Push the strip in until it stops.

## Note



The Reader's built-in meter turns off after 2 minutes of inactivity. If you don't see a reading after 2 minutes, turn the Reader on again. You can also turn the Reader off if you are not using the sensor or not. Be sure to read the test strip instructions for use.



◀ Previous

Next ▶

# Checking Glucose with a Test Strip

## How to do it:



Use your lancing device to obtain a blood drop and apply blood to the white area at the end of the test strip. Refer to your lancing device instructions for use if you need help using your lancing device. If sounds are turned on, the Reader beeps once to let you know you have applied enough blood.

You will see a butterfly on the screen while you wait for your result. Do not remove the test strip while the butterfly is on the screen. If sounds are turned on, the Reader beeps once when your result is ready.



If the butterfly does not appear, you may not have applied enough blood to the test strip. Apply a second drop of blood to the test strip within 5 seconds of the first drop. If the butterfly still does not appear or if more than 5 seconds have passed, discard the test strip. Turn off the Reader and repeat the steps in this section with a new test strip.



◀ Previous

Next ▶



# Checking Glucose

## How to do it:

Use your lancing device to obtain blood to the white area at the end of the test strip. Follow the lancing device instructions for your lancing device. If sounds like a needle, it is time to let you know you have applied enough blood.

You will see a butterfly on the screen when you apply blood to your result. Do not remove the test strip until you see the butterfly is on the screen. If so, place the test strip into the Reader beeps once when you hear the beep.

If the butterfly does not appear after 5 seconds, apply more blood to the test strip. Place the first drop of blood to the test strip within 5 seconds of the first drop. If the butterfly still does not appear after 5 seconds, than 5 seconds have passed, turn off the Reader and repeat the process with a new test strip.

## Caution




Test on your fingers in accordance with the Intended Use.

### Intended Use

The FreeStyle Libre Reader's built-in meter is for use outside the body only (in vitro diagnostic use) in the quantitative measurement of glucose in fresh whole blood for self testing by lay users from the fingers. It is not intended to be used for testing neonatal blood samples or for the diagnosis or screening of diabetes.

The FreeStyle Libre Reader's built-in meter is indicated for the home (lay) user in the management of patients with diabetes. It is intended to be used by a single person and should not be shared.

The FreeStyle Precision Neo Blood Glucose Test Strips are for use with the FreeStyle Libre Reader's built-in meter to quantitatively measure glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips.



# Checking Glucose

## How to do it:

Use your lancing device to obtain blood to the white area at the end of the test strip. Follow the steps in your lancing device instruction manual or refer to the instructions for your lancing device. If sounds are turned on, the Reader beeps once to let you know you have applied enough blood.

You will see a butterfly on the screen when the Reader is ready for your result. Do not remove the test strip until the butterfly disappears. If the butterfly is on the screen. If sounds are turned on, the Reader beeps once when your result is ready.


If the butterfly does not appear, you may not have applied enough blood to the test strip. Apply a second drop of blood to the test strip within 5 seconds of the first drop. If the butterfly still does not appear or if more than 5 seconds have passed, discard the test strip. Turn off the Reader and repeat the steps in this section with a new test strip.

## Note



- E-3 means the blood drop is too small, or incorrect test procedure, or there may be a problem with the test strip.
- E-4 means the blood glucose level may be too high to be read by the System or there may be a problem with the test strip.

See Troubleshooting section of the Reader Kit User's Manual for more information.



◀ Previous

Next ▶



# Checking Glucose with a Test Strip

## How to do it:



After reviewing your result, remove and discard the used test strip according to local regulations.



Blood glucose results are marked on the results screen and in the Reader's Logbook with the  symbol.



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Checking Glucose

## How to do it:

After reviewing your result, remove the used test strip according to regulations.

Blood glucose results are marked on the screen and in the Reader's Log symbol.

## Caution




Test on your fingers in accordance with the Intended Use.

### Intended Use

The FreeStyle Libre Reader's built-in meter is for use outside the body only (in vitro diagnostic use) in the quantitative measurement of glucose in fresh whole blood for self testing by lay users from the fingers. It is not intended to be used for testing neonatal blood samples or for the diagnosis or screening of diabetes.

The FreeStyle Libre Reader's built-in meter is indicated for the home (lay) user in the management of patients with diabetes. It is intended to be used by a single person and should not be shared.

The FreeStyle Precision Neo Blood Glucose Test Strips are for use with the FreeStyle Libre Reader's built-in meter to quantitatively measure glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips.



◀ Previous

Next ▶





## Important

### Checklist

#### How to do it

After reviewing the used regulations

Blood glucose screen area symbol.

FreeStyle  
Li  
FLASH GLUCOSE MONITOR

- Use only FreeStyle Precision Neo test strips. Other test strips may produce inaccurate results.
- Read all the instructions in this section. Failure to follow instructions may cause incorrect blood glucose results. Practice the testing procedures before using the Reader's built-in meter.
- Read the test strip instructions for use before performing your first blood glucose test as they contain important information. They also let you know how to store and handle the test strips and give you information about sample types.
- The Reader's built-in meter is not for use on people who are dehydrated, hypotensive, in shock, or for individuals in hyperglycemic-hyperosmolar state, with or without ketosis.
- The Reader's built-in meter is not for use on neonates, in critically-ill patients, or for diagnosis or screening of diabetes.
- Follow your health care professional's advice when testing blood glucose levels.
- Observe caution when using around children. Small parts may constitute a choking hazard.
- You should clean and disinfect the Reader once per week. The Reader should also be cleaned and disinfected prior to being handled by any person providing testing assistance to the user. Refer to the Maintenance and Disposal section of the Reader Kit User's manual for instructions.
- The Reader is for use by a single person. It must not be used on more than one person including other family members due to the risk of spreading infection. All parts of the Reader are considered biohazardous and can potentially transmit infectious diseases, even after performing the cleaning and disinfection procedure.
- Use the Reader's built-in meter within the test strip operating temperature range or you will see Error Message E-1.
- Use a test strip immediately after removing from its foil packet.
- Only use a test strip once.
- Do not put urine on the test strip.
- Do not use expired test strips as they may cause inaccurate results.
- Do not use a wet, bent, scratched, or damaged test strip.
- Do not use the test strip if the foil packet has a hole or is torn.
- Results from the built-in meter are shown only in your Reader's Logbook and not in other history options.
- Refer to your lancing device instructions for use for how to use your lancing device.
- After performing a blood glucose test, wash your hands with soap and water and thoroughly dry them.
- The built-in meter displays results from 20 - 500 mg/dL. Low or high blood glucose results can indicate a potentially serious medical condition.



# Checking Glucose


## How to do it:

After reviewing your result, remove the used test strip according to local regulations.

Blood glucose results are marked on the results screen and in the Reader's Logbook with the  symbol.

### Note

If you are using the App, you can enter your blood glucose result into the App's Logbook.



[Previous](#) | [Next](#)

# Adding Notes

Both the Reader and the App let you track food, insulin and exercise.

## How to do it with the Reader:

Press the Home Button to turn on the Reader. Check your glucose.

From the Glucose Reading screen, add notes by touching the symbol.

Select the checkbox next to the note you would like to add. Touch the down arrow to view additional note options.

After checking the box for food and insulin notes, the + symbol appears to the right of the note. Touch + to add more specific information to the note. Then touch OK.

Touch OK to save notes.



FreeStyle Libre



FreeStyle  
Libre

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



## Adding Notes

Both the Reader and the App have the ability to add notes.

### How to do it with the Reader:

Press the Home Button to turn on the Reader and view your glucose.

From the Glucose Reading screen, add notes by touching the pencil symbol.

Select the checkbox next to the note you would like to add. Touch the down arrow to view additional note options.

After checking the box for food and insulin notes, the + symbol appears to the right of the note. Touch + to add more specific information to the note. Then touch OK.

Touch OK to save notes.

### Note



- You can add a note at the time of your glucose reading or within 15 minutes after your reading was obtained.
- Food and rapid-acting insulin notes are shown on your glucose graphs and in your Logbook as symbols.



FreeStyle  
Libre

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



## Adding Notes


### How to do it with the App:

Tap the symbol on the My Glucose screen.

Select the checkbox next to the note you would like to add.

After you check the box, you can add more specific information to your note.

Tap **DONE** to save your note.



◀ Previous | Next ▶



## Adding Notes

### How to do it with the App:

Tap the symbol on the My

Select the checkbox next to the note you would like to

After you check the box, you can add more specific information to your note.


Tap **DONE** to save your note.

### Note



You cannot add a blood glucose result note to a glucose reading.

See Reviewing History section for how to enter blood glucose results in the App's Logbook.



FreeStyle  
*Libre*

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Adding Notes

## How to do it with the App:

Tap the symbol on the My

Select the checkbox next to the note you would like to

After you check the box, you can add more specific information to your note.

Tap **DONE** to save your note.


## Note

Notes are shown on your glucose graphs and in your Logbook as symbols. See App Symbols section for more information.



8:06 AM

115 ↑  
mg/dL



◀ Previous

Next ▶


# Reviewing History

Reviewing and understanding your glucose history can be an important tool for improving your glucose control. Both the Reader and App store about 90 days of information and have several ways to review your past glucose readings, notes, and other information.

## How to do it with the Reader:

Touch the **Review History** icon on the Home Screen.

Use the arrows to view the available options.



Previous | Next





## Reviewing History

Reviewing and understanding your glucose history is an important tool for improving your diabetes management.

Reader and App store about 9 days until the sensor needs to be replaced.

several ways to review your past glucose readings, notes, and other information.


### Important

Work with your health care professional to understand your glucose history.

#### How to do it with the Reader:

Touch the **Review History** icon on the Home Screen.

Use the arrows to view the available options.






# Reviewing History

## How to do it with the App:

From the Main Menu, tap Logbook to view the **Logbook** or tap on one of the other history options under **Reports**.



◀ Previous | Next ▶



# Reviewing History

## How to do it with the App:


From the Main Menu, tap Logs

or tap on one of the other history options under Reports.

### Important



Work with your health care professional to understand your glucose history.



**FreeStyle  
Libre**

FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Reviewing History

Summary of History Options - Reader

	What data is used?	What does this show?	Example
<b>Logbook</b>	Sensor scan results and blood glucose test results from each day.	Shows entries for each time you scanned your Sensor or performed a blood glucose test. If you entered notes with a glucose reading, the  symbol appears. Touch the entry to review detailed information including any notes you entered.	
Daily Graph			
Average Glucose			
Daily Patterns			
Time In Target			
Low Glucose Events			
Sensor Usage			




Previous | Next



# Reviewing History

Summary of History Options - Reader

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings from each day.	Shows a graph of your Sensor glucose readings by day. The blue bar indicates your Target Glucose Range. Symbols indicate any food or rapid-acting insulin notes you have entered.	
<b>Daily Graph</b>			
Average Glucose			
Daily Patterns			
Time In Target			
Low Glucose Events			
Sensor Usage			




◀ Previous | Next ▶



# Reviewing History

Summary of History Options - Reader

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows information about the average of your Sensor glucose readings. It includes the overall average and the average for four different 6-hour periods of the day.	
Daily Graph			
<b>Average Glucose</b>			
Daily Patterns			
Time In Target			
Low Glucose Events			
Sensor Usage			




 Previous | Next 



# Reviewing History

Summary of History Options - Reader

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows the pattern and variability of your Sensor glucose over a typical day.	
Daily Graph			
Average Glucose			
<b>Daily Patterns</b>			
Time In Target			
Low Glucose Events			
Sensor Usage			




◀ Previous | Next ▶



# Reviewing History

Summary of History Options - Reader

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows the percentage of time your Sensor glucose readings were above, below, or within your Target Glucose Range.	
Daily Graph			
Average Glucose			
Daily Patterns			
<b>Time In Target</b>			
Low Glucose Events			
Sensor Usage			




◀ Previous | Next ▶



# Reviewing History

Summary of History Options - Reader

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows the number of low glucose events measured by your Sensor. It includes the total number of events and the events in four different 6-hour periods of the day.	
Daily Graph			
Average Glucose			
Daily Patterns			
Time In Target			
<b>Low Glucose Events</b>			
Sensor Usage			




◀ Previous | Next ▶



# Reviewing History

Summary of History Options - Reader

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows how often you scan your Sensor. It includes an average of how many times you scanned your Sensor each day, and the percentage of possible Sensor data the Reader recorded from your scans.	
Daily Graph			
Average Glucose			
Daily Patterns			
Time In Target			
Low Glucose Events			
<b>Sensor Usage</b>			




 Previous | Next 



# Reviewing History

Summary of History Options - App

	What data is used?	What does this show?	Example
<b>Logbook</b>	Sensor scan results from each day. You can also use the Logbook to manually enter your blood glucose test results.	Shows entries for each time you scanned your Sensor as well as notes you added. The Logbook also lets you record a blood glucose test you performed. To do this, tap the  symbol and enter your result.	
Daily Graph			
Average Glucose			
Daily Patterns			
Time In Target			
Low Glucose Events			
Sensor Usage			



 Previous | Next 



# Reviewing History

Summary of History Options - App

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings from each day.	Shows a graph of your Sensor glucose readings by day. The green bar indicates your Target Glucose Range. Symbols indicate any notes you have entered.	
<b>Daily Graph</b>			
Average Glucose			
Daily Patterns			
Time In Target			
Low Glucose Events			
Sensor Usage			



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM


◀ Previous

Next ▶



# Reviewing History

Summary of History Options - App

	What data is used?	What does this show?	Example																				
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows information about the average of your Sensor glucose readings. It includes the overall average and the average for different periods of the day.	 <p>AVERAGE GLUCOSE ▾</p> <p>April 29–July 27, 2017</p> <table border="1"><thead><tr><th>Time</th><th>Glucose Level</th></tr></thead><tbody><tr><td>12am</td><td>149</td></tr><tr><td>3am</td><td>122</td></tr><tr><td>6am</td><td>112</td></tr><tr><td>9am</td><td>142</td></tr><tr><td>12pm</td><td>158</td></tr><tr><td>3pm</td><td>153</td></tr><tr><td>6pm</td><td>156</td></tr><tr><td>9pm</td><td>163</td></tr><tr><td>12am</td><td>163</td></tr></tbody></table> <p>Average: 144 mg/dL</p> <p>Data available for 90 of 90 days</p> <p>7 DAYS   14 DAYS   30 DAYS   90 DAYS</p>	Time	Glucose Level	12am	149	3am	122	6am	112	9am	142	12pm	158	3pm	153	6pm	156	9pm	163	12am	163
Time	Glucose Level																						
12am	149																						
3am	122																						
6am	112																						
9am	142																						
12pm	158																						
3pm	153																						
6pm	156																						
9pm	163																						
12am	163																						
Daily Graph																							
<b>Average Glucose</b>																							
Daily Patterns																							
Time In Target																							
Low Glucose Events																							
Sensor Usage																							




◀ Previous | Next ▶



# Reviewing History

Summary of History Options - App

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows the pattern and variability of your Sensor glucose over a typical day.	 <p>DAILY PATTERNS ▾</p> <p>June 28–July 27, 2017</p> <p>350 300 250 200 150 100 50</p> <p>12am 3am 6am 9am 12pm 3pm 6pm 9pm 12am</p> <p>8.3 mmol/L</p> <p>Data available for 30 of 30 days</p> <p>7 DAYS 14 DAYS 30 DAYS 90 DAYS</p>
Daily Graph			
Average Glucose			
<b>Daily Patterns</b>			
Time In Target			
Low Glucose Events			
Sensor Usage			



**FreeStyle**  
**Libre**  
FLASH GLUCOSE MONITORING SYSTEM

◀ Previous

Next ▶



# Reviewing History

Summary of History Options - App

	What data is used?	What does this show?	Example												
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows the percentage of time your Sensor glucose readings were above, below, or within your Target Glucose Range.	<p>TIME IN TARGET ▾</p> <p>July 14–27, 2017</p> <p>mg/dL</p> <table border="1"> <thead> <tr> <th>Glucose Range (mg/dL)</th> <th>Percentage (%)</th> </tr> </thead> <tbody> <tr> <td>&gt; 240</td> <td>3%</td> </tr> <tr> <td>141–240</td> <td>38%</td> </tr> <tr> <td>100–140</td> <td>32%</td> </tr> <tr> <td>70–99</td> <td>22%</td> </tr> <tr> <td>&lt; 70</td> <td>5%</td> </tr> </tbody> </table> <p>Target Range: 100 - 140 mg/dL Data available for 14 of 14 days</p> <p>7 DAYS 14 DAYS 30 DAYS 90 DAYS</p>	Glucose Range (mg/dL)	Percentage (%)	> 240	3%	141–240	38%	100–140	32%	70–99	22%	< 70	5%
Glucose Range (mg/dL)	Percentage (%)														
> 240	3%														
141–240	38%														
100–140	32%														
70–99	22%														
< 70	5%														
Daily Graph															
Average Glucose															
Daily Patterns															
<b>Time In Target</b>															
Low Glucose Events															
Sensor Usage															



◀ Previous

Next ▶



# Reviewing History

Summary of History Options - App

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows the number of low glucose events measured by your Sensor. It includes the total number of events and the events in different periods of the day.	
Daily Graph			
Average Glucose			
Daily Patterns			
Time In Target			
<b>Low Glucose Events</b>			
Sensor Usage			




◀ Previous | Next ▶



# Reviewing History

Summary of History Options - App

	What data is used?	What does this show?	Example
Logbook	Sensor glucose readings collected in the last 7, 14, 30, and 90 days.	Shows how often you scan your Sensor. It includes the total number of scans, an average of how many times you scanned your Sensor each day, and the percentage of possible Sensor data the App recorded from your scans.	
Daily Graph			
Average Glucose			
Daily Patterns			
Time In Target			
Low Glucose Events			
<b>Sensor Usage</b>			



◀ Previous | Next ▶



## More Information

[Important Safety Information](#) >

[Reader Symbols](#) >

[App Symbols](#) >

[System Specifications](#) >

several sequential readings over time. The System is intended for single patient use and requires a prescription.

## Contraindications



**MRI/CT/Diathermy:** The System must be removed prior to Magnetic Resonance Imaging (MRI), Computed Tomography (CT) scan, or high-frequency electrical heat (diathermy) treatment. The effect of MRI, CT scans, or diathermy on the performance of the System has not been evaluated. The exposure may damage the Sensor and may impact proper function of the device which could cause incorrect readings.

## WARNINGS:

- **Do not ignore symptoms that may be due to low or high blood glucose:** if you are experiencing symptoms that are not consistent with your glucose readings, consult your health care professional.
- **Checking Sensor glucose readings with a blood glucose meter:** Under the following conditions, Sensor glucose readings may not be accurate and you should conduct a fingerstick test using a blood glucose meter. You should not use Sensor glucose readings to make a diabetes treatment decision:
  - If you suspect that your reading may be inaccurate for any reason
  - When you are experiencing symptoms that may be due to low or high blood glucose
  - When you are experiencing symptoms that do not match the Sensor glucose readings
  - During times of rapidly changing glucose (more than 2 mg/dL per minute), when interstitial fluid glucose levels as measured by the Sensor may not accurately reflect blood glucose levels
  - When the Sensor glucose reading does not include a Current Glucose number or Glucose Trend Arrow
  - In order to confirm hypoglycemia or impending hypoglycemia as reported by the Sensor

- When you see the  symbol, you must check your blood glucose with a blood glucose meter before making any treatment decisions. Sensor readings may not accurately reflect blood glucose levels.
- If you are using the FreeStyle LibreLink app, you must also have access to a blood glucose monitoring system as the App does not provide one.
- **Hypoglycemic unawareness:** The System has not been evaluated for use in patients with hypoglycemic unawareness and will not automatically alert you of a hypoglycemic event without you scanning your Sensor.
- **No alarms without a Sensor scan:** The System does not have alarms that will automatically notify you when you are having a severe low (hypoglycemic) or high (hyperglycemic) glucose event unless you scan your Sensor. For example, the System does not have an alarm that can alert or wake you when you are sleeping in the case of low or high glucose.
- **Choking hazard:** The System contains small parts that may be dangerous if swallowed.

## Cautions and Limitations

Below are important cautions and limitations to keep in mind so you can use the System safely. They are grouped into categories for easy reference.



### What to know about Alarms/Alerts:

- There are NO alarms or alerts unless you scan the Sensor.



### What to know before using the System:

- Review all product information before use.
- Take standard precautions for transmission of blood borne pathogens to avoid contamination.



### Who should not use the System:

- **Do not use the System in people less than 18 years of age.** The System is not approved for use in people under 18 years of age and Sensor readings in this population may be inaccurate. In general, continuous glucose monitoring systems are recognized to be less accurate in children than in adults.
- **Do not use the System in critically ill patients.** The System is not approved for use in these patients. It is not known how different conditions or medications common to the critically ill population may affect performance of the System. Sensor glucose readings may be inaccurate in critically ill patients.
- **Do not use the System in pregnant women or persons on dialysis.** The System is not approved for use in pregnant women or persons on dialysis and has not been evaluated in these populations.
- Performance of the System when used with other implanted medical devices, such as pacemakers, has not been evaluated.



### What should you know about wearing a Sensor:

- After the start-up period, the Sensor can be worn for up to the wear duration specified in your Sensor Kit's product insert.
- Some individuals may be sensitive to the adhesive that keeps the Sensor attached to the skin. If you notice significant skin irritation around or under your Sensor, remove the Sensor and stop using the System. Contact your health care professional before continuing to use the System.
- Intense exercise may cause your Sensor to loosen due to sweat or movement of the Sensor. Remove and replace your Sensor if it starts to loosen and follow the instructions to select an appropriate application site.
- Do not reuse Sensors. The Sensor and Sensor Applicator are designed for single use. Reuse may result in no glucose readings and infection. Not suitable for re-sterilization. Further exposure to irradiation may cause inaccurate results.
- If a Sensor breaks inside your body, call your health care professional.



### **How to Store the Sensor Kit:**

- Store the Sensor Kit between 39°F and 77°F. Storage outside of this range may cause inaccurate Sensor glucose readings. While you don't need to keep your Sensor Kit in a refrigerator, you can as long as the refrigerator is between 39°F and 77°F. Do not freeze.
- Store the Sensor Kit between 10-90% non-condensing humidity.




### **When not to use the System:**

- Do NOT use if the Sensor Kit package, Sensor Pack or Sensor Applicator appear to be damaged or already opened due to risk of no results and/or infection.
- Do NOT use if Sensor Kit contents are past expiration date.
- Do NOT use if the Reader appears to be damaged due to risk of electric shock and/or no results.



### **What to know before you Apply the Sensor:**

- The Sensor Pack and Sensor Applicator are packaged as a set (separately from the Reader) and have the same Sensor code. Check that the Sensor codes match before using your Sensor Pack and Sensor Applicator. Do not use Sensor Packs and Sensor Applicators with different Sensor codes together as this will result in incorrect glucose readings.
- Clean the application site and ensure that it is dry prior to Sensor insertion. This helps the Sensor stay attached to your body.
- Clean hands prior to Sensor handling/insertion to help prevent infection.
- Change the application site for the next Sensor application to prevent discomfort or skin irritation.
- Sensor placement is not approved for sites other than the back of the arm. If placed in other areas, the Sensor may not function properly.



- Select an appropriate Sensor site to help the Sensor stay attached to the body and prevent discomfort or skin irritation. Avoid areas with scars, moles, stretch marks, or lumps. Select an area of skin that generally stays flat during normal daily activities (no bending or folding). Choose a site that is at least 1 inch away from an insulin injection site.



#### **When is Sensor Glucose different from Blood Glucose:**

- Physiological differences between the interstitial fluid and capillary blood may result in differences in glucose readings between the System and results from a fingerstick test using a blood glucose meter. Differences in glucose readings between interstitial fluid and capillary blood may be observed during times of rapid change in blood glucose, such as after eating, dosing insulin, or exercising.



#### **What to know about interfering substances such as Vitamin C and Aspirin:**

- Taking ascorbic acid (vitamin C) while wearing the Sensor may falsely raise Sensor glucose readings. Taking salicylic acid (used in some pain relievers such as aspirin and some skin care products) may slightly lower Sensor glucose readings. The level of inaccuracy depends on the amount of the interfering substance active in the body.
- Test results did not indicate interference for methyldopa (used in some drugs to treat high blood pressure) or tolbutamide (infrequently used in some drugs to treat diabetes in the US) at maximum circulating levels. However, concentrations of potential interferents in interstitial fluid are unknown compared to circulating blood.



#### **What to know about X-Rays:**

- The Sensor should be removed prior to exposing it to an X-ray machine. The effect of X-rays on the performance of the System has not been evaluated. The exposure may damage the Sensor and may impact proper function of the device to detect trends and track patterns in glucose values during the wear period.



#### **When to remove the Sensor:**

- If the Sensor is becoming loose or if the Sensor tip is coming out of your skin, you may get no readings or unreliable readings, which may not match how you feel. Check to make sure your Sensor has not come loose. If it has come loose, remove it and apply a new one.
- If you believe your glucose readings are not correct or are inconsistent with how you feel, perform a blood glucose test on your finger to confirm your glucose. If the problem continues, remove the current Sensor and apply a new one.



#### **What to do if you are dehydrated:**

- Severe dehydration and excessive water loss may cause inaccurate Sensor glucose readings. If you believe you are suffering from dehydration, consult your health care professional immediately.



#### **What to know about the Reader's Built-in Meter:**

- The FreeStyle Libre Reader has a built-in blood glucose meter that is designed to be used only with FreeStyle Precision Neo blood glucose test strips and MediSense Glucose and Ketone Control Solution. Using other test strips with the Reader's built-in meter will produce an error or cause the Reader's built-in meter to not turn on or start a test. The Reader's built-in meter does not have ketone testing functionality.
- The Reader's built-in meter is not for use on people who are dehydrated, hypotensive, in shock, or for individuals in hyperglycemic-hyperosmolar state, with or without ketosis.
- The Reader's built-in meter is not for use on neonates, in critically-ill patients, or for diagnosis or screening of diabetes.
- See Using the Reader's Built-in meter section of the Reader Kit User's Manual for additional important information on the use of the Reader's built-in meter.



#### **Where to charge your Reader:**

- Be sure to select a location for charging that allows the power adapter to be easily unplugged. Do NOT block access to the charger due to the potential risk of electrical shock.



#### **What to know about FreeStyle LibreLink:**

- FreeStyle LibreLink installed on a smartphone is intended for use by a single person. It must not be used by more than one person due to the risk of misinterpreting glucose information.
- FreeStyle LibreLink and FreeStyle Libre Readers do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.



◀ Previous

Next ▶

# Important information about the FreeStyle Libre System

## Important Safety Information

### Indications for Use

The FreeStyle Libre Flash Glucose Monitoring System is a continuous glucose monitoring (CGM) device indicated for the management of diabetes in persons age 18 and older. It is designed to replace blood glucose testing for diabetes treatment decisions.

The System detects trends and tracks patterns aiding in the detection of episodes of hyperglycemia and hypoglycemia, facilitating both acute and long-term therapy adjustments. Interpretation of the System readings should be based on the glucose trends and several sequential readings over time. The System is intended for single patient use and requires a prescription.

### Contraindications

**MRI/CT/Diathermy:** The System must be removed prior to Magnetic Resonance Imaging (MRI), Computed Tomography (CT) scan, or high-frequency electrical heat (diathermy) treatment. The effect of MRI, CT scans, or diathermy on the performance of the System has not been evaluated. The exposure may damage the Sensor and may impact proper function of the device which could cause incorrect readings.

### WARNINGS:

- **Do not ignore symptoms that may be due to low or high blood glucose:** if you are experiencing symptoms that are not consistent with your glucose readings, consult your health care professional.
- **Checking Sensor glucose readings with a blood glucose meter:** Under the following conditions, Sensor glucose readings may not be accurate and you should conduct a fingerstick test using a blood glucose meter. You should not use Sensor glucose readings to make a diabetes treatment decision:
  - If you suspect that your reading may be inaccurate for any reason
  - When you are experiencing symptoms that may be due to low or high blood glucose
  - When you are experiencing symptoms that do not match the Sensor glucose readings
  - During times of rapidly changing glucose (more than 2 mg/dL per minute), when interstitial fluid glucose levels as measured by the Sensor may not accurately reflect blood glucose levels
  - When the Sensor glucose reading does not include a Current Glucose number or Glucose Trend Arrow
  - In order to confirm hypoglycemia or impending hypoglycemia as reported by the Sensor

- When you see the  symbol, you must check your blood glucose with a blood glucose meter before making any treatment decisions. Sensor readings may not accurately reflect blood glucose levels.
- If you are using the FreeStyle LibreLink app, you must also have access to a blood glucose monitoring system as the App does not provide one.
- **Hypoglycemic unawareness:** The System has not been evaluated for use in patients with hypoglycemic unawareness and will not automatically alert you of a hypoglycemic event without you scanning your Sensor.
- **No alarms without a Sensor scan:** The System does not have alarms that will automatically notify you when you are having a severe low (hypoglycemic) or high (hyperglycemic) glucose event unless you scan your Sensor. For example, the System does not have an alarm that can alert or wake you when you are sleeping in the case of low or high glucose.
- **Choking hazard:** The System contains small parts that may be dangerous if swallowed.

## Cautions and Limitations

Below are important cautions and limitations to keep in mind so you can use the System safely. They are grouped into categories for easy reference.



### What to know about Alarms/Alerts:

- There are NO alarms or alerts unless you scan the Sensor.



### What to know before using the System:

- Review all product information before use.
- Take standard precautions for transmission of blood borne pathogens to avoid contamination.



### Who should not use the System:

- **Do not use the System in people less than 18 years of age.** The System is not approved for use in people under 18 years of age and Sensor readings in this population may be inaccurate. In general, continuous glucose monitoring systems are recognized to be less accurate in children than in adults.
- **Do not use the System in critically ill patients.** The System is not approved for use in these patients. It is not known how different conditions or medications common to the critically ill population may affect performance of the System. Sensor glucose readings may be inaccurate in critically ill patients.
- **Do not use the System in pregnant women or persons on dialysis.** The System is not approved for use in pregnant women or persons on dialysis and has not been evaluated in these populations.
- Performance of the System when used with other implanted medical devices, such as pacemakers, has not been evaluated.



### What should you know about wearing a Sensor:

- After the start-up period, the Sensor can be worn for up to the wear duration specified in your Sensor Kit's product insert.
- Some individuals may be sensitive to the adhesive that keeps the Sensor attached to the skin. If you notice significant skin irritation around or under your Sensor, remove the Sensor and stop using the System. Contact your health care professional before continuing to use the System.
- Intense exercise may cause your Sensor to loosen due to sweat or movement of the Sensor. Remove and replace your Sensor if it starts to loosen and follow the instructions to select an appropriate application site.
- Do not reuse Sensors. The Sensor and Sensor Applicator are designed for single use. Reuse may result in no glucose readings and infection. Not suitable for re-sterilization. Further exposure to irradiation may cause inaccurate results.
- If a Sensor breaks inside your body, call your health care professional.



### **How to Store the Sensor Kit:**

- Store the Sensor Kit between 39°F and 77°F. Storage outside of this range may cause inaccurate Sensor glucose readings. While you don't need to keep your Sensor Kit in a refrigerator, you can as long as the refrigerator is between 39°F and 77°F. Do not freeze.
- Store the Sensor Kit between 10-90% non-condensing humidity.




### **When not to use the System:**

- Do NOT use if the Sensor Kit package, Sensor Pack or Sensor Applicator appear to be damaged or already opened due to risk of no results and/or infection.
- Do NOT use if Sensor Kit contents are past expiration date.
- Do NOT use if the Reader appears to be damaged due to risk of electric shock and/or no results.



### **What to know before you Apply the Sensor:**

- The Sensor Pack and Sensor Applicator are packaged as a set (separately from the Reader) and have the same Sensor code. Check that the Sensor codes match before using your Sensor Pack and Sensor Applicator. Do not use Sensor Packs and Sensor Applicators with different Sensor codes together as this will result in incorrect glucose readings.
- Clean the application site and ensure that it is dry prior to Sensor insertion. This helps the Sensor stay attached to your body.
- Clean hands prior to Sensor handling/insertion to help prevent infection.
- Change the application site for the next Sensor application to prevent discomfort or skin irritation.
- Sensor placement is not approved for sites other than the back of the arm. If placed in other areas, the Sensor may not function properly.



- Select an appropriate Sensor site to help the Sensor stay attached to the body and prevent discomfort or skin irritation. Avoid areas with scars, moles, stretch marks, or lumps. Select an area of skin that generally stays flat during normal daily activities (no bending or folding). Choose a site that is at least 1 inch away from an insulin injection site.



#### **When is Sensor Glucose different from Blood Glucose:**

- Physiological differences between the interstitial fluid and capillary blood may result in differences in glucose readings between the System and results from a fingerstick test using a blood glucose meter. Differences in glucose readings between interstitial fluid and capillary blood may be observed during times of rapid change in blood glucose, such as after eating, dosing insulin, or exercising.



#### **What to know about interfering substances such as Vitamin C and Aspirin:**

- Taking ascorbic acid (vitamin C) while wearing the Sensor may falsely raise Sensor glucose readings. Taking salicylic acid (used in some pain relievers such as aspirin and some skin care products) may slightly lower Sensor glucose readings. The level of inaccuracy depends on the amount of the interfering substance active in the body.
- Test results did not indicate interference for methyldopa (used in some drugs to treat high blood pressure) or tolbutamide (infrequently used in some drugs to treat diabetes in the US) at maximum circulating levels. However, concentrations of potential interferents in interstitial fluid are unknown compared to circulating blood.



#### **What to know about X-Rays:**

- The Sensor should be removed prior to exposing it to an X-ray machine. The effect of X-rays on the performance of the System has not been evaluated. The exposure may damage the Sensor and may impact proper function of the device to detect trends and track patterns in glucose values during the wear period.



#### **When to remove the Sensor:**

- If the Sensor is becoming loose or if the Sensor tip is coming out of your skin, you may get no readings or unreliable readings, which may not match how you feel. Check to make sure your Sensor has not come loose. If it has come loose, remove it and apply a new one.
- If you believe your glucose readings are not correct or are inconsistent with how you feel, perform a blood glucose test on your finger to confirm your glucose. If the problem continues, remove the current Sensor and apply a new one.



#### **What to do if you are dehydrated:**

- Severe dehydration and excessive water loss may cause inaccurate Sensor glucose readings. If you believe you are suffering from dehydration, consult your health care professional immediately.



#### **What to know about the Reader's Built-in Meter:**

- The FreeStyle Libre Reader has a built-in blood glucose meter that is designed to be used only with FreeStyle Precision Neo blood glucose test strips and MediSense Glucose and Ketone Control Solution. Using other test strips with the Reader's built-in meter will produce an error or cause the Reader's built-in meter to not turn on or start a test. The Reader's built-in meter does not have ketone testing functionality.
- The Reader's built-in meter is not for use on people who are dehydrated, hypotensive, in shock, or for individuals in hyperglycemic-hyperosmolar state, with or without ketosis.
- The Reader's built-in meter is not for use on neonates, in critically-ill patients, or for diagnosis or screening of diabetes.
- See Using the Reader's Built-in meter section of the Reader Kit User's Manual for additional important information on the use of the Reader's built-in meter.



#### **Where to charge your Reader:**

- Be sure to select a location for charging that allows the power adapter to be easily unplugged. Do NOT block access to the charger due to the potential risk of electrical shock.



#### **What to know about FreeStyle LibreLink:**

- FreeStyle LibreLink installed on a smartphone is intended for use by a single person. It must not be used by more than one person due to the risk of misinterpreting glucose information.
- FreeStyle LibreLink and FreeStyle Libre Readers do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.



◀ Previous

Next ▶



# Reader Symbols



Sensor may be inaccurate. Check blood glucose with a test strip before making any treatment decisions



Sensor

↑ ↗ → ↘ ↓ Direction glucose is going



Caution

◀ ▼ ▲ ▶ View previous/next screen



Notes



Add more information to notes



Food note



FreeStyle  
*Libre*

FLASH GLUCOSE MONITORING SYSTEM



Rapid-acting insulin note



Time changed on Reader



Reminders



Blood glucose test



Settings



Control solution test result



Low battery



Battery charging



Sensor too cold



Sensor too hot

◀ Previous

Next ▶

# App Symbols



App icon



Sensor may be inaccurate. Check blood glucose with a test strip before making any treatment decisions



Scan button



Direction glucose is going



Caution



Add/edit notes



Manually entered blood glucose result note



Add blood glucose result note



Food note



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM



Insulin (Rapid or Long-acting) note



Exercise note



Time change



Main menu



Multiple/Custom notes



Share report



Additional information



Calendar



Sensor too cold



Sensor too hot

Previous

Next

# System Specifications

[Sensor Specifications](#)[Reader Specifications](#)

Sensor glucose assay method	Amperometric electrochemical sensor
Sensor glucose reading range	40 to 500 mg/dL
Sensor size	5 mm height and 35 mm diameter
Sensor weight	5 grams
Sensor power source	One silver oxide battery
Sensor data	Up to the wear duration specified in your Sensor Kit's product insert
Sensor memory	8 hours (glucose readings stored every 15 minutes)
Operating temperature	50 °F to 113 °F
Sensor Applicator and Sensor Pack storage temperature	39 °F to 77 °F
Operating and storage relative humidity	10-90%, non-condensing
Sensor water resistance	IP27: Can withstand immersion into 3 ft (one meter) of water for up to 30 minutes. Protected against insertion of objects > 12 mm diameter
Operating and storage altitude	-1,250 ft (-381 meters) to 10,000 ft (3,048 meters)

Please refer to the Reader Kit User's Manual for updates to the information.



**FreeStyle  
Libre**  
FLASH GLUCOSE MONITORING SYSTEM

[◀ Previous](#)[Next ▶](#)



# System Specifications

Sensor Specifications

Reader Specifications

<b>Blood glucose assay range</b>	20 to 500 mg/dL	<b>Reader display timeout</b>	60 seconds (120 seconds when test strip is inserted)
<b>Reader size</b>	95 mm x 60 mm x 16 mm	<b>Radio Frequency</b>	Near Field Communication* (13.56 MHz RFID); ASK Modulation; 124 dBuV/m; 1.5 inch communication range
<b>Reader weight</b>	65 grams	<b>Data port</b>	Micro USB
<b>Reader power source</b>	One lithium-ion rechargeable battery	<b>Minimum Computer Requirements</b>	System must only be used with EN60950-1 rated computers
<b>Reader battery life</b>	7 days of typical use	<b>Mean service life</b>	3 years of typical use
<b>Reader memory</b>	90 days of typical use	<b>Reader cleaning and disinfection</b>	The Reader has a mean use life of 3 years, which is 156 cleaning and disinfection cycles (1 cycle per week for 3 years).
<b>Reader operating temperature</b>	50 °F to 113 °F	<b>Power Adapter</b>	Abbott Diabetes Care PRT25611 Operating temperature: 50 °F to 104 °F
<b>Reader storage temperature</b>	-4 °F to 140 °F	<b>USB Cable</b>	Abbott Diabetes Care PRT21373 Length: 37 inches (94 cm)
<b>Operating and storage relative humidity</b>	10-90%, non-condensing		
<b>Reader moisture protection</b>	Keep dry		
<b>Operating and storage altitude</b>	-1,250 ft (-381 meters) to 10,000 ft (3,048 meters)		

\*Security measures: The communication between Reader and Sensor is a short range near field communication method making it difficult to interfere with or intercept data that is being transferred. The Sensor and Reader are protected by proprietary data format, memory mapping, and cyclic redundancy check (CRC) generation and verification of data.

Quality of Service (QoS): QoS for the FreeStyle Libre Reader and Sensor wireless communications using the near field communications is assured within the effective range of 4 cm between the Sensor and Reader that is specified to occur within 15 seconds.



◀ Previous

| Home

▶