FreeStyle Libre Flash Glucose Monitoring System

NOTE: Please consult your health care team on how to use the information in this section.

Performance Characteristics

Performance of the Sensor was evaluated in a controlled clinical study. The study was conducted in 5 centres and a total of 146 subjects with diabetes were included in the effectiveness analysis. Each subject wore up to two Sensors for up to 14 days on the back of the upper arm. During the study, subjects had their venous blood glucose analysed over three separate visits to the clinical centre using the Yellow Springs Instrument Life Sciences 2300 STAT Plus™. Three lots of Sensors were evaluated in the study.


Fig 1. Comparison of the Sensors vs. YSI reference.

Slope	0,97
Intercept	-1,3 mg/dL (-0,1 mmol/L)
Correlation	0,98
N	18926
Range	37 - 479 mg/dL (2,0 – 26,6 mmol/L)
Overall mean bias	-5,6 mg/dL (-0,3 mmol/L)
Mean Absolute Relative Difference (MARD)	9,2%

Table 2. Sensor accuracy for all results vs. YSI reference

Sensor accuracy results for glucose concentrations <80 mg/dL (4,4 mmol/L)	Within ±15 mg/dL (within ±0,83 mmol/L)	Within ±20 mg/dL (within ±1,11 mmol/L)	Within ±30 mg/dL (within ±1,67 mmol/L)
	4199 / 4595 (91,4%)	4482 / 4595 (97,5%)	4583 / 4595 (99,7%)
Sensor accuracy results for glucose concentrations ≥80 mg/dL (4,4 mmol/L)	Within ±15%	Within ±20%	Within ±30%
	12143 / 14331 (84,7%)	13153 / 14331 (91,8%)	14012 / 14331 (97,8%)
Sensor accuracy for all results	Within ±20 mg/dL (±1,11 mmol/L) and within ±20% of reference		
	17635 / 18926 (93,2%)		

Table 3. Sensor performance relative to YSI reference at different glucose levels

Glucose	Mean Absolute Relative Difference
≤50 mg/dL (2,8 mmol/L)	9,1 mg/dL (0,5 mmol/L)*
51-80 mg/dL (2,8-4,4 mmol/L)	7,0 mg/dL (0,4 mmol/L)*
81-180 mg/dL (4,5-10,0 mmol/L)	10,1%
181-300 mg/dL (10,0-16,7 mmol/L)	7,5%
301-400 mg/dL (16,7-22,2 mmol/L)	7,1%
>400 mg/dL (22,2 mmol/L)	10,2%

* For glucose ≤80 mg/dL (4,4 mmol/L), the differences in mg/dL (mmol/L) are presented instead of relative differences (%).

Table 4. Sensor accuracy over wear duration vs. YSI reference

	Beginning	Early Middle	Late Middle	End
Within ±20 mg/dL (±1,11 mmol/L) and within ±20% of reference	91,2%	95,1%	94,2%	93,7%
Mean Absolute Relative Difference (%)	10,0	8,5	8,8	9,1

Skin Interaction

Based on the examination of 146 study participants, the following incidence of skin issues were observed. Four occurrences of erythema were reported to be moderate in intensity. All other skin issues were reported to be mild in intensity.

Bleeding – 0,7% of the time
Bruising – 0,7% of the time
Erythema – 2,7% of the time
Pain – 0,7% of the time
Scabbing – 2,7% of the time

Expected Clinical Benefits

Complications as a result of diabetes mellitus (including, but not limited to: diabetic retinopathy, diabetic nephropathy) are well documented.¹ Self-monitoring of blood glucose (SMBG) by patients has revolutionized management of diabetes.² Using glucose monitoring devices patients with diabetes can work to achieve and maintain specific glycaemic goals. Given the results of the Diabetes Control and Complications Trial (DCCT)³ and other studies, there is broad consensus on the health benefits of normal or near-normal blood glucose levels and on the importance, especially in insulin-treated patients, of glucose monitoring devices in treatment efforts designed to achieve these glycaemic goals. Based principally on the DCCT results, experts recommend that most individuals with diabetes should attempt to achieve and maintain blood glucose levels as close to normal as is safely possible. Most patients with diabetes, especially insulin treated patients, can achieve this goal only by using glucose monitoring devices.

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