FreeStyle  
*LibreLink*

English

Deutsch

Français

Italiano

Nederlands

Svenska

Norsk

Suomi

X

**App Symbols****Important Information**

Indications for Use

**FreeStyle LibreLink Overview**

Home Screen

Sensor Kit

**App Setup****Applying Your Sensor****Starting Your Sensor****Checking Your Glucose****Understanding Your Glucose Readings****Alarms with a FreeStyle Libre 2 Sensor****Setting Alarms****Using Alarms****Adding Notes****Reviewing Your History**

Logbook

Other History Options

**Removing Your Sensor****Replacing Your Sensor****Setting Reminders****App Settings and Other Options****Using the Connected Apps Option****Living with your Sensor**

Activities

Maintenance

Disposal

**Troubleshooting**

Problems at the Sensor Application Site

Problems Starting Your Sensor or Receiving  
Sensor Readings

Problems Receiving Glucose Alarms

**Customer Service****Labelling Symbols and Definitions****Electromagnetic Compatibility****Performance Characteristics**

## User's Manual

# App Symbols

Direction your glucose is going. See [Understanding Your Glucose Readings](#) for more information.

Caution

Add/edit notes

Food note

Insulin (Rapid or Long-acting) note

Exercise note

Time change

Alarms you have turned on are unavailable

Sensor too cold

Sensor too hot

App icon

Multiple/Custom notes

Share report

---

Additional information

---

Main menu

---

Overflow menu

---

Calendar

---

## Important Information

### Indications for Use

FreeStyle LibreLink app ("App") when used with a FreeStyle Libre or FreeStyle Libre 2 Flash Glucose Monitoring System Sensor ("Sensor") is indicated for measuring interstitial fluid glucose levels in people (age 4 and older) with diabetes mellitus, including pregnant women. The App and Sensor are designed to replace blood glucose testing in the self-management of diabetes, including dosing of insulin.

The indication for children (age 4 - 12) is limited to those who are supervised by a caregiver who is at least 18 years of age. The caregiver is responsible for managing or assisting the child to manage the App and Sensor and also for interpreting or assisting the child to interpret Sensor glucose readings.

**Note:** Not all Sensors are available in all countries.

**WARNING:** If you are using FreeStyle LibreLink, you must also have access to a blood glucose monitoring system as the App does not provide one.

### CAUTION:

- FreeStyle LibreLink installed on a smartphone is intended for use by a single person. It must not be used by more than one person due to the risk of misinterpreting glucose information.
- If you are using a FreeStyle Libre Sensor with the App or you started your

FreeStyle Libre 2 Sensor with the Reader, you will not receive alarms from the App.

### No Alarms from the App

 You are using a FreeStyle Libre Sensor.

 You started a FreeStyle Libre 2 Sensor with the Reader before using it with the App.

### Alarms from the App

 You started a FreeStyle Libre 2 Sensor with the App.

- You will only get alarms from the App if you use the App to start a FreeStyle Libre 2 Sensor. For you to receive alarms, make sure to:
  - Turn alarms **ON** and ensure that your smartphone is within 6 metres (20 ft) of you at all times. The transmission range is 6 metres (20 ft) unobstructed. If you are out of range, you may not receive glucose alarms.
  - Do not force close the App.
  - Check to make sure that you have the correct settings and permissions enabled on your phone to receive alarms.
    - Enable Bluetooth and Notifications for the App. Be sure not to turn on any features or modify phone settings that may interrupt the presentation of notifications.
    - Turn on Lock Screen notifications, Notification sounds, Channel notifications, Sound and Pop-up notifications, and general phone sounds or vibration.
    - Turn off Do Not Disturb mode or manage the settings for your alarms to override Do Not Disturb. Override Do Not Disturb functionality depends on the phone model and Android operating system version you are using.
    - Turn on Location Permissions. Android 6.0 and above require apps to have Location Permissions enabled to connect with Bluetooth devices. If

you wish to receive alarms, Location Permissions must be turned on.

- Turn off Battery Optimisation for the App to allow the App to run in the background and ensure you receive alarms, even when your phone has low battery.
- You may need to add the App to the list of apps that will not be put to sleep.
- Remember that alarm settings will follow your smartphone's sound and vibration settings so these should be at a level you can hear to prevent missed alarms.
- You should disconnect headphones when you are not using them as you may not receive audio with alarms.
- If you are using peripheral devices connected to your phone, such as wireless headphones or a smartwatch, you may receive alarms on only one device or peripheral, not all.
- Keep your smartphone well charged and turned on.

## Additional Safety Information

FreeStyle LibreLink and FreeStyle Libre or FreeStyle Libre 2 Readers ("Readers") do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.

## Security Information

- You are responsible for properly securing and managing your smartphone. If you suspect an adverse cybersecurity event related to FreeStyle LibreLink, contact Customer Service.
- FreeStyle LibreLink is not intended for use on a smartphone that has been altered or customised to remove, replace or circumvent the manufacturer's approved configuration or use restriction, or that otherwise violates the manufacturer's warranty.

**The following Contraindication, Warnings, and other safety information apply to the Sensor, when used with FreeStyle LibreLink.**

**CONTRAINDICATION:** The Sensor must be removed prior to Magnetic Resonance Imaging (MRI).

**WARNING:**

- The Sensor contains small parts that may be dangerous if swallowed.
- Do not ignore symptoms that may be due to low or high blood glucose. If you have symptoms that do not match the Sensor glucose reading, or suspect your reading may be inaccurate, check the reading by conducting a fingerstick test using a blood glucose meter. If you are experiencing symptoms that are not consistent with your glucose readings, consult your health care professional.
- The FreeStyle Libre 2 Sensor can be used with the FreeStyle Libre Reader but the FreeStyle Libre Reader will NOT issue alarms.

**CAUTION:**

- On rare occasions, you may get inaccurate Sensor glucose readings. If you believe your readings are not correct or are inconsistent with how you feel, perform a blood glucose test on your finger to confirm your glucose and check to make sure your Sensor has not come loose. If the problem continues or if your Sensor is coming loose, remove the current Sensor and apply a new one.
- Intense exercise may cause your Sensor to loosen due to sweat or movement of the Sensor. If your Sensor comes loose, you may get no readings or unreliable readings, which may not match how you feel. Follow the instructions to select an appropriate application site.
- The Sensor uses all available glucose data to give you readings so you should scan your Sensor at least once every 8 hours for the most accurate performance. Scanning less frequently may result in decreased performance. If you are using both the App and Reader with the same Sensor, be sure to scan frequently with both devices.
- Some individuals may be sensitive to the adhesive that keeps the Sensor attached to the skin. If you notice significant skin irritation around or under your Sensor, remove the Sensor and stop using the Sensor. Contact your health care professional before continuing to use the Sensor.
- Performance of the Sensor when used with other implanted medical devices,

such as pacemakers, has not been evaluated.

- Do not reuse Sensors. The Sensor and Sensor Applicator are designed for single use. Reuse may result in no glucose readings and infection. Not suitable for re-sterilisation. Further exposure to irradiation may cause inaccurate results.
- The Sensor Pack and Sensor Applicator are packaged as a set and have the same Sensor code. Check that the Sensor codes match before using your Sensor Pack and Sensor Applicator. Sensor Packs and Sensor Applicators with the same Sensor code should be used together or your Sensor glucose readings may be incorrect.

### **Additional Safety Information**

- Physiological differences between the interstitial fluid and capillary blood may result in differences in glucose readings. Differences in Sensor glucose readings between interstitial fluid and capillary blood may be observed during times of rapid change in blood glucose, such as after eating, dosing insulin, or exercising.
- Store the Sensor Kit between 4°C and 25°C. While you don't need to keep your Sensor Kit in a refrigerator, you can as long as the refrigerator is between 4°C and 25°C.
- If you have a medical appointment that includes strong magnetic or electromagnetic radiation, for example an X-ray, MRI (Magnetic Resonance Imaging), or CT (Computed Tomography) scan, remove the Sensor you are wearing and apply a new one after the appointment. The effect of these types of procedures on the performance of the Sensor has not been evaluated.
- The Sensor has not been evaluated for use in persons on dialysis or people less than 4 years of age.
- The Sensor Pack is sterile unless opened or damaged.
- Your Sensor has been tested to withstand immersion into one metre (3 ft) of water for up to 30 minutes. It is also protected against insertion of objects > 12mm diameter. (IP27)
- Do not freeze the Sensor. Do not use if expiry date has passed.

## **FreeStyle LibreLink Overview**

**IMPORTANT:** Read all of the information in this User's Manual before using FreeStyle LibreLink with a Sensor. Refer to your smartphone's instructions for use for how to use your smartphone. If you are using a Reader, refer to the User's Manual in the Reader Kit.


FreeStyle LibreLink is available for download from the Google Play Store. When you're ready to start using FreeStyle LibreLink, you'll prepare and apply a Sensor to the back of your upper arm. You can then use the App to get glucose readings from the Sensor and store your glucose history and notes. The App can be used with either the FreeStyle Libre Sensor or FreeStyle Libre 2 Sensor. Each Sensor comes in a [Sensor Kit](#) and can be worn on your body for up to 14 days.

### Note:

- Not all Sensors are available in all countries.
- Go to [www.FreeStyleLibre.com](http://www.FreeStyleLibre.com) for smartphone requirements and compatibility. Please keep in mind that the ease of scanning a Sensor may vary between devices.

## Home Screen

The Home Screen gives you access to information about your glucose and the App. To return to the Home Screen from another screen, go to the Main Menu and tap **Home**.

**Main Menu** - Tap to access the Home Screen, Logbook, other history options, and Connected Apps.

**Glucose Graph** - Graph of your stored Sensor glucose readings.

**Overflow Menu** - Tap to change App settings and view App information.

**Scan Prompt** - Indicates if the App is ready to scan a Sensor.

**Glucose Information** - Your Time In Target, information about your last scan, and average glucose for the last 24 hours.

**Target Glucose Range** - The graph shows your target glucose range. This is not related to glucose alarm levels.

**High Glucose Alarm Level** - Your High Glucose Alarm level displays only when you have used the App to start a FreeStyle Libre 2 Sensor and have turned the alarm **ON**.

**Low Glucose Alarm Level** - Your Low Glucose Alarm level displays only when you have used the App to start a FreeStyle Libre 2 Sensor and have turned the alarm **ON**.

## Sensor Kit

The Sensor Kit includes:

- Sensor Pack
- Sensor Applicator
- Product insert

When opening your kit, check that the contents are undamaged and that you have all parts listed. If any parts are missing or damaged, contact Customer Service. The Sensor (only visible after applied) is initially in two parts: one part is in the Sensor Pack and the other part is in the Sensor Applicator. Once prepared and applied to your body, the Sensor measures your glucose using a small, flexible tip that inserts just under the skin.

**Sensor Pack.** Used with the Sensor Applicator to prepare the Sensor for use.

**Sensor Applicator.** Applies the Sensor to your body.

## App Setup

Before using the App for the first time, you must complete the setup.

1. Check that your smartphone is connected to a network (WiFi or mobile). You can then install FreeStyle LibreLink from the Google Play store. Tap the App icon to open the App.

**Note:** You only need to be connected to a network for setup, using LibreView, and sharing with other apps. You do not need to be connected to scan a Sensor, add notes, or review your history in the App.

2. Swipe left to view some helpful tips or tap **GET STARTED NOW** at any point.
3. Confirm your country and tap **NEXT**.
4. You need a LibreView account to use the App. Follow onscreen instructions to review legal information and create a new account or login to your existing account.

LibreView Data Management Software is developed by Newyu, Inc. Use of FreeStyle LibreLink requires registration with LibreView, a service provided by Abbott and Newyu, Inc.

5. Confirm your glucose unit of measure and tap **NEXT**.
6. Select how you count carbohydrates (in grams or portions) and tap **NEXT**. The carbohydrate unit will be used in any food notes you enter in the App.
7. Select whether you want sound and vibration OR vibration only when you scan your Sensor. Tap **NEXT**.

**Note:** This setting does not affect alarms.

8. The App now displays some useful information. Tap **NEXT** to review each screen.
9. Apply a new Sensor and then tap **NEXT**. Go to [Starting Your Sensor](#).


**Note:** If you need help applying your Sensor, tap **HOW TO APPLY A SENSOR** or go to [Applying Your Sensor](#).

## Applying Your Sensor

### CAUTION:


- The Sensor Pack and Sensor Applicator are packaged together as a set and have the same Sensor code. Check that the Sensor codes match before using your Sensor Pack and Sensor Applicator. Sensor Packs and Sensor Applicators with the same Sensor code should be used together or your Sensor glucose readings may be incorrect.

- Intense exercise may cause your Sensor to loosen due to sweat or movement of the Sensor. If your Sensor comes loose, you may get no readings or unreliable readings, which may not match how you feel. Follow the instructions to select an appropriate application site.
- 

1. Apply Sensors only on the back of your upper arm. Avoid areas with scars, moles, stretch marks or lumps. Select an area of skin that generally stays flat during your normal daily activities (no bending or folding). Choose a site that is at least 2.5 cm (1 inch) away from an insulin injection site. To prevent discomfort or skin irritation, you should select a different site other than the one most recently used.

2. Wash application site using a plain soap, dry, and then clean with an alcohol wipe. This will help remove any oily residue that may prevent the Sensor from sticking properly. Allow site to air dry before proceeding.

**Note:** The area **MUST** be clean and dry, or the Sensor may not stick to the site.


3. Open the Sensor Pack by peeling the lid off completely. Unscrew the cap from the Sensor Applicator and set the cap aside.

**CAUTION:** Do NOT use if the Sensor Pack or Sensor Applicator seem to be damaged or already opened. Do NOT use if past expiry date.


4. Line up the dark mark on the Sensor Applicator with the dark mark on the Sensor Pack. On a hard surface, press down firmly on the Sensor Applicator until it comes to a stop.

5. Lift the Sensor Applicator out of the Sensor Pack.


6. The Sensor Applicator is prepared and ready to apply the Sensor.

**CAUTION:** The Sensor Applicator now contains a needle. Do NOT touch inside the Sensor Applicator or put it back into the Sensor Pack.


7. Place the Sensor Applicator over the prepared site and push down firmly to apply the Sensor to your body.

**CAUTION:** Do NOT push down on Sensor Applicator until placed over prepared site to prevent unintended results or injury.

8. Gently pull the Sensor Applicator away from your body. The Sensor should now be attached to your skin.

**Note:** Applying the Sensor may cause bruising or bleeding. If there is bleeding that does not stop, remove the Sensor, and apply a new one at a different site.

9. Make sure Sensor is secure after application. Put the cap back on the Sensor Applicator. Discard the used Sensor Applicator and Sensor Pack. See [Disposal](#).

**Note:** You can access an in-app tutorial on applying a Sensor. Tap  at the top of your screen and then tap **Help**.

## Starting Your Sensor

### IMPORTANT:

- The App requires that your smartphone be set to network time. This is the default time setting for most smartphones.
- When using the App, you should keep your smartphone well charged and be sure you have access to a blood glucose meter.
- Make sure your smartphone's sound is set so you can hear scan sounds or

alarms if you turned those on.

- Keep in mind that the ease of scanning a Sensor may vary between devices. Once you determine the location of the NFC (Near Field Communication) antenna on your smartphone, you should be able to reliably scan your Sensor by holding that area near your Sensor. You may need to adjust your scan distance based on what clothing you are wearing. In addition to proximity and orientation, other factors can affect NFC performance. For example, a bulky or metallic case can interfere with the NFC signal.
- 

1. Hold the back of your smartphone near the Sensor (this can be done over clothing) and position the NFC antenna over the Sensor. Do not move the smartphone until you hear the first tone and/or feel a vibration. This indicates your smartphone and Sensor have established an NFC connection.
2. Continue holding your smartphone near the Sensor until you hear a second tone and/or feel a vibration. This completes the scan.

#### Note:

- If you need help, tap **HOW TO SCAN A SENSOR** to view an in-app tutorial. You can also access this later by tapping  at the top of your screen and then tapping **Help**.
- If your Sensor is not successfully scanned, you may receive one of these Scan Errors:
  - The smartphone was unable to scan the Sensor. Scan the Sensor again. Make sure to hold the back of the smartphone near the Sensor. When you hear the first tone or vibration, hold the smartphone still until you hear the second tone or vibration.
  - Open the app to scan your Sensor because another app that uses NFC was detected.

See [Troubleshooting](#) for additional error messages.

3. The Sensor can be used to check your glucose after 60 minutes. While the Sensor is starting up you can navigate away from the App. You will see a notification when the Sensor is ready.


## Note:

- If you like, you can use a Sensor with both the App and the Reader. To do this, you must start the Sensor with the Reader first and then scan with the App. If you start a FreeStyle Libre 2 Sensor with the FreeStyle Libre 2 Reader, remember that you will only get alarms from the Reader. The App can only issue alarms if you use it to start a FreeStyle Libre 2 Sensor.
- Performance between the Reader and App may vary depending on your Reader's software version. Please refer to the performance data insert that's included in the Reader Kit for Reader performance information.
- Remember that FreeStyle LibreLink and Readers do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.

## Checking Your Glucose

1. Open the App and hold the back of your smartphone near the Sensor. If scan sounds are turned on, you will get two separate tones along with the vibrations when your Sensor has been scanned.
2. The My Glucose screen now displays your glucose reading. It includes your Current Glucose, a Glucose Trend Arrow indicating which way your glucose is going, and a graph of your current and stored glucose readings.

**Message** - Tap for more information.

**Back** - Tap to return to Home screen.

**Current Glucose** - Glucose value from your latest scan.

**Add Note** - Tap to add notes to the glucose reading.

**Glucose Trend Arrow** - Direction your glucose is going.

**Note Symbol** - Tap to review notes you've entered.

**Glucose Graph** - Graph of your current and stored glucose readings.

**Target Glucose Range** - The graph shows your target glucose range. This is not related to glucose alarm levels.

**High Glucose Alarm Level** - Your High Glucose Alarm level displays only when you have used the App to start a FreeStyle Libre 2 Sensor and have turned the alarm **ON**.

**Low Glucose Alarm Level** - Your Low Glucose Alarm level displays only when you have used the App to start a FreeStyle Libre 2 Sensor and have turned the alarm **ON**.

### Note:

- A Sensor can store up to 8 hours of glucose data, so scan it at least once every 8 hours to capture all of your available glucose data.
- The graph will scale to 27.8 mmol/L to accommodate glucose readings above 21 mmol/L.
- The ⌚ symbol may appear, indicating the smartphone's time was changed. Gaps in the graph may result or glucose readings may be hidden.
- All available glucose data is used to make your graph so you can expect to see some differences between the graph line and previous current glucose readings.
- Your current glucose value determines the background colour on the My Glucose screen:

**Orange** - High glucose (above 13.3 mmol/L)

**Yellow** - Between the Target Glucose Range and high or low glucose level

**Green** - Within the Target Glucose Range

**Red** - Low glucose (below 3.9 mmol/L)

## Understanding Your Glucose Readings

### Glucose Trend Arrow

The Glucose Trend Arrow gives you an indication of the direction your glucose is going.

 Glucose is rising quickly  
(more than 0.1 mmol/L per minute)

 Glucose is rising (between 0.06 and 0.1 mmol/L per minute)

 Glucose is changing slowly  
(less than 0.06 mmol/L per minute)


 Glucose is falling (between 0.06 and 0.1 mmol/L per minute)

 Glucose is falling quickly  
(more than 0.1 mmol/L per minute)


### Messages

Below are messages you may see with your glucose readings.

**LO | HI:** If **LO** appears, your reading is lower than 2.2 mmol/L. If **HI** appears, your reading is higher than 27.8 mmol/L. You can touch  for more information. Check your blood glucose on your finger with a test strip. If you get a second **LO** or **HI** result, contact your health care professional **immediately**.


**Low Glucose | High Glucose:** If your glucose is higher than 13.3 mmol/L or lower than 3.9 mmol/L, you will see a message on the screen. You can touch  for more information and set a reminder to check your glucose.

---

**Glucose Going Low | Glucose Going High:** If your glucose is projected to be higher than 13.3 mmol/L or lower than 3.9 mmol/L, within 15 minutes, you will see a message on the screen. The background colour corresponds to your current glucose value. You can touch  for more information and set a reminder to check your glucose.

### Note:

- If you are not sure about a message or reading, contact your health care professional for information.
- Messages you receive with your glucose readings are not related to glucose alarm settings.

## Alarms with a FreeStyle Libre 2 Sensor

If you use the App to start a FreeStyle Libre 2 Sensor, you can get Low and High Glucose Alarms from the Sensor if you turn them **ON**. These alarms are turned **OFF** by default.

This section explains how to turn on and set alarms as well as how to use them. Please read all the information in this section before setting and using alarms.

### CAUTION:

- If you are using a FreeStyle Libre Sensor with the App or you started your FreeStyle Libre 2 Sensor with the Reader, you will not receive alarms from the App.

### No Alarms from the App

You are using a FreeStyle Libre Sensor.

You started a FreeStyle Libre 2 Sensor with the Reader before using it with the App.

## Alarms from the App

You started a FreeStyle Libre 2 Sensor with the App.

- You will only get alarms from the App if you use the App to start a FreeStyle Libre 2 Sensor. For you to receive alarms, make sure to:
  - Turn alarms **ON** and ensure that your smartphone is within 6 metres (20 ft) of you at all times. The transmission range is 6 metres (20 ft) unobstructed. If you are out of range, you may not receive glucose alarms.
  - Do not force close the App.
  - Check to make sure that you have the correct settings and permissions enabled on your phone to receive alarms.
    - Enable Bluetooth and Notifications for the App. Be sure not to turn on any features or modify phone settings that may interrupt the presentation of notifications.
    - Turn on Lock Screen notifications, Notification sounds, Channel notifications, Sound and Pop-up notifications, and general phone sounds or vibration.
    - Turn off Do Not Disturb mode or manage the settings for your alarms to override Do Not Disturb. Override Do Not Disturb functionality depends on the phone model and Android operating system version you are using.
    - Turn on Location Permissions. Android 6.0 and above require apps to have Location Permissions enabled to connect with Bluetooth devices. If you wish to receive alarms, Location Permissions must be turned on.
    - Turn off Battery Optimisation for the App to allow the App to run in the background and ensure you receive alarms, even when your phone has

low battery.

- You may need to add the App to the list of apps that will not be put to sleep.
- Remember that alarm settings will follow your smartphone's sounds and vibration settings so these should be at a level you can hear to prevent missed alarms.
- You should disconnect headphones when you are not using them as you may not receive audio with alarms.
- If you are using peripheral devices connected to your phone, such as wireless headphones or a smartwatch, you may receive alarms on only one device or peripheral, not all.
- Keep your smartphone well charged and turned on.

### **IMPORTANT:**

- Scan your Sensor often to check your glucose. If you get a Low or High Glucose Alarm, you must obtain a glucose result to determine what to do next.
- The Low and High Glucose Alarms should not be used exclusively to detect low or high glucose conditions. The glucose alarms should always be used along with your current glucose, glucose trend arrow, and glucose graph.
- Low and High Glucose Alarm levels are different from your Target Glucose Range values. Low and High Glucose Alarms tell you when your glucose has passed the level you set in the alarm. Your Target Glucose Range is displayed on glucose graphs in the App and used to calculate your Time In Target.
- Make sure your smartphone is near you. The Sensor itself will not issue alarms.
- If the Sensor is not communicating with the App, you will not receive glucose alarms, and you may miss detecting low glucose or high glucose episodes. You will see the 📌 or 📌 symbol on the screen when the Sensor is not communicating with the App. Make sure the Signal Loss Alarm is on so you will be notified if your Sensor has not communicated with the App for 20 minutes.
- If you see the 📌 or 📌 symbol, this means you are not getting glucose alarms because of one or more of the following:
  - Bluetooth is **OFF**


- App notifications are **OFF**
- Sensor is not communicating with the App
- Lock Screen notifications or Notification sounds are **OFF**
- Channel notifications or Sound and Pop-up notifications are **OFF**
- Battery optimisation is **ON**
- Location Permissions are **OFF**

## Setting Alarms

To set or turn on alarms, go to the Main Menu and tap **Alarms**. Select the alarm you want to turn on and set.

### Low Glucose Alarm


1. The Low Glucose Alarm is off by default. Tap the slider to turn the alarm on.
2. If the alarm is on, you will be notified when your glucose falls below the alarm level, which is initially set to 3.9 mmol/L. Tap to change this value between 3.3 mmol/L and 5.6 mmol/L. Tap **SAVE**.
3. Choose the sound for this alarm. Volume and vibration will match your smartphone settings.
4. Select whether to manage the settings for this alarm to Override Do Not Disturb. Override Do Not Disturb functionality depends on the phone model and Android operating system version you are using.
5. Tap the back button to return to the main alarm settings screen.

## High Glucose Alarm


1. The High Glucose Alarm is off by default. Tap the slider to turn the alarm on.
2. If the alarm is on, you will be notified when your glucose rises above the alarm level, which is initially set to 13.3 mmol/L. Tap to change this value between 6.7 mmol/L and 22.2 mmol/L. Tap **SAVE**.
3. Choose the sound for this alarm. Volume and vibration will match your smartphone settings.
4. Select whether to manage the settings for this alarm to Override Do Not Disturb. Override Do Not Disturb functionality depends on the phone model and Android operating system version you are using.
5. Tap the back button to return to the main alarm settings screen.


## Signal Loss Alarm

1. Tap the slider to turn the alarm on. If the alarm is on, you will be notified when your Sensor has not communicated with the App for 20 minutes and you are not receiving Low or High Glucose Alarms.  
**Note:** The Signal Loss Alarm automatically turns on the first time you turn the Low or High Glucose Alarm on.
2. Choose the sound for this alarm. Volume and vibration will match your smartphone settings.
3. Tap the back button to return to the main alarm settings screen.


## Using Alarms


**Low Glucose Alarm** notifies you if your glucose drops below the level you set. The alarm does not include your glucose reading, so you need to scan your Sensor to check your glucose. Swipe or tap to dismiss the alarm and check your glucose. You will only receive one alarm per low glucose episode.

**High Glucose Alarm** notifies you if your glucose rises above the level you set. The alarm does not include your glucose reading, so you need to scan your Sensor to check your glucose. Swipe or tap to dismiss the alarm and check your glucose. You will only receive one alarm per high glucose episode.

**Signal Loss Alarm** notifies you if your Sensor has not communicated with the App for 20 minutes and you are not receiving Low or High Glucose Alarms. Signal loss could be caused by the Sensor being too far away from your smartphone (over 6 metres (20 ft)) or another issue such as an error or problem with your Sensor. Swipe or tap to dismiss the alarm.

### Note:

- If you ignore an alarm, you will receive it again in 5 minutes if the condition still exists.
- Only your most recent alarms will display on your screen.

## Adding Notes

Notes can be saved with your glucose readings to help you track food, insulin, and exercise. You can also add your own comment.

1. Tap  on the My Glucose screen.
2. Select the tickbox next to the notes you would like to add. After you tick the box, you can add more specific information to your note.
  - Food notes: Enter meal type and grams or portion information
  - Insulin notes: Enter the number of units taken
  - Exercise notes: Enter intensity and duration
3. Tap **DONE** to save your note.

Notes you add are shown on your glucose graph and in your Logbook as symbols. You can review a note by tapping its symbol on your glucose graph or by going to the Logbook. See [Reviewing Your History](#) for more information about the Logbook. To edit a note from the glucose graph, tap the symbol and

then tap on the information you would like to change. Tap **DONE** when you are finished.

Food

---

Insulin (Rapid or Long-acting)

---

Exercise

---

Food + insulin

---

Multiple/Custom notes – indicates different types of notes entered together or notes entered within a short period of time. A numbered badge next to the symbol indicates the number of notes.

---

## Reviewing Your History

Reviewing and understanding your glucose history can be an important tool for improving your glucose control. The App stores about 90 days of information and has several ways to review your past glucose readings and notes. From the Main Menu, tap **Logbook** to view the Logbook or tap on one of the other history options under **Reports**.

### **IMPORTANT:**

- Work with your health care professional to understand your glucose history.
  - Remember that FreeStyle LibreLink and Readers do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.
-

## Logbook

The Logbook contains entries for each time you scanned your Sensor as well as notes you added. If you would like to view a different day, tap the  symbol or use the arrows. To add a note to a Logbook entry, tap on the entry and then tap . Select your note information and tap **DONE**.

To add a note that is independent of a Logbook entry, tap  on the main Logbook screen. Tap  if you want to add a note on a different date.

## Other History Options

**Daily Patterns:** A graph showing the pattern and variability of your Sensor glucose readings over a typical day. The thick black line shows the median (midpoint) of your glucose readings. The light blue shading represents the 10th - 90th percentile range of your glucose readings. Dark blue shading represents the 25th - 75th percentile range.

**Note:** Daily Patterns needs at least 5 days of glucose data.

**Time In Target:** A graph showing the percentage of time your Sensor glucose readings were above, below, or within your Target Glucose Range.

**Low Glucose Events:** Information about the number of low glucose events measured by your Sensor. A low glucose event is recorded when your Sensor glucose reading is lower than 3.9 mmol/L for longer than 15 minutes. The total number of events is displayed below the graph. The bar graph displays the low glucose events in different periods of the day.

**Average Glucose:** Information about the average of your Sensor glucose readings. The overall average for the selected time period is displayed below the graph. The average is also shown for different periods of the day. Readings above or below your Target Glucose Range are yellow, orange, or red. Readings in range are green.

**Daily Graph:** A graph of your Sensor glucose readings by day. The graph shows your Target Glucose Range and symbols for notes you have entered.

- The graph will scale to 27.8 mmol/L to accommodate glucose readings above 21 mmol/L.

- You might see gaps in the graph during times when you have not scanned at least once in 8 hours.
- The  symbol may appear indicating a time change. Gaps in the graph may result or glucose readings may be hidden.

**Estimated A1c:** Your estimated A1c level (also called HbA1c) is based on available Sensor glucose data from the last 90 days. The more data available, the better your estimation will be. However, the estimated level may not match your A1c measured in a laboratory \*. A1c can be used as an indicator of how well your glucose levels have been controlled and may be used to monitor your diabetes treatment regimen.

\* The formula is based on the published reference, which compared average Sensor glucose and laboratory-measured A1c:

$$A1c_{\%} = (\text{Avg } SG_{\text{mg/dL}} + 46.7)/28.7$$

$$A1c_{\%} = (\text{Avg } SG_{\text{mmol/L}} + 2.59)/1.59$$

Reference: Nathan DM, Kuenen J, Borg R, Zheng H, Schoenfeld D, Heine RJ for the A1c-Derived Average Glucose (ADAG). Study Group: Translating the hemoglobin A1c assay into estimated average glucose values. Diabetes Care 2008, 31:1473-8.

**Sensor Usage:** Information about how often you scan your Sensor. This includes the total number of scans, an average of how many times you scanned your Sensor each day, and the percentage of possible Sensor data recorded from your scans.

### Note:

- Tap the  symbol on any report to share a screenshot of the report.
- Tap the  symbol to view a description of the report.
- From any report screen, swipe left or right to view the next or previous report.
- On all reports except the Daily Graph and Estimated A1c, you can select to show information about your last 7, 14, 30, or 90 days.

## Removing Your Sensor

1. Pull up the edge of the adhesive that keeps your Sensor attached to your skin. Slowly peel away from your skin in one motion.

**Note:** Any remaining adhesive residue on the skin can be removed with warm soapy water or isopropyl alcohol.

2. Discard the used Sensor. See [Disposal](#). When you are ready to apply a new Sensor, follow the instructions in [Applying Your Sensor](#) and [Starting Your Sensor](#). If you removed your last Sensor before 14 days of use, you will be prompted to confirm that you would like to start a new Sensor when you first scan it.

## Replacing Your Sensor

Your Sensor automatically stops working after 14 days of wear and must be replaced. You should also replace your Sensor if you notice any irritation or discomfort at the application site or if the App reports a problem with the Sensor currently in use. Taking action early can keep small problems from turning into larger ones.

**CAUTION:** If the glucose readings from the Sensor do NOT seem to match with how you feel, check to make sure that your Sensor has not come loose. If the Sensor tip has come out of your skin, or your Sensor is coming loose, remove the Sensor and apply a new one.

## Setting Reminders

You can create single or repeating reminders to help you remember things like checking your glucose or taking insulin. There is one default reminder to help you remember to scan your Sensor. This Scan Sensor reminder can be changed or disabled but cannot be deleted.

**Note:** If you want to receive a sound/vibration with your reminder, ensure that sound/vibration on your smartphone is turned on, sound is set at a level you can

hear, and your smartphone's Do Not Disturb feature is turned off. If Do Not Disturb is on, you will only see your reminder on the screen.

1. To add a new reminder, go to the Main Menu and tap **Reminders**. Tap **ADD REMINDER**.
2. Name your reminder.
3. Tap the time fields to set the time for the reminder.

The screenshot shows a reminder creation interface for an event titled "Exercise". The time is set to 16:00. A "Repeating" toggle switch is turned on. Below this, there are two columns of checkboxes, each with a blue checkmark, indicating that all days of the week are selected for the reminder to repeat. The days listed are All, Monday, Tuesday, Wednesday, Thursday, Friday, Saturday, and Sunday. At the bottom of the screen, there are two buttons: "CANCEL" on the left and "DONE" on the right.

**Note:** If you would like the reminder to repeat, tap the slider to the right. You can also select which days you would like to receive the reminder.

4. Tap **DONE**. You will now see your reminder on the list along with the time you will receive it.

**Note:**

- To turn off a reminder, tap the slider to the left.
- To delete a reminder, swipe the reminder to the right and tap the  symbol. The Scan Sensor reminder cannot be deleted.
- Your reminders will be received as notifications that you can swipe or tap to dismiss.

## App Settings and Other Options

You can go to the Overflow Menu to change App settings and also view information about the App.

1. Tap the Overflow menu  at the top right of your screen.
2. Tap **App Settings** and select your setting. Tap **SAVE** when you are done.

**Unit of Measurement** - View the glucose unit of measure used in the App.

**Report Settings** - Work with your health care professional to set your Target Glucose Range, which is displayed on glucose graphs in the App and used to calculate your Time In Target. The Target Glucose Range setting will not set glucose alarm levels.

**Carbohydrate Units** - Choose grams or portions for food notes that you enter.

**Scan Sounds** - Select whether you would like to hear a sound in addition to a vibration when you scan the Sensor. Remember that Scan Sounds inherit the volume settings on your smartphone. If your smartphone volume is turned off, you will not hear a scan sound. The Scan Sounds setting does not affect alarms.

**Text to Speech** - Turn on Text to Speech to have the glucose reading read aloud when you scan the Sensor. You will hear only your current glucose value and trend arrow direction. Additional information, such as the glucose graph and any message, are available on your My Glucose screen. Always review your My Glucose screen to get complete information. Remember that Text to Speech inherits the volume settings on your smartphone. If your smartphone volume is turned off, you will not hear the glucose reading read aloud.

### Other options:

**Account Settings:** View/change your LibreView account information.

**Account Password:** Change your LibreView account password.

**Help:** View in-app tutorials, access this User's Manual, and review the App's legal information. You can also view the Event Log, which is a list of events recorded by the App. This may be used by Customer Service to help

troubleshoot.

**About:** View App software version and other information.

## Using the Connected Apps Option

The **Connected Apps** option in the Main menu opens a web browser within the App. It lists different apps you can connect with to share your data. The available apps may vary based on your country. To connect your data with apps listed in the **Connected Apps** option, select them from the list of apps, and follow the onscreen instructions.

## Living with your Sensor

### Activities

**Bathing, Showering, and Swimming:** Your Sensor is water-resistant and can be worn while bathing, showering, or swimming. Do NOT take your Sensor deeper than 3 feet (1 metre) or immerse it longer than 30 minutes in water.

**Sleeping:** Your Sensor should not interfere with your sleep. It is recommended that you scan your Sensor before going to sleep and when you wake up because your Sensor holds 8 hours of data at a time. If you have reminders set to go off while you are sleeping or have glucose alarms set, place your smartphone nearby.

### Travelling by Air:

- You may use your Sensor while on an aircraft, following any requests from the flight crew. After you have put your smartphone in airplane mode, you can turn NFC back on to continue to get Sensor glucose readings.

**IMPORTANT:** Glucose alarms (if available) will not be issued while your smartphone is in airplane mode unless you enable Bluetooth.

- Some airport full-body scanners include x-ray or millimetre radio-wave, which you cannot expose your Sensor to. The effect of these scanners has not been evaluated and the exposure may damage the Sensor or cause inaccurate results. To avoid removing your Sensor, you may request another type of screening. If you do choose to go through a full-body scanner, you must remove your Sensor.

- The Sensor can be exposed to common electrostatic (ESD) and electromagnetic interference (EMI), including airport metal detectors.

**Note:** Changing the time affects the graphs and statistics. The ⌚ symbol may appear on your glucose graph indicating a time change. Gaps in the graph may result or glucose readings may be hidden.

## Maintenance

The Sensor has no serviceable parts.

## Disposal

### Reader and Sensor:

These devices must not be disposed of via municipal waste collection. Separate collection for electrical and electronic equipment waste per Directive 2012/19/EC in the European Union is required. Contact the manufacturer for details.

As Readers and Sensors may have been exposed to bodily fluids, you may wipe prior to disposing, such as by using a cloth dampened with a mixture of one part household bleach and nine parts water.

**Note:** Readers and Sensors contain non-removeable batteries and must not be incinerated. Batteries may explode upon incineration.

### Sensor Applicator:

Please consult your local waste management authority for instructions on how to dispose Sensor Applicators at a designated sharps collection site. Ensure the cap is on the Sensor Applicator as it contains a needle.

### Sensor Pack:

Used Sensor Packs may be disposed of via municipal waste collection.

## Troubleshooting

This section lists problems that you may experience, the possible cause(s), and recommended actions. If there is an error, a message will appear on the screen with directions to resolve the error.

**IMPORTANT:** If you are having issues with the App, please keep in mind that uninstalling the App and/or clearing data will cause you to lose all historical data

and end the Sensor currently in use. Please call Customer Service if you have any questions.

## Problems at the Sensor Application Site

**Problem: The Sensor is not sticking to your skin.**

What it may mean: The site is not free of dirt, oil, hair, or sweat.

What to do: 1. Remove the Sensor. 2. Clean the site with a plain soap and water and consider shaving. 3. Follow the instructions in [Applying Your Sensor](#) and [Starting Your Sensor](#).

---

**Problem: Skin irritation at the Sensor application site.**

What it may mean: Seams or other constrictive clothing or accessories causing friction at the site **OR** you may be sensitive to the adhesive material.

What to do: Ensure that nothing rubs on the site. If the irritation is where the adhesive touches skin, contact your health care professional to identify the best solution.

## Problems Starting Your Sensor or Receiving Sensor Readings

**Display: Sensor Starting Up**

What it may mean: Sensor is not ready to read glucose.

What to do: Wait until the 60 minute Sensor start-up period has completed.

---

**Display: Signal Loss Alarm**

What it may mean: Sensor has not automatically communicated with the App in the last 20 minutes.

What to do: Make sure your smartphone is within 6 metres (20 ft) of the Sensor. Try scanning the Sensor to get a glucose reading. If the Signal Loss Alarm shows again after scanning your Sensor, contact Customer Service.

---

Display: **Sensor Ended**

What it may mean: The Sensor life has ended.

What to do: Apply and start a new Sensor.

---

Display: **New Sensor Found**

What it may mean: You scanned a new Sensor before your previous Sensor ended.

What to do: Your smartphone can only be used with one Sensor at a time. If you start a new Sensor, you will no longer be able to scan your old Sensor. If you would like to begin using the new Sensor, select "Yes".

---

Display: **Sensor Error**

What it may mean: The Sensor is unable to provide a glucose reading.

What to do: Scan again in 10 minutes.

---

Display: **Glucose Reading Is Unavailable**

What it may mean: The Sensor is unable to provide a glucose reading.

What to do: Scan again in 10 minutes.

---

Display: **Sensor Too Hot**

What it may mean: Your Sensor is too hot to provide a glucose reading.

What to do: Move to a location where the temperature is appropriate and scan again in a few minutes.

---

**Display: Sensor Too Cold**

What it may mean: Your Sensor is too cold to provide a glucose reading.

What to do: Move to a location where the temperature is appropriate and scan again in a few minutes.

---

**Display: Check Sensor**

What it may mean: The Sensor tip may not be under your skin.

What to do: Try to start your Sensor again. If you see "Check Sensor" again on the screen, your Sensor was not applied properly. Apply and start a new Sensor.

---

**Display: Replace Sensor**

What it may mean: The App has detected a problem with your Sensor.

What to do: Apply and start a new Sensor.

---

**Display: Unexpected Application Error**

What it may mean: The App has detected an unexpected error.

What to do: Shut down the App completely and restart.

---

**Display: Incompatible Sensor**

What it may mean: The Sensor cannot be used with the App.

What to do: Call Customer Service.

---

**Display: Scan Error**

What it may mean: The smartphone was unable to scan the Sensor **OR** another NFC app is competing for your smartphone's NFC.

What to do: Try scanning the Sensor again. Make sure to hold your smartphone still once you get the first tone and/or vibration. Wait until you get the second tone and/or vibration to move it away from the Sensor. Make sure you are not

touching any buttons on the smartphone or the screen **OR** open the app to scan your Sensor because another app that uses NFC was detected.

## Problems Receiving Glucose Alarms

What it may mean: You have not turned glucose alarms on.

What to do: Go to the main menu and then select **Alarms**. Choose the alarm you want to turn on and set.

---

What it may mean: You are using a FreeStyle Libre Sensor or you started a FreeStyle Libre 2 Sensor with the Reader before using it with the App.

What to do: Start a new FreeStyle Libre 2 Sensor with the App.

---

What it may mean: The Sensor is not communicating with the App or there may be a problem with the Sensor.

What to do: The Sensor must be within range (6 metres (20 ft)) of your smartphone for you to receive alarms. Make sure that you are within this range. You will see the  or  symbol when your Sensor has not communicated with the App in 5 minutes. If the Signal Loss Alarm is on, you will be notified if there has been no communication for 20 minutes. Try scanning your Sensor. If the Signal Loss Alarm is on and shows again after scanning your Sensor, contact Customer Service.

---

What it may mean: One or more of the following is turned off: Bluetooth, Notifications, Lock Screen notifications, Notification sounds, Location Permissions, Channel notifications, Sound and Pop-up notifications, or general phone sounds or vibration. Or, you have turned on Battery Optimisation. Or, you have turned on Do Not Disturb mode without managing the settings for your alarms to override Do Not Disturb.

What to do: Check to make sure that you have the correct settings and permissions enabled on your phone to receive alarms. Go to [Setting Alarms](#) for

more information.

---

What it may mean: You may have set an alarm level that is higher or lower than you intended.

What to do: Confirm your alarm settings are appropriate.

---

What it may mean: If you are using peripherals such as wireless headphones or a smartwatch you may receive alarms on only one device or peripheral, not all.

What to do: Disconnect headphones or peripherals when you are not using them.

---

What it may mean: The App was put to sleep by the phone operating system.

What to do: Put the App on the list of apps that will not be put to sleep.

---

What it may mean: You have already dismissed this type of alarm.

What to do: You will receive another alarm when a new low or high glucose episode starts.

---

What it may mean: You have closed the App.

What to do: Make sure the App is always open in the background.

---

What it may mean: Your Sensor has ended.

What to do: Replace your Sensor with a new one.

---

## Customer Service

Customer Service is available to answer any questions you may have about FreeStyle LibreLink. Please go to [www.FreeStyleLibre.com](http://www.FreeStyleLibre.com) or refer to the product

insert in your Sensor Kit for your Customer Service phone number. A printed copy of this User's Manual is available upon request.

## Reporting of Serious Incidents

If a serious incident has occurred in relation to this device it should be reported to Abbott Diabetes Care. Please go to [www.FreeStyleLibre.com](http://www.FreeStyleLibre.com) or refer to the product insert in your Sensor Kit for your Customer Service phone number.

In European Union Member States, serious incidents should also be reported to the competent authority, (the government department responsible for medical devices) in your country. Please refer to your government website for details of how to contact your competent authority.

A 'serious incident' means any incident that directly or indirectly led, might have led or might lead to:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health

## Sensor Specifications

**Sensor glucose assay method:** Amperometric electrochemical sensor

**Sensor glucose reading range:** 2.2 to 27.8 mmol/L

**Sensor size:** 5mm height and 35mm diameter

**Sensor weight:** 5 grams

**Sensor power source:** One silver oxide battery

**Sensor life:** Up to 14 days

**Sensor memory:** 8 hours (glucose readings stored every 15 minutes)

**Operating temperature:** 10°C to 45°C

**Sensor Applicator and Sensor Pack storage temperature:** 4°C to 25°C

**Operating and storage relative humidity:** 10-90%, non-condensing

**Sensor water resistance and ingress protection:** IP27: Can withstand immersion into one metre (3 ft) of water for up to 30 minutes. Protected against insertion of objects > 12mm diameter.

**Operating and storage altitude:** -381 metres (-1,250 ft) to 3,048 metres (10,000 ft)

**Radio Frequency (FreeStyle Libre 2 Sensor):** 2.402-2.480 GHz BLE; GFSK; 0dBm EIRP

**Sensor transmission range (FreeStyle Libre 2 Sensor):** 6 metres (20 ft) unobstructed

## Labelling Symbols and Definitions

	Consult instructions for use
	Temperature limit
	Manufacturer
	Date of Manufacture
	CE Mark
	Authorised Representative in the European Community
	Single sterile barrier system
	Batch code
	Type BF applied part
<b>CODE</b>	Sensor code
	Do not re-use
	Use-by date
	Catalogue number

**SN** Serial number

---

 Caution

---

**STERILE R** Sterilised by irradiation

---

Sterile Barrier. Refer to Instructions for Use if opened or damaged.

---

 Humidity limitation

---

Do not use if package is damaged.

 **For Sterile Barrier:** Do not use if the product sterile barrier system or its packaging is compromised.

---

 This product must not be disposed of via municipal waste collection. Separate collection for electrical and electronic equipment waste per Directive 2012/19/EC in the European Union is required. Contact the manufacturer for details.

---

## Electromagnetic Compatibility

- The Sensor needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual.
- Portable and mobile RF communications equipment can affect the Sensor.
- The use of accessories, transducers and cables other than those specified by Abbott Diabetes Care

may result in increased EMISSIONS or decreased IMMUNITY of the Sensor.

- The Sensor should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the Sensor should be observed to verify normal operation in the configuration in which it will be used.

## **Guidance and manufacturer's declaration – electromagnetic emissions**

The Sensor is intended for use in the electromagnetic environment specified below. The customer or the user of the Sensor should assure that it is used in such an environment.

Emissions test: RF emissions; CISPR 11

Compliance: Group 1

Electromagnetic environment – guidance: The Sensor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Emissions test: RF emissions; CISPR 11

Compliance: Class B

Electromagnetic environment – guidance: The Sensor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.

## **Guidance and manufacturer's declaration – electromagnetic immunity**

The Sensor is intended for use in the electromagnetic environment specified below. The customer or the user of the Sensor should assure that it is used in such an environment.

Immunity test: Electrostatic discharge (ESD); IEC 61000-4-2

IEC 60601 test level:  $\pm 8$  kV contact;  $\pm 15$  kV air

Compliance level:  $\pm 8$  kV contact;  $\pm 15$  kV air

Electromagnetic environment – guidance: Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Immunity test: Power frequency (50/60 Hz); magnetic field; IEC 61000-4-8

IEC 60601 test level: 30 A/m

Compliance level: 30 A/m

Electromagnetic environment – guidance: Power frequency magnetic fields should be at levels characteristic of a typical location in a typical domestic, commercial, or hospital environment.

Immunity test: Radiated RF; IEC 61000-4-3

IEC 60601 test level: 10 V/m; 80 MHz to 2.7 GHz

Compliance level: 10 V/m

Electromagnetic environment – guidance:

Recommended separation distance

$$d = 1.2 \sqrt{P}$$

80 MHz to 800 MHz

$$d = 2.3 \sqrt{P}$$

800 MHz to 2.5 GHz

$P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and  $d$  is the recommended separation distance in metres (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup>

Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (mobile/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Sensor is used exceeds the applicable RF compliance level above, the Sensor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Sensor.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the Sensor

The Sensor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Sensor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Sensor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.


## Performance Characteristics

**Note:** Please consult your health care team on how to use the information in this section.

### Performance Characteristics

Performance of the Sensor was evaluated in a controlled clinical study. The study was conducted in 5 centres and a total of 146 subjects with diabetes were included in the effectiveness analysis. Each subject wore up to two Sensors for up to 14 days on the back of the upper arm. During the study, subjects had their venous blood glucose analysed over three separate visits to the clinical centre using the Yellow Springs Instrument Life Sciences 2300 STAT Plus.™ Three lots of Sensors were evaluated in the study.

**Fig 1.** Comparison of the Sensors vs. YSI reference.

**Table 1.** Regression analysis of the Sensors vs. YSI reference

<b>Slope</b>	0.97
<b>Intercept</b>	-1.3 mg/dL (-0.1 mmol/L)
<b>Correlation</b>	0.98
<b>N</b>	18926
<b>Range</b>	37 - 479 mg/dL (2.0 - 26.6 mmol/L)
<b>Overall mean bias</b>	-5.6 mg/dL (-0.3 mmol/L)
<b>Mean Absolute Relative Difference (MARD)</b>	9.2%

**Table 2.** Sensor accuracy for all results vs. YSI reference

Sensor accuracy results for glucose concentrations	Within $\pm 15$ mg/dL (within $\pm 0.83$ mmol/L)	Within $\pm 20$ mg/dL (within $\pm 1.11$ mmol/L)	Within $\pm 30$ mg/dL (within $\pm 1.67$ mmol/L)
Sensor accuracy results for glucose concentrations <80 mg/dL (4.4 mmol/L)	4199 / 4595 (91.4%)	4482 / 4595 (97.5%)	4583 / 4595 (99.7%)
Sensor accuracy results for glucose concentrations $\geq 80$ mg/dL (4.4 mmol/L)	Within $\pm 15\%$	Within $\pm 20\%$	Within $\pm 30\%$
	12143 / 14331 (84.7%)	13153 / 14331 (91.8%)	14012 / 14331 (97.8%)
Sensor accuracy for all results	Within $\pm 20$ mg/dL ( $\pm 1.11$ mmol/L) and within $\pm 20\%$ of reference		
	17635 / 18926 (93.2%)		

**Table 3.** Sensor performance relative to YSI reference at different glucose levels

Glucose	Mean Absolute Relative Difference
≤50 mg/dL (2.8 mmol/L)	9.1 mg/dL (0.5 mmol/L)*
51-80 mg/dL (2.8-4.4 mmol/L)	7.0 mg/dL (0.4 mmol/L)*
81-180 mg/dL (4.5-10.0 mmol/L)	10.1%
181-300 mg/dL (10.0-16.7 mmol/L)	7.5%
301-400 mg/dL (16.7-22.2 mmol/L)	7.1%
>400 mg/dL (22.2 mmol/L)	10.2%

\* For glucose ≤80 mg/dL (4.4 mmol/L), the differences in mg/dL (mmol/L) are presented instead of relative differences (%).

**Table 4.** Sensor accuracy over wear duration vs. YSI reference

	Beginning	Early Middle	Late Middle	End
<b>Within ±20 mg/dL (±1.11 mmol/L) and within ±20% of reference</b>	91.2%	95.1%	94.2%	93.7%
<b>Mean Absolute Relative Difference (%)</b>	10.0	8.5	8.8	9.1

## Skin Interaction

Based on the examination of 146 study participants, the following incidence of skin issues were observed. Four occurrences of erythema were reported to be moderate in intensity. All other skin issues were reported to be mild in intensity.

Bleeding – 0.7% of the time

Bruising – 0.7% of the time

Erythema – 2.7% of the time

Pain – 0.7% of the time

Scabbing – 2.7% of the time

## Expected Clinical Benefits

Complications as a result of diabetes mellitus (including, but not limited to: diabetic retinopathy, diabetic nephropathy) are well documented.<sup>1</sup> Self-monitoring of blood glucose (SMBG) by patients has revolutionised management of diabetes.<sup>2</sup> Using glucose monitoring devices patients with diabetes can work to achieve and maintain specific glycaemic goals. Given the results of the Diabetes Control and Complications Trial (DCCT)<sup>3</sup> and other studies, there is broad consensus on the health benefits of normal or near-normal blood glucose levels and on the importance, especially in insulin-treated patients, of glucose monitoring devices in treatment efforts designed to achieve these glycaemic goals. Based principally on the DCCT results, experts recommend that most individuals with diabetes should attempt to achieve and maintain blood glucose levels as close to normal as is safely possible. Most patients with diabetes, especially insulin treated patients, can achieve this goal only by using glucose monitoring devices.

<sup>1</sup>Textbook of Diabetes, Volumes 1 & 2; Pickup and Williams, 1999.

<sup>2</sup>ADA Position Statement. Test of glycemia in diabetes. Diabetes Care 2003; 26(Suppl.1) S106-108.

<sup>3</sup>Diabetes Control and Complications Trial Research Group (DCCT): The effect of intensive treatment of diabetes on the development and progression of long term complications in insulin dependent diabetes mellitus. New Engl J Med, 329: 977-86; 1993.

Customer Service: [www.FreeStyleLibre.com](http://www.FreeStyleLibre.com)

Patent: <https://www.abbott.com/patents>

FreeStyle, Libre, and related brand marks are marks of Abbott. Other trademarks are the property of their respective owners.

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©2021 Abbott
Abbott B.V.
Abbott Diabetes Care Ltd.