


TIPS FOR KIDS





Remember these important things:



The Sensor has to go on the back of my upper arm.



When I see  during the first 12 hours of wearing a Sensor: I must do a blood glucose test with my meter before taking action like eating carbs or dosing insulin.



I must do a blood glucose test with my meter if I think the number doesn't match how I feel.



If I take insulin for a meal, I should not take any more insulin for 2 hours. I need to wait and watch.




Warning!

Don't share your System with another person to avoid confusing glucose information.

Important!


- While using a FreeStyle Libre 2 Sensor or Libre 3 Sensor, don't take more than 500 mg of Vitamin C per day.
- While using a FreeStyle Libre 2 Plus Sensor or Libre 3 Plus Sensor, don't take more than 1000 mg of Vitamin C per day.
- Taking more than these levels may affect your Sensor readings. Some supplements like Airborne® and Emergen-C® contain 1000 mg of Vitamin C.



My Reader

This arrow shows me which way my glucose is going

Graph of my glucose results



My App

