FreeStyle *LibreLink*

User's Manual
Manual del usuario

English | **Español**

Important Safety Information

Online Labeling

Indications for Use

Contraindications

Warnings

Cautions and Limitations

App Symbols

FreeStyle LibreLink Overview

Home Screen

Sensor Kit

App Setup

Applying Your Sensor

Starting Your Sensor

Checking Your Glucose

Understanding Your Glucose Readings

Making Treatment Decisions

Adding Notes

Reviewing Your History

Logbook

Other History Options

Removing Your Sensor

Replacing Your Sensor

Setting Reminders

Settings and Other Options in the Main Menu

Activities

Maintenance and Disposal

Troubleshooting

Problems at the Sensor Application Site

Problems Starting Your Sensor or Receiving Sensor Readings

Customer Service

Labeling Symbols and Definitions

Performance Characteristics

Electromagnetic Compatibility

User's Manual

Important Safety Information

Online Labeling

The latest version of the User's Manual, including performance data, can always be accessed at www.FreeStyleLibre.us/support/overview.html

You can also order a free printed copy from Customer Service: 1-855-632-8658, 7 Days a Week from 8AM to 8PM Eastern Time; excluding holidays.

Indications for Use

The FreeStyle Libre and FreeStyle Libre 14 day Flash Glucose Monitoring Systems are continuous glucose monitoring (CGM) devices indicated for the management of diabetes in persons age 18 and older. They are designed to replace blood glucose testing for diabetes treatment decisions.

Each System detects trends and track patterns aiding in the detection of episodes of hyperglycemia and hypoglycemia, facilitating both acute and long-term therapy adjustments. Interpretation of the Sensor glucose readings should be based on the glucose trends and several sequential readings over time. The Systems are intended for single patient use and require a prescription.


Contraindications

 **MRI/CT/Diathermy:** The System must be removed prior to Magnetic Resonance Imaging (MRI), Computed Tomography (CT) scan, or high-frequency electrical heat (diathermy) treatment. The effect of MRI, CT scans, or diathermy on the performance of the System has not been evaluated. The exposure may damage the Sensor and may impact proper function of the device which could cause incorrect readings.

Warnings

- **Do not ignore symptoms that may be due to low or high blood glucose:** If you are experiencing symptoms that are not consistent with your glucose readings, consult your health care professional.

- **Check Sensor glucose readings by conducting a fingerstick test with a blood glucose meter under the following conditions, when Sensor glucose readings may not be accurate and should not be used to make a diabetes treatment decision:**
 - If you suspect that your reading may be inaccurate for any reason
 - When you are experiencing symptoms that may be due to low or high blood glucose
 - When you are experiencing symptoms that do not match the Sensor glucose readings
 - During the first 12 hours of wearing a FreeStyle Libre 14 day Sensor
 - During times of rapidly changing glucose (more than 2 mg/dL per minute)
 - When the Sensor glucose reading does not include a Current Glucose number or Glucose Trend Arrow
 - In order to confirm hypoglycemia or impending hypoglycemia as reported by the Sensor
 - When you see the  symbol, you must check your blood glucose with a blood glucose meter before making any treatment decisions. Sensor readings may not accurately reflect blood glucose levels.

- If you are using FreeStyle LibreLink app, you must have access to a blood glucose monitoring system as the App does not provide one.
- **Hypoglycemic unawareness:** The System has not been evaluated for use in patients with

hypoglycemic unawareness and will not automatically alert you of a hypoglycemic event without you scanning your Sensor.

- **No alarms without a Sensor scan:** The System does not have alarms that will automatically notify you when you are having a severe low (hypoglycemic) or high (hyperglycemic) glucose event unless you scan your Sensor. For example, the System does not have an alarm that can alert or wake you when you are sleeping in the case of low or high glucose.
- **Choking hazard:** The System contains small parts that may be dangerous if swallowed.

Cautions and Limitations

Below are important cautions and limitations to keep in mind so you can use the System safely. They are grouped into categories for easy reference.

What to know about Alarms/Alerts:

- There are NO alarms or alerts unless you scan the Sensor.

What to know before using the System:

- Review all product information before use.
- Take standard precautions for transmission of blood borne pathogens to avoid contamination.

Who should not use the System:

- **Do not use the System in people younger than 18 years of age.** The System is not approved for use in people under 18 years of age and Sensor readings in this population may be inaccurate. In general, continuous glucose monitoring systems are recognized to be less accurate in children than in adults.
- **Do not use the System in critically ill patients.** The System is not approved for use in these patients. It is not known how different conditions or medications common to the critically ill population may affect performance of the System. Sensor glucose readings may be inaccurate in critically ill patients.
- **Do not use the System in pregnant women or persons on dialysis.** The System is not approved for use in pregnant women or persons on dialysis and has not been evaluated in these populations.
- Performance of the System when used with other implanted medical devices, such as pacemakers, has not been evaluated.

What should you know about wearing a Sensor:

- After the start-up period, the Sensor can be worn for up to the wear duration specified in your Sensor Kit's product insert.
- Some individuals may be sensitive to the adhesive that keeps the Sensor attached to the skin. If you notice significant skin irritation around or under your Sensor, remove the Sensor and stop using the System. Contact your health care professional before continuing to use the System.
- Intense exercise may cause your Sensor to loosen due to sweat or movement of the Sensor. Remove and replace your Sensor if it starts to loosen and follow the instructions to select an appropriate application site.
- If you are wearing a FreeStyle Libre 14 day Sensor, all available glucose data is used to give you readings so you should scan your Sensor at least once every 8 hours for the most accurate performance. Scanning less frequently may result in decreased performance. If you are using both the App and Reader with the same FreeStyle Libre 14 day Sensor, be sure to scan frequently with both devices.
- Do not reuse Sensors. The Sensor and Sensor Applicator are designed for single use. Reuse may result in no glucose readings and infection. Not suitable for re-sterilization. Further exposure to irradiation may cause inaccurate results.
- If a Sensor breaks inside your body, call your health care professional.

How to store the Sensor Kit:

- Store the Sensor Kit between 39°F and 77°F. Storage outside of this range may cause inaccurate Sensor glucose readings. While you don't need to keep your Sensor Kit in a refrigerator, you can as long as the refrigerator is between 39°F and 77°F. Do not freeze.
- Store the Sensor Kit between 10-90% non-condensing humidity.

How to store the Reader:


- Store the Reader between -4°F and 140°F. Storage in temperatures outside of this range, such as in a parked car on a hot day, may cause the Reader to not function properly.

When not to use the System:

- Do NOT use if the Sensor Kit package, Sensor Pack, or Sensor Applicator appear to be damaged or already opened due to risk of no results and/or infection.
- Do NOT use if Sensor Kit contents are past expiration date.
- Do NOT use if the Reader appears to be damaged due to risk of electric shock and/or no results.

What to know before you apply the Sensor:

- The Sensor Pack and Sensor Applicator are packaged as a set (separately from the Reader) and have the same Sensor code. Check that the Sensor codes match before using your Sensor Pack and Sensor Applicator. Do not use Sensor Packs and Sensor Applicators with different Sensor codes together as this will result in incorrect glucose readings.

- Clean the application site and ensure that it is dry prior to Sensor insertion. This helps the Sensor stay attached to your body.
- Clean hands prior to Sensor handling/insertion to help prevent infection.
- Change the application site for the next Sensor application to prevent discomfort or skin irritation.
- Sensor placement is not approved for sites other than the back of the arm. If placed in other areas, the Sensor may not function properly.
- Select an appropriate Sensor site to help the Sensor stay attached to the body and prevent discomfort or skin irritation. Avoid areas with scars, moles, stretch marks, or lumps. Select an area of skin that generally stays flat during normal daily activities (no bending or folding). Choose a site that is at least 1 inch away from an insulin injection site.

When is Sensor Glucose different from Blood Glucose:

- Physiological differences between the interstitial fluid and capillary blood may result in differences in glucose readings between the System and results from a fingerstick test using a blood glucose meter. Differences in glucose readings between interstitial fluid and capillary blood may be observed during times of rapid change in blood glucose, such as after eating, dosing insulin, or exercising.

What to know about interfering substances such as Vitamin C and Aspirin:

- Taking ascorbic acid (vitamin C) while wearing the Sensor may falsely raise Sensor glucose readings. Taking salicylic acid (used in some pain relievers such as aspirin and some skin care products) may slightly lower Sensor glucose readings. The level of inaccuracy depends on the amount of the interfering substance active in the body.

- Test results did not indicate interference for methyldopa (used in some drugs to treat high blood pressure) or tolbutamide (infrequently used in some drugs to treat diabetes in the US) at maximum circulating levels. However, concentrations of potential interferents in interstitial fluid are unknown compared to circulating blood.

What to know about X-Rays:

- The Sensor should be removed prior to exposing it to an X-ray machine. The effect of X-rays on the performance of the System has not been evaluated. The exposure may damage the Sensor and may impact proper function of the device to detect trends and track patterns in glucose values during the wear period.

When to remove the Sensor:

- If the Sensor is becoming loose or if the Sensor tip is coming out of your skin, you may get no readings or unreliable readings, which may not match how you feel. Check to make sure your Sensor has not come loose. If it has come loose, remove it and apply a new one.
- If you believe your glucose readings are not correct or are inconsistent with how you feel, perform a blood glucose test on your finger to confirm your glucose. If the problem continues, remove the current Sensor and apply a new one.

What to do if you are dehydrated:

- Severe dehydration and excessive water loss may cause inaccurate Sensor glucose readings. If you believe you are suffering from dehydration, consult your health care professional immediately.

What to know about the Reader:

- Do NOT place the Reader in water or other liquids as this may cause it to not function properly and may lead to **risk of fire or burns**.
- The FreeStyle Libre and FreeStyle Libre 14 day Readers have a built-in blood glucose meter that is designed to be used only with FreeStyle Precision Neo blood glucose test strips and MediSense Glucose and Ketone Control Solution. Using other test strips with the Reader's built-in meter will produce an error or cause the Reader's built-in meter to not turn on or start a test. The Reader's built-in meter does not have ketone testing functionality.
- The Reader's built-in meter is not for use on people who are dehydrated, hypotensive, in shock, or for individuals in hyperglycemic-hyperosmolar state, with or without ketosis.
- The Reader's built-in meter is not for use on neonates, in critically-ill patients, or for diagnosis or screening of diabetes.
- See *Using the Reader's Built-in Meter* section of the Reader Kit User's Manual for additional

important information on the use of the Reader's built-in meter.

What to know about charging your Reader:

- Always use the Abbott provided power adapter and yellow USB cable that came with your Reader to minimize the **risk of fire or burns**. Take care when plugging and unplugging your USB cable. Do not force or bend the end of the USB cable into the Reader's USB port.
- Choose a location for charging where you can easily access the power adapter and quickly disconnect to prevent the potential risk of electrical shock.
- The maximum surface temperature of the Reader and/or the power adapter could go as warm as 120°F when it's charging or 118°F during normal use. Under these conditions, do not hold the Reader or the power adapter for five minutes or more. People with disorders of peripheral circulation or sensation should use caution at this temperature.
- Do NOT expose the USB cable or power adapter to water or other liquids as this may cause them to not function properly and may lead to **risk of fire or burns**.

What to know about FreeStyle LibreLink:

- FreeStyle LibreLink installed on a smartphone is intended for use by a single person. It must not be used by more than one person due to the risk of misinterpreting glucose information.
- FreeStyle LibreLink and FreeStyle Libre or FreeStyle Libre 14 day Readers do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.
- Disable your phone's automatic operating system (OS) updates. Prior to updating your phone's OS or updating the App, you should check the Mobile Device and OS Compatibility Guide to determine if the FreeStyle LibreLink App is compatible with your OS and your phone. The OS Compatibility Guide is available in the Help Section of the App or on www.FreeStyleLibre.com. You should check the OS Compatibility Guide periodically to make sure that your OS and your phone continue to be compatible with the App.
- In the event that an App or OS update causes your previously compatible phone to become incompatible, you may be notified ahead of time via e-mail or through the App. Make sure that your LibreView account has your current e-mail address to receive important information.
- After an OS update, open your App and check your device settings to make sure it's working properly.

App Symbols

App icon

Sensor may be inaccurate. Check blood glucose with a test strip before making any treatment decisions

Direction your glucose is going. See [Understanding Your Glucose Readings](#) for more information.

Caution

Add/edit notes

Manually entered blood glucose result note

Add blood glucose result note

Food note

Insulin note

Food + insulin note

Exercise note

Time change

Main menu

Multiple/Custom notes

Share report

Additional information

Calendar

Sensor too cold

Sensor too hot

FreeStyle LibreLink Overview

IMPORTANT:

- Before you use your System, review all the product instructions and the Interactive Tutorial. You can access the Interactive Tutorial at www.FreeStyleLibre.com. The Quick Reference Guide and Interactive Tutorial give you quick access to important aspects and limitations of the System. The User's Manual includes all safety information and instructions for use. Refer to your smartphone instructions for use for how to use your smartphone.
- During the first 12 hours after insertion of a FreeStyle Libre 14 day Sensor, Sensor readings will be accompanied by the  symbol. Whenever  is displayed, a blood glucose test should be performed to confirm the Sensor reading prior to treatment.
- Talk to your health care professional about how you should use your Sensor glucose information to help manage your diabetes.
- FreeStyle LibreLink is available for download from the Google Play Store. When you're ready to start using FreeStyle LibreLink, you'll prepare and apply a Sensor to the back of your upper arm. You can then use the App to get glucose readings from the Sensor and store your glucose history and notes.

Note:


- The Sensor comes in a FreeStyle Libre or FreeStyle Libre 14 day Sensor kit. See [Sensor Kit](#).
- The App is not compatible with all smartphones. Before upgrading your smartphone or its operating system, check www.FreeStyleLibre.com.

Security Information

- You are responsible for properly securing and managing your smartphone. If you suspect an adverse cybersecurity event related to FreeStyle LibreLink, contact Customer Service.
- Make sure that your phone and Sensor kit are kept in a safe place, under your control. This is important to help prevent anyone from accessing or tampering with the System.
- FreeStyle LibreLink is not intended for use on a smartphone that has been altered or customized to remove, replace or circumvent the manufacturer's approved configuration or use restriction, or that otherwise violates the manufacturer's warranty.

Home Screen

The Home Screen gives you access to information about your glucose and the App. To return to the Home Screen from another screen, go to the Main Menu and tap **Home**.

Main Menu - Tap to access the Home Screen, Logbook, other history options, and Connected Apps. You can also access Settings, Help, and other information.

Glucose Graph - Graph of your stored Sensor glucose readings.

Scan Prompt - Indicates if the App is ready to scan a Sensor.

Glucose Information - Your Time In Range, information about your Last Reading, and average glucose for the last 24 hours.

Sensor Kit

The FreeStyle Libre or FreeStyle Libre 14 day Sensor Kit includes:


- Sensor Pack
- Sensor Applicator
- Product insert

When opening your kit, check that the contents are undamaged and that you have all parts listed. If any parts are missing or damaged, contact Customer Service. The Sensor (only visible after applied) is initially in two parts: one part is in the Sensor Pack and the other part is in the Sensor Applicator. Once prepared and applied to your body, the Sensor measures your glucose using a small, flexible tip that inserts just under the skin. After the start-up period, the Sensor can be worn for up to the wear duration specified in the Sensor Kit's product insert.

Sensor Pack. Used with the Sensor Applicator to prepare the Sensor for use.

Sensor Applicator. Applies the Sensor to your body.

App Setup

Before using the App for the first time, you must complete the setup.

1. Check that your smartphone is connected to a network (WiFi or cellular). You can then install FreeStyle LibreLink from the Google Play Store. Tap the App icon to open the App.
Note: You only need to be connected to a network for setup, using LibreView, and sharing with other apps. You do not need to be connected to scan a Sensor, add notes, or review your history in the App.
2. Swipe left to view some helpful tips or tap **GET STARTED NOW** at any point.
3. Confirm your country/region and tap **NEXT**.
4. You need a LibreView account to use the App. Follow onscreen instructions to review legal information and create a new account or login to your existing account.
5. Confirm your glucose unit of measure and tap **NEXT**.
6. Set your Target Glucose Range and tap **NEXT**. Work with your health care professional to determine your Target Glucose Range. Your Target Glucose Range is displayed on glucose graphs in the App and is used to calculate your Time In Range.
7. Select how you count carbohydrates (in grams or servings) and tap **NEXT**. The carbohydrate unit will be used in any food notes you enter in the App.
8. Select whether you want sound and vibration OR vibration only when you scan your Sensor. Tap **NEXT**.

9. The App now displays important information about:

- How to understand the My Glucose screen
- When to do a blood glucose test
- Where to apply the Sensor

Tap **NEXT** after reviewing each screen.

10. Accept required permissions.


11. Apply a new Sensor and then tap **NEXT**. Go to [Starting Your Sensor](#).

Note: If you need help applying your Sensor, tap **HOW TO APPLY A SENSOR** or go to [Applying Your Sensor](#).

Applying Your Sensor

CAUTION:

- The Sensor Pack and Sensor Applicator are packaged together as a set and have the same Sensor code. Check that the Sensor codes match before using your Sensor Pack and Sensor Applicator. Do not use Sensor Packs and Sensor Applicators with different Sensor codes together as this will result in incorrect glucose readings.


- Intense exercise may cause your Sensor to loosen due to sweat or movement of the Sensor. Remove and replace your Sensor if it starts to loosen and follow the instructions to select an appropriate application site.

IMPORTANT: Before using your Sensor Pack and Sensor Applicator, make sure you have an alcohol wipe (70% isopropyl alcohol) on hand to prepare the application site. This is not included in the Sensor Kit.

1. Apply Sensors only on the back of your upper arm. If placed in other areas, the Sensor

may not function properly and could give inaccurate readings. The application of the Sensor is not approved for other sites. Avoid areas with scars, moles, stretch marks or lumps. Select an area of skin that generally stays flat during your normal daily activities (no bending or folding). Choose a site that is at least 1 inch (2.5 cm) away from an insulin injection site. To prevent discomfort or skin irritation, you should select a different site other than the one most recently used.


2. Clean application site with an alcohol wipe and allow site to dry before proceeding. This helps the Sensor stay attached to your body.

Note: The area **MUST** be clean and dry, or the Sensor may not stick to the site.


3. Open the Sensor Pack by peeling the lid off completely. Unscrew the cap from the Sensor Applicator and set the cap aside.

CAUTION: Do NOT use if the Sensor Pack or Sensor Applicator seem to be damaged or already opened. Do NOT use if past expiration date.

4. Line up the dark mark on the Sensor Applicator with the dark mark on the Sensor Pack. On a hard surface, press down firmly on the Sensor Applicator until it comes to a stop.

5. Lift the Sensor Applicator out of the Sensor Pack.


6. The Sensor Applicator is prepared and ready to apply the Sensor.

CAUTION: The Sensor Applicator now contains a needle. Do NOT touch inside the Sensor Applicator or put it back into the Sensor Pack.


7. Place the Sensor Applicator over the prepared site and push down firmly to apply the Sensor to your body.

CAUTION: Do NOT push down on Sensor Applicator until placed over prepared site to prevent unintended results or injury.

8. Gently pull the Sensor Applicator away from your body. The Sensor should now be attached to your skin.

Note: Applying the Sensor may cause bruising or bleeding. If there is bleeding that does not stop, remove the Sensor, and contact your health care professional.

9. Make sure Sensor is secure after application. Put the cap back on the Sensor Applicator. Discard the used Sensor Applicator and Sensor Pack according to local regulations.

Note: Tap **Help** in the Main Menu to access an in-app tutorial on applying a Sensor.

Starting Your Sensor

IMPORTANT:

- Keep in mind that the ease of scanning a Sensor may vary between devices. Once you determine the location of the NFC (Near Field Communication) antenna on your smartphone, you should be able to reliably scan your Sensor by holding that area near your Sensor. You may need to adjust your scan distance based on what clothing you are wearing. In addition to proximity and orientation, other factors can affect NFC performance. For example, a bulky or metallic case can interfere with the NFC signal.
- The App requires that your smartphone be set to network time. This is the default time setting for most smartphones.
- When using the App, you should keep your smartphone well charged and be sure you have access to a blood glucose meter.
- Make sure your smartphone's sound is set so you can hear scan sounds if you turned those on.

-
1. Hold the back of your smartphone near the Sensor (this can be done over clothing) and position the NFC antenna over the Sensor. Do not move the smartphone until you hear the first tone and/or feel a vibration. This indicates your smartphone and Sensor have established an NFC connection.
 2. Continue holding your smartphone near the Sensor until you hear a second tone and/or feel a vibration. This completes the scan.

Note:

- If you need help, tap **HOW TO SCAN A SENSOR** to view an in-app tutorial. You can also access this later by going to the Main Menu and then tapping **Help**.

- If your Sensor is not successfully scanned, you may receive a scan error.

See [Troubleshooting](#) for additional error messages.


3. The Sensor can be used to check your glucose after the start-up period. While the Sensor is starting up you can navigate away from the App. You will see a notification when the Sensor is ready.

Note:

- If you like, you can use a Sensor with both the App and the Reader. To do this, you must start the Sensor with the Reader first and then scan with the App.
- Remember that FreeStyle LibreLink and Readers do not share data. For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.

Checking Your Glucose

1. Open the App and hold the back of your smartphone near the Sensor. If scan sounds are turned on, you will get two separate tones along with the vibrations when your Sensor has been scanned.
2. The My Glucose screen now displays your glucose reading. It includes your Current Glucose, a Glucose Trend Arrow indicating which way your glucose is going, and a graph of your current and stored glucose readings.

Message - Tap for more information.

Back - Tap to return to Home screen.

Current Glucose - Glucose value from your latest scan.

Check Blood Glucose - Tap for more information. When you see this symbol, do a blood glucose test before making treatment decisions.

Add Note - Tap to add notes to the glucose reading.

Glucose Trend Arrow - Direction your glucose is going.

Note Symbol - Tap to review notes you've entered.

Glucose Graph - Graph of your current and stored glucose readings.

Target Glucose Range - The graph shows your target glucose range.

Note:

- A Sensor can store up to 8 hours of glucose data, so scan it at least once every 8 hours to capture all of your available glucose data.
- The graph will scale to 500 mg/dL to accommodate glucose readings above 350 mg/dL.
- The ⌚ symbol may appear, indicating the smartphone's time was changed. Gaps in the graph may result or glucose readings may be hidden.
- Your current glucose value determines the background color on the My Glucose screen:

Orange - High glucose (above 240 mg/dL)

Yellow - Between the Target Glucose Range and high or low glucose level

Green - Within the Target Glucose Range

Red - Low glucose (below 70 mg/dL)

Understanding Your Glucose Readings

Glucose Trend Arrow


The Glucose Trend Arrow gives you an indication of the direction your glucose is going.

Glucose is rising quickly (more than 2 mg/dL per minute)

Glucose is rising (between 1 and 2 mg/dL per minute)

Glucose is changing slowly (less than 1 mg/dL per minute)

Glucose is falling (between 1 and 2 mg/dL per minute)


Glucose is falling quickly (more than 2 mg/dL per minute)

Note: The Glucose Trend Arrow may not always appear with your reading.


Messages

Below are messages you may see with your glucose readings.

LO | HI: If **LO** appears, your reading is lower than 40 mg/dL. If **HI** appears, your reading is higher than 500 mg/dL. You can touch the  symbol for more information. Check your blood glucose on your finger with a test strip. If you get a second **LO** or **HI** result, contact your health care professional **immediately**.


Low Glucose | High Glucose: If your glucose is higher than 240 mg/dL or lower than 70 mg/dL, you will see a message on the screen. You can touch the  symbol for more information and set a reminder to check your glucose.


Glucose Going Low | Glucose Going High: If your glucose is projected to be higher than 240 mg/dL or lower than 70 mg/dL within 15 minutes, you will see a message on the screen.

The background color corresponds to your current glucose value. You can touch the  symbol for more information and set a reminder to check your glucose.

You will see this symbol  if your glucose reading is less than 70 mg/dL, projected to be less than 70 mg/dL, rapidly changing, has no number or trend arrow, or it is in the first 12 hours of wearing a FreeStyle Libre 14 day Sensor. You can touch the symbol for more information. Check your blood glucose on your finger with a test strip before making treatment decisions.


Note: If you are not sure about a message or reading, contact your health care professional before you do anything.

Making Treatment Decisions

Work with your health care professional to put together a plan for managing your diabetes that includes when to use the System information for making treatment decisions.

WARNING: The System can replace blood glucose testing except in a few situations. These are the times when you need to do a blood glucose test before deciding what to do or what treatment decision to make as Sensor readings may not accurately reflect blood glucose levels:

Do a blood glucose test if you see the Check Blood Glucose  symbol. The  symbol means your Sensor glucose reading may not be accurate. For example, there may be times when you get a low glucose reading but you do not actually have low glucose.

Do a blood glucose test if you think your glucose readings are not correct or do not match how you feel. Do not ignore symptoms that may be due to low or high glucose.

Note: The  symbol will **NOT** appear in this situation.

Making Treatment Decisions – Getting Started

Before you start using the System for treatment decisions, make sure you have a good understanding of how the System works for your body. **Continue to use your blood glucose meter for treatment decisions until you are comfortable with the information you receive from your System.** This includes understanding that: Sensor performance can vary in between Sensors, within a Sensor wear period, and in different situations.

Getting familiar with the System could take days, weeks, or even months. The more you check readings from the System with a blood glucose meter, the better you will understand how the System works for you.

Work with your health care professional to put together a plan for managing your diabetes that includes when to use the System information for making treatment decisions.


Helpful Tips

- Confirm your Sensor glucose readings with a blood glucose meter until you understand:
 - Sensor accuracy may vary between Sensors.
 - Sensor accuracy may vary during a Sensor wear session.
 - Sensor accuracy may vary in different situations (meals, exercise, first day of use, etc.).
- Scan your Sensor often to see how carbs, medication, exercise, illness, or stress levels impact your Sensor glucose readings. The information you get can help you figure out why your glucose sometimes goes too high or too low, and how to prevent it from doing so in the future.
- Talk to your health care professional about how your insulin works. The more you understand about your insulin, including how long it takes to start working and how long it lasts in your body, the more likely you will be to make better treatment decisions.
- Making a treatment decision does not just mean taking insulin. Treatment decisions can also include things like taking fast-acting carbs, eating, or even doing nothing and scanning again later.
- Your health care professional can also help you to understand when doing nothing and scanning again later is the right treatment decision. For example, if your glucose is high and going up, your first instinct may be to take more insulin to lower your glucose, however, depending on when you last took insulin or your recent activity, the right treatment decision may be to do nothing and scan again later. Avoid "insulin stacking".

When not to use Sensor Glucose readings for treatment decisions

Glucose is Falling Quickly or Rising Quickly

Sensor glucose values, which are based on interstitial fluid glucose levels, can be different from blood glucose levels (fingersticks), particularly during times when your blood glucose is changing quickly. For example after eating, taking insulin, or exercising. When glucose levels are falling quickly, glucose readings from the Sensor may be higher than blood glucose levels. On the other hand, when glucose levels are rising quickly, glucose readings from the Sensor may be lower than blood glucose levels. If glucose is rising quickly or falling quickly, you will see the  symbol. Whenever you see the  symbol, do a blood glucose test and treat based on that result.

Low Glucose or Glucose Going Low message

The System lets you know about hypoglycemia or impending hypoglycemia with a Low Glucose or Glucose Going Low message. These messages may not accurately reflect blood glucose. When there is a Low Glucose or Glucose Going Low message, you will also see the 🩸🔍 symbol. Whenever you see the 🩸🔍 symbol, do a blood glucose test and treat based on that result.

No Glucose Trend Arrow

When there is no Glucose Trend Arrow, the System can't tell if your glucose is rising quickly or falling quickly and will display the 🩸🔍 symbol. Whenever you see the 🩸🔍 symbol, you should do a blood glucose test and treat based on that result.

No Current Glucose Number

When there is no Current Glucose number, such as when you receive an error message or a LO or HI result, you don't have enough information to make a treatment decision. When there is no Current Glucose you will see the 🩸🔍 symbol. Whenever you see the 🩸🔍 symbol, do a blood glucose test and treat based on that result.

During the first 12 hours of wearing a FreeStyle Libre 14 day Sensor

During the first 12 hours of wearing a FreeStyle Libre 14 day Sensor, you should not use glucose readings from the Sensor to make treatment decisions. You will see the  symbol displayed during this time, and whenever you see the , do a blood glucose test and treat based on that result.

Think Your Readings are Incorrect?

Don't trust Sensor glucose readings that you think may be incorrect or that don't match what you would expect based on your recent activity. For example, if you ate dinner but forgot to take insulin before eating, you would expect your glucose to be high. If your glucose reading is low, then it doesn't match your recent activity, so don't use it to make treatment decisions. Don't make treatment decisions if you think your Sensor glucose readings are incorrect. Do a blood glucose test and treat based on that result.

You Have Low or High Blood Glucose Symptoms

Don't ignore symptoms that may be due to low or high blood glucose. Do a blood glucose test and treat based on that result.

Symptoms Don't Match Readings

There may be times when your symptoms don't match your Sensor glucose readings. For example, you are feeling shaky, sweaty, and dizzy – symptoms you generally get when you have low glucose, but your glucose reading is within your target range. When symptoms don't match readings, do a blood glucose test and treat based on that result. Don't ignore symptoms that may be due to low or high blood glucose.

If you're the caregiver, pay attention to times when the symptoms of the one you're caring for don't match their Sensor glucose readings. When symptoms don't match readings, do a blood glucose test and treat based on that result.

Note: The  symbol will **NOT** appear in these situations.


When to do Nothing and Scan Again Later

Your health care professional can help you understand when doing nothing and scanning again later is the right treatment decision. For example, if your glucose is high and going up, your first instinct may be to take more insulin to lower your glucose, however depending on when you last took insulin or your recent activity, the right treatment decision may be to do nothing and scan again later.

Don't take a correction dose within 2 hours of your meal dose. This may result in "insulin stacking" and low glucose.

Making Treatment Decisions - Advanced

After you scan your Sensor, use all of the information on the screen when deciding what to do or what treatment decision to make.

This table provides some information on how you can factor the Glucose Trend Arrow into your treatment decisions. Remember that you should never make a treatment decision based on the Glucose Trend Arrow alone.

Treatment Decision Considerations for: No Arrow or No Number

You will see the symbol. Do not treat based on Sensor glucose reading. Do a blood

glucose test.

Treatment Decision Considerations for Glucose Trend Arrow: ↑

You will see the  symbol. Do not treat based on Sensor glucose reading. Do a blood glucose test.

Treatment Decision Considerations for Glucose Trend Arrow: ↗

Low Glucose (< 70 mg/dL): You will see the  symbol. Do not treat based on Sensor glucose reading. Do a blood glucose test.

Glucose in Target Range: If you are about to eat, take insulin to cover your meal. Consider taking a little more since glucose is rising.

If you have taken insulin recently, do nothing and scan again later.

Avoid "insulin stacking".

High Glucose (> 240 mg/dL): If you are about to eat, take insulin to cover your meal. Consider taking a little more since glucose is high and rising.

If this is between meals, consider taking an insulin correction dose, unless you have taken insulin recently. If you have taken insulin recently, do nothing and scan again later.

Avoid "insulin stacking".

Treatment Decision Considerations for Glucose Trend Arrow: →

Low Glucose (< 70 mg/dL): You will see the  symbol. Do not treat based on Sensor glucose reading. Do a blood glucose test.

Glucose in Target Range: If you are about to eat, take insulin to cover your meal.

If this is between meals, do nothing and scan again later.

High Glucose (> 240 mg/dL): If you are about to eat, take insulin to cover your meal. Consider taking a little more since glucose is high.

If this is between meals, consider taking an insulin correction dose, unless you have taken insulin recently. If you have taken insulin recently, do nothing and scan again later.

Avoid "insulin stacking".

Treatment Decision Considerations for Glucose Trend Arrow: ↘

Low Glucose (< 70 mg/dL): You will see the  symbol. Do not treat based on Sensor glucose reading. Do a blood glucose test.

Glucose in Target Range: If you are about to eat, take insulin to cover your meal. Consider taking a little less since glucose is falling.

If this is between meals, consider eating a snack or fast-acting carbohydrates to stay within target and scan again later.

High Glucose (> 240 mg/dL): If you are about to eat, take insulin to cover your meal. Consider taking a little less since glucose is falling.

If this is between meals, consider doing nothing and scan again later.

Avoid "insulin stacking".


Treatment Decision Considerations for Glucose Trend Arrow: ↓

You will see the  symbol. Do not treat based on Sensor glucose reading. Do a blood glucose test.

Example Scenarios

Here are some example scenarios to help you understand how to use the information on your screen. Always use all of the information on the screen before deciding what to do or treatment decision to make. If you are not sure about what to do, consult your health care professional.

What you see - *When you wake-up*


What it means

When you wake up, your current glucose is 65 mg/dL and the trend arrow shows it is changing slowly →. There is also a  message at the top of the screen and the  symbol.

Anytime you see the  symbol, you should do a blood glucose test before deciding what to do.

What you see - *Before breakfast*


What it means

Before breakfast, your current glucose is 115 mg/dL. The graph shows that your glucose is going up and so does the trend arrow [↗].

Consider what might be causing your glucose to go up and what you might do to prevent a high glucose. For example:

- How much insulin should you take before your meal?
- Since you see [↗], should you consider taking a little more insulin?

What you see - *After breakfast*


What it means

After breakfast, your current glucose is 108 mg/dL. The trend arrow shows it is going down quickly ↓. There is also a  GLUCOSE GOING LOW message at the top of the screen and the  symbol.

Anytime you see the  symbol, you should do a blood glucose test before deciding what to do.


What you see - *Before lunch*

What it means

When you checked your glucose before lunch, it was 90 mg/dL and rising. Before eating lunch, you took enough insulin to cover the meal and a little more since your trend arrow was $\nearrow$.

What you see - *After lunch*

What it means


90 minutes later, your current glucose is 225 mg/dL. The graph shows that your glucose is still going up, and so does the trend arrow ↗.

Don't take a correction dose within 2 hours of your meal dose. This may result in "insulin stacking" and low glucose.

Consider what might be causing your glucose to go up and what you might do to prevent a high glucose. For example:

- Has the insulin you took for your meal reached its full effect?
- Scan your Sensor again later.

What you see - *After exercising*


What it means

After exercising, you are feeling shaky, sweaty, and dizzy – symptoms you generally get when you have low glucose. But, your current glucose is 204 mg/dL.

Anytime you get a reading that doesn't match how you feel, do a blood glucose test.

Note: The Check Blood Glucose  symbol will **NOT** appear in this situation.

What you see - *Before dinner*


What it means

Before dinner, your current glucose is 134 mg/dL. The graph shows that your glucose is going down and so does the trend arrow $\searrow$.

Consider what might be causing your glucose to go down and what you might do to prevent a low glucose. For example:

- How much insulin should you take to cover your meal?
- Since you see $\searrow$, should you think about taking a little less insulin?

What you see - *After dinner*

What it means

After dinner, your current glucose is 215 mg/dL but there is no trend arrow. There is also the  symbol on the screen.

Anytime you see the  symbol, you should do a blood glucose test before deciding what to do.

Other considerations

Deciding how much rapid-acting insulin to take for different meals and situations can be difficult. Work with your health care professional to discuss different situations and what might work best for you. Here are some questions to consider:

Meal dosing

- What do you do if your before meal glucose is high?
- What do you do if your before meal glucose is low?
- How much time do you wait to eat after taking your meal insulin?
- Do you adjust the amount of meal insulin based on the number of carbs or how much you are planning to eat?

- Do you adjust your meal insulin dose for high fat foods such as pizza?
- Do you know how to adjust your insulin doses when drinking alcoholic beverages?

High glucose corrections

- Do you take extra insulin if your glucose is high?
- How do you decide how much insulin to take for a high glucose?
- How long do you wait between insulin doses to avoid insulin stacking?

Bedtime

- How often do you check your glucose before bed?
- What do you consider a safe bedtime glucose?
- What do you do if your bedtime glucose is high?
- What do you do if your bedtime glucose is low?
- When should you eat a bedtime snack?
- What do you do if your before meal glucose is high?
- What do you do if your before meal glucose is low?

Other factors

- How do you adjust your insulin dose based on the Glucose Trend Arrow?
- How do you adjust your insulin dose for different types of exercise or activities?
- How do you adjust your insulin doses for stress?
- How do you adjust your insulin doses for illness?

Adding Notes

Notes can be saved with your glucose readings to help you track food, insulin, and exercise. You can also add your own comment.

1. Tap the  symbol on your glucose reading screen.
2. Select the checkbox next to the notes you would like to add. After you check the box, you can add more specific information to your note.
 - Food notes: Enter meal type and grams or serving information

- Insulin notes: Enter the number of units taken
- Exercise notes: Enter intensity and duration

Note: You cannot add a blood glucose result note to a glucose reading. See [Reviewing Your History](#) for how to enter blood glucose results in the App's Logbook.

3. Tap **DONE** to save your note.

Notes you add are shown on your glucose graph and in your Logbook as symbols. You can review a note by tapping its symbol on your glucose graph or by going to the Logbook. See [Reviewing Your History](#) for more information about the Logbook. To edit a note from the glucose graph, tap the symbol and then tap on the information you would like to change. Tap **DONE** when you are finished.

Food

Insulin

Exercise

Food + insulin

Multiple/Custom notes – indicates different types of notes entered together or notes entered within a short period of time. A numbered badge next to the symbol indicates the number of notes.

Reviewing Your History

Reviewing and understanding your glucose history can be an important tool for improving your glucose control. The App stores about 90 days of information and has several ways to review your past glucose readings and notes. From the Main Menu, tap **Logbook** to view the Logbook or tap on one of the other history options under **Reports**.

IMPORTANT:

- Work with your health care professional to understand your glucose history.
- FreeStyle LibreLink does not share data with the FreeStyle Libre Reader.

- For complete information on a device, be sure to scan your Sensor every 8 hours with that device; otherwise, your reports will not include all your data.

Logbook

The **Logbook** contains entries for each time you scanned your Sensor as well as notes you added. If you would like to view a different day, tap the  symbol or use the arrows. To add a note to a Logbook entry, tap on the entry and then tap the  symbol. Select your note information and tap **DONE**.

To add a note that is independent of a Logbook entry, tap the  symbol on the main Logbook screen. Tap the  symbol if you want to add a note on a different date.

The logbook also allows you to record a blood glucose test you performed. To do this, tap the  symbol and enter your result. You can enter a result between 20-500 mg/dL but cannot enter LO or HI.

Other History Options

Daily Patterns: A graph showing the pattern and variability of your Sensor glucose readings over a typical day. The thick black line shows the median (midpoint) of your glucose readings. The light blue shading represents the 10th - 90th percentile range of your glucose readings. Dark blue shading represents the 25th - 75th percentile range.

Note: Daily Patterns needs at least 5 days of glucose data.

Time In Range: A graph showing the percentage of time your Sensor glucose readings were above, below, or within your Target Glucose Range.

Low Glucose Events: Information about the number of low glucose events measured by your Sensor. A low glucose event is recorded when your Sensor glucose reading is lower than 70 mg/dL for 15 minutes or longer. The total number of events is displayed below the graph. The bar graph displays the low glucose events in different periods of the day.

Average Glucose: Information about the average of your Sensor glucose readings. The overall average for the selected time period is displayed below the graph. The average is also shown for different periods of the day. Readings above or below your Target Glucose Range are yellow, orange, or red. Readings in range are green.

Daily Graph: A graph of your Sensor glucose readings by day. The graph shows your Target Glucose Range and symbols for notes you have entered.

- The graph will scale to 500 mg/dL to accommodate glucose readings above 350 mg/dL.
- You might see gaps in the graph during times when you have not scanned at least once in 8 hours.
- The ⌚ symbol may appear indicating a time change. Gaps in the graph may result or glucose readings may be hidden.

Sensor Usage: Information about how often you scanned your Sensor or viewed your Sensor glucose readings in the App and how much information has been captured from your Sensor.

Note:

- Tap the ➦ symbol on any report to share a screenshot of the report.
- Tap the ⓘ symbol to view a description of the report.
- From any report screen, swipe left or right to view the next or previous report.
- On all reports except the **Daily Graph**, you can select to show information about your last 7, 14, 30, or 90 days.

Removing Your Sensor

1. Pull up the edge of the adhesive that keeps your Sensor attached to your skin. Slowly peel away from your skin in one motion.

Note: Any remaining adhesive residue on the skin can be removed with warm soapy water or isopropyl alcohol.

2. Discard the used Sensor according to local regulations. See [Maintenance and Disposal](#). When you are ready to apply a new Sensor, follow the instructions in [Applying Your Sensor](#) and [Starting Your Sensor](#).

Replacing Your Sensor

Your Sensor automatically stops working after the wear duration specified in your Sensor Kit's product insert and must be replaced. You should also replace your Sensor if you notice any irritation or discomfort at the application site or if the App reports a problem with the Sensor currently in use. Taking action early can keep small problems from turning into larger ones.

CAUTION: If the Sensor is becoming loose or if the Sensor tip is coming out of your skin, you may get no readings or unreliable readings, which may not match how you feel. Check to make sure your Sensor has not come loose. If it has come loose, remove it and apply a new one.

Setting Reminders

You can create single or repeating reminders to help you remember things like checking your glucose or taking insulin. There is one default reminder to help you remember to check your glucose. This reminder can be changed or disabled but cannot be deleted.

Note:

- If you want to receive a sound/vibration with your reminder, ensure that sound/vibration on your smartphone is turned on, sound is set at a level you can hear, and your smartphone's Do Not Disturb feature is turned off. If Do Not Disturb is on, you will only see your reminder on the screen.
 - Allow Alarms and Reminders permission for the App. This permission is required for Android 12 and above.
1. To add a new reminder, go to the Main Menu and tap **Reminders**. Tap **ADD REMINDER**.
 2. Name your reminder.
 3. Tap the time fields to set the time for the reminder.
 4. Tap **DONE**. You will now see your reminder on the list along with the time you will receive it.

Note:

- To delete a reminder, swipe the reminder and tap the  symbol. The reminder to check your glucose cannot be deleted.

- Your reminders will be received as notifications that you can swipe or tap to dismiss.

Settings and Other Options in the Main Menu

Settings

App Settings:

Unit of Measurement - View the glucose unit of measure used in the App.

Report Settings - Set the target range you want to be displayed on the App glucose graph. It is also used to calculate your Time In Range. Tap **SAVE** when you are done.

Carbohydrate Units - Choose grams or servings for food notes that you enter. Tap **SAVE** when you are done.

Scan Sounds - Select whether you would like to hear a sound in addition to a vibration when you scan the Sensor. Remember that Scan Sounds inherit the volume settings on your smartphone. If your smartphone volume is turned off, you will not hear a scan sound. Tap **SAVE** when you are done.

Text to Speech - Turn on to have the glucose reading read aloud when you scan the Sensor. You will hear only your current glucose value and trend arrow direction along with whether you need to do a blood glucose test before making treatment decisions. Additional information, such as the glucose graph and any message, are available on your glucose reading screen. Always review your glucose reading to get complete information. Remember that this feature inherits the volume settings on your phone. If your phone volume is turned off, you will not hear the glucose reading read aloud. Tap **SAVE** when you are done.

Account Settings:

Account Details - View/change your LibreView account information.

Account Password - Change your LibreView account password.

Connected Apps

The **Connected Apps** option in the Main menu opens a web browser within the App. It lists different apps you can connect with to share your data. To connect your data with apps listed in the Connected Apps option, select them from the list of apps, and follow the onscreen instructions.

Help

View in-app tutorials, access the App's product labeling, and review the App's legal information. You can also view the Event Log, which is a list of events recorded by the App. This may be used by Customer Service to help troubleshoot.

About

View App software version and other information.

Activities

Bathing, Showering, and Swimming: Your Sensor is water-resistant and can be worn while bathing, showering, or swimming. Do NOT take your Sensor deeper than 3 ft (1 meter) or immerse it longer than 30 minutes in water.

Sleeping: Your Sensor should not interfere with your sleep. It is recommended that you scan your Sensor before going to sleep and when you wake up because your Sensor holds only 8 hours of data at a time. For example, if you sleep for 9 hours without scanning your Sensor, 1 hour of data will not be collected and a gap will appear on your glucose graph. If you have reminders set to go off while you are sleeping, place your smartphone nearby.

Traveling by Air: You can safely use your System at all times while on an aircraft:

- You can continue to get Sensor glucose readings after you put your smartphone in airplane mode and turn NFC back on. Please note you must comply with requests from the flight crew to not scan your Sensor.
- Some airport full-body scanners include X-ray or millimeter radio-wave, which you cannot expose your Sensor to. The effect of these scanners has not been evaluated and the exposure may damage the Sensor or cause inaccurate results. To avoid removing your Sensor, you may request another type of screening. If you do choose to go through a full-body scanner, you must remove your Sensor.
- The Sensor can be exposed to common electrostatic (ESD) and electromagnetic interference (EMI), including airport metal detectors.

Note: Changing the time affects the graphs, statistics, and settings programmed by time of day since these are dependent on the time setting being accurate. The ⌚ symbol may appear on your glucose graph indicating a time change. Gaps in the graph may result or glucose readings may be hidden.

Maintenance and Disposal

Maintenance: The Sensor has no serviceable parts.

Disposal: The product contains electronic equipment, batteries, sharps and materials that may contact bodily fluids during use. Dispose of product in accordance with all applicable local regulations. Contact Customer Service for further information on the appropriate disposal of components.

Troubleshooting

This section lists problems that you may experience, the possible cause(s), and recommended actions. If there is an error, a message will appear on the screen with directions to resolve the error.

IMPORTANT: If you are having issues with the App, please keep in mind that uninstalling the App and/or clearing data will cause you to lose all historical data and end the Sensor currently in use. Please call Customer Service if you have any questions.

Problems at the Sensor Application Site

Problem: **The Sensor is not sticking to your skin.**

What it may mean: The site is not free of dirt, oil, hair, or sweat.

What to do: 1. Remove the Sensor. 2. Consider shaving and/or cleaning the site with soap and water. 3. Follow the instructions in [Applying Your Sensor](#) and [Starting Your Sensor](#).

Problem: **Skin irritation at the Sensor application site.**

What it may mean: Seams or other constrictive clothing or accessories causing friction at the site **OR** you may be sensitive to the adhesive material.

What to do: Ensure that nothing rubs on the site. If the irritation is where the adhesive touches skin, contact your health care professional to identify the best solution.

Problems Starting Your Sensor or Receiving Sensor Readings

Display: **Sensor Starting Up**

What it may mean: Sensor is not ready to read glucose.

What to do: Wait until the start-up period has completed.

Display: **Sensor Ended**

What it may mean: The Sensor life has ended.

What to do: Apply and start a new Sensor.

Display: New Sensor Found

What it may mean: You scanned a new Sensor before your previous Sensor ended.

What to do: Your smartphone can only be used with one Sensor at a time. If you start a new Sensor, you will no longer be able to scan your old Sensor. If you would like to begin using the new Sensor, select **Yes** and scan again.

Display: Sensor Error

What it may mean: The Sensor is unable to provide a glucose reading.

What to do: Scan again after the duration specified in the message.

Note: If you are using a FreeStyle Libre 14 day Sensor and receive this error during your first 12 hours of wearing a Sensor, it may mean that your body is still adjusting to the Sensor. Use a blood glucose meter to check your glucose while you wait. You do not need to remove your Sensor.

Display: Glucose Reading Is Unavailable

What it may mean: The Sensor is unable to provide a glucose reading.

What to do: Scan again in 10 minutes.

Display: Sensor Too Hot

What it may mean: Your Sensor is too hot to provide a glucose reading.

What to do: Move to a location where the temperature is appropriate and scan again in a few minutes.

Display: Sensor Too Cold

What it may mean: Your Sensor is too cold to provide a glucose reading.

What to do: Move to a location where the temperature is appropriate and scan again in a few minutes.

Display: Check Sensor

What it may mean: The Sensor tip may not be under your skin.

What to do: Try to start your Sensor again. If you see **Check Sensor** again on the screen, your Sensor was not applied properly. Apply and start a new Sensor.

Display: Replace Sensor

What it may mean: The App has detected a problem with your Sensor.

What to do: Apply and start a new Sensor.

Display: Unexpected Application Error

What it may mean: The App has detected an unexpected error.

What to do: Shut down the App completely and restart.

Display: Incompatible Sensor

What it may mean: The Sensor cannot be used with the App. Check that you have installed the app that is compatible with your Sensor. You may need to download a different app if your Sensor is not compatible.

What to do: Tap Learn More to find out what Sensors can be used. If you still have questions, call Customer Service.

Display: Scan Error

What it may mean: The smartphone was unable to scan the Sensor **OR** another NFC app is competing for your smartphone's NFC.

What to do: Try scanning the Sensor again. Make sure to hold your smartphone still once you get the first tone and/or vibration. Wait until you get the second tone and/or vibration to move it away from the Sensor. Make sure you are not touching any buttons on the smartphone or the screen **OR** open the App to scan your Sensor because another app that uses NFC was detected.

Display: App Permission Required

What it may mean: A required App permission is turned off.

What to do: Follow the instructions on the screen to turn the permission on.

Customer Service

Customer Service is available at 1-855-632-8658 7 Days a Week from 8AM to 8PM Eastern Time; excluding holidays. A printed copy of this User's Manual is available upon request. The latest version of the User's Manual is available at

www.FreeStyleLibre.us/support/overview.html

Sensor Specifications

Sensor glucose assay method: Amperometric electrochemical sensor

Sensor glucose reading range: 40 to 500 mg/dL

Sensor size: 5mm height and 35mm diameter

Sensor weight: 5 grams

Sensor power source: One silver oxide battery

Sensor life: Up to the wear duration specified in your Sensor Kit's product insert

Sensor memory: 8 hours (glucose readings stored every 15 minutes)

Operating temperature: 50°F to 113°F

Sensor Applicator and Sensor Pack storage temperature: 39°F to 77°F

Operating and storage relative humidity: 10-90%, non-condensing

Sensor water resistance and ingress protection: IP27: Can withstand immersion into 3 ft (1 meter) of water for up to 30 minutes. Protected against insertion of objects > 12mm diameter.

Operating and storage altitude: -1250 ft (-381 meters) to 10,000 ft (3,048 meters)

Labeling Symbols and Definitions

Consult instructions for use

Temperature limit

Manufacturer

Batch code

Unique Device Identifier

Type BF applied part

CODE Sensor code

Do not re-use

MR unsafe

FCC Declaration of Conformity mark

Use-by date

Catalog number

Serial number

Caution

Sterilized by irradiation

Humidity limitation

Do not use if package is damaged.

CAUTION: Federal law restricts this device to sale by or on the order of a physician.

This product contains electronic equipment, batteries, sharps and materials that may contact bodily fluids during use.

Dispose of product in accordance with all applicable local regulations.

Performance Characteristics - FreeStyle Libre Sensor

Clinical Study Overview

Performance of the FreeStyle Libre Flash Glucose Monitoring System (the System) was evaluated in a clinical study. The study was conducted at 4 centers with a total of 48 subjects with diabetes (95.8% Type 1, 4.2% Type 2). All subjects were aged eighteen and older. Subjects in the study required insulin to manage their diabetes. Each subject wore up to two

FreeStyle Libre Sensors on the back of the upper arm. During the study, subjects tested their blood glucose using fingerstick capillary samples at least eight times during each day of the study. Subjects used the blood glucose meter built into the FreeStyle Libre Reader. Additionally, subjects had their venous blood glucose analyzed up to 128 times over four separate visits to the clinical center. Venous blood was analyzed using the Yellow Springs Instrument Life Sciences 2300 STAT Plus™ Glucose & Lactate Analyzer (YSI). YSI is a laboratory glucose and lactate analyzer of whole blood and plasma and is a widely recognized standard in laboratory analysis of blood glucose. Glucose readings obtained from the System were compared to glucose readings obtained from the YSI to evaluate the performance of the System. Three lots of Sensors were evaluated in the study.

Agreement with YSI Levels

Agreement between FreeStyle Libre Glucose Measurement (CGM) and venous blood was characterized by using paired CGM and Yellow Springs Instrument measurements (YSI). The accuracy of CGM versus YSI reference was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for reference values 80 mg/dL and above, and 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL when glucose levels are assigned using the YSI values. Overall 91.1% of results were within ± 20 mg/dL / 20% of YSI reference.

Agreement with CGM Glucose Levels

Agreement between CGM and venous blood was characterized by using paired CGM and Yellow Springs Instrument measurements (YSI). The accuracy of CGM versus YSI reference was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for reference values 80 mg/dL and above, and 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL. The results are presented in [Table 1](#) for YSI reference. Overall 91.0% of results were within ± 20 mg/dL / 20% of YSI reference.

Table 1. Number and Percent of Results within YSI Reference

CGM Glucose Level (mg/dL)	Number of CGM-Reference Pairs	Within $\pm 15\%$ / ± 15 mg/dL	Within $\pm 20\%$ / ± 20 mg/dL	Within $\pm 30\%$ / ± 30 mg/dL	Within $\pm 40\%$ / ± 40 mg/dL	Outside $\pm 40\%$ / ± 40 mg/dL
Overall	5772	82.1	91.0	97.8	99.3	0.7
40-50	38	44.7	57.9	81.6	94.7	5.3
51-80	461	72.2	81.1	92.0	97.6	2.4
81-180	3236	82.9	91.2	97.9	99.3	0.7
181-300	1799	84.9	93.6	99.2	99.7	0.3
301-400	226	77.0	95.1	99.6	99.6	0.4
401-500	12	58.3	75.0	100.0	100.0	0.0

Agreement on Day 1 against YSI Reference

The accuracy of CGM versus YSI reference on the first day of Sensor wear was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for reference values 80 mg/dL and above, and within 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL by hourly intervals. The results are presented in [Table 2](#).

Table 2. Number and Percent of Results within YSI Reference

Time Interval (hours)	Number of CGM-Reference Pairs	Within $\pm 15\%$ / $\pm 15\text{mg/dL}$	Within $\pm 20\%$ / $\pm 20\text{mg/dL}$	Within $\pm 30\%$ / $\pm 30\text{mg/dL}$	Within $\pm 40\%$ / $\pm 40\text{mg/dL}$	Outside $\pm 40\%$ / $\pm 40\text{mg/dL}$
(0-2)	81	69.1	87.7	100.0	100.0	0.0
(2-4)	318	73.9	84.6	97.2	99.7	0.3
(4-6)	374	76.7	88.0	97.3	99.7	0.3
(6-8)	369	79.9	90.8	99.2	100.0	0.0

Overall Accuracy against YSI reference

Accuracy was measured by comparing the absolute relative difference between the System and reference YSI glucose values. The absolute relative difference measures the level of disagreement between the System and the reference value, but does not tell you whether the System glucose value was, on average, higher or lower than the reference glucose value. The Mean Absolute Relative Difference gives an indication of the average percent disagreement between the CGM and the reference. [Table 3](#) shows the overall absolute difference measure. Overall the Mean Absolute Relative Difference was 9.7% for the comparison with YSI reference. The Median Absolute Relative Difference shows that half of the time the System was within 7.7% of the YSI reference.

Table 3. Difference Measures with YSI Reference

Number of CGM-Reference Pairs	Median Absolute Relative Difference	Mean Absolute Relative Difference
5772	7.7%	9.7%

Agreement with BG Levels

Agreement between the System and capillary blood glucose values (BG) as measured by the Reader's built-in meter was characterized by using paired System CGM and BG value. The accuracy of CGM versus BG value was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for BG values 80 mg/dL and above, and within 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL. The results are presented in [Table 4](#) for BG values. Overall 84.3% of results were within $\pm 20\text{ mg/dL}$ / 20% of BG values.

Table 4. Number and Percent of Results within BG Values*

--	--	--	--	--	--	--

CGM Glucose Level (mg/dL)	Number of CGM-Reference Pairs	Within $\pm 15\%$ / $\pm 15\text{mg/dL}$	Within $\pm 20\%$ / $\pm 20\text{mg/dL}$	Within $\pm 30\%$ / $\pm 30\text{mg/dL}$	Within $\pm 40\%$ / $\pm 40\text{mg/dL}$	Outside $\pm 40\%$ / $\pm 40\text{mg/dL}$
Overall	3680	72.8	84.3	95.0	98.0	2.0
40-50	23	47.8	73.9	87.0	95.7	4.3
51-80	288	65.6	77.1	89.6	97.6	2.4
81-180	1722	71.7	82.9	94.7	97.4	2.6
181-300	1193	75.0	87.3	96.2	98.4	1.6
301-400	362	76.8	87.8	97.2	99.2	0.8
401-500	92	78.3	85.9	95.7	98.9	1.1

* Comparison to BG was performed using the FreeStyle Libre Reader's built-in blood glucose meter. Different performance may be expected when compared to other models of blood glucose meters.

Overall Accuracy against BG values

Accuracy was measured by comparing the absolute relative difference between the System and BG values. The absolute relative difference measures the level of disagreement between the System and the BG value, but does not tell you whether the System glucose value was, on average, higher or lower than the BG glucose value. The Mean Absolute Relative Difference gives an indication of the average percent disagreement between the CGM and the BG value. [Table 5](#) shows the overall absolute difference measure. Overall the Mean Absolute Relative Difference was 12.1% for the comparison with BG value. The Median Absolute Relative Difference shows that half of the time the System was within 9.4% of the BG value.

Table 5. Difference Measures with BG Value*

Number of CGM-Reference Pairs	Median Absolute Relative Difference	Mean Absolute Relative Difference
3680	9.4%	12.1%

* Comparison to BG was performed using the FreeStyle Libre Reader's built-in blood glucose meter. Different performance may be expected when compared to other models of blood glucose meters.

Concurrence of System and Reference (CGM vs. YSI)

The percentage of concurring glucose values (CGM vs. YSI) in each glucose reference range is presented for each CGM range in [Table 6](#). For example, in the clinical study, when the System glucose results were within the 81 to 120 mg/dL range, actual blood glucose values were less than 40 mg/dL 0% of the time, between 40 and 60 mg/dL 0.2% of the time, between 61 and 80 mg/dL 5.6% of the time, between 81 and 120 mg/dL 75.9% of the time, between 121 and 160 mg/dL 17.6% of the time, between 161 and 200 mg/dL 0.6% of the time, between 201 and 250 mg/dL 0.1% of the time and above 250 mg/dL 0% of the time.

Table 6. Concurrence Analysis by Glucose Level

CGM (mg/dL)	YSI Glucose Level (mg/dL)												N
	<40*	40-60	61-80	81-120	121-160	161-200	201-250	251-300	301-350	351-400	401-500	>500*	
<40	0.0	19.0	61.9	19.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	21
40-60	0.7	25.2	58.7	15.4	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	143
61-80	0.0	7.3	45.8	46.3	0.6	0.0	0.0	0.0	0.0	0.0	0.0	0.0	356
81-120	0.0	0.2	5.6	75.9	17.6	0.6	0.1	0.0	0.0	0.0	0.0	0.0	1222
121-160	0.0	0.0	0.1	13.2	72.0	14.3	0.4	0.0	0.0	0.0	0.0	0.0	1435
161-200	0.0	0.0	0.0	0.3	21.9	67.5	10.2	0.1	0.0	0.0	0.0	0.0	1087
201-250	0.0	0.0	0.0	0.0	0.9	28.8	64.5	5.7	0.0	0.0	0.0	0.0	905
251-300	0.0	0.0	0.0	0.0	0.0	0.3	41.2	53.4	5.2	0.0	0.0	0.0	386
301-350	0.0	0.0	0.0	0.0	0.0	0.0	1.2	55.3	40.6	2.9	0.0	0.0	170
351-400	0.0	0.0	0.0	0.0	0.0	0.0	0.0	3.6	66.1	30.4	0.0	0.0	56
401-500	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	25.0	33.3	41.7	0.0	12
>500	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0

* Levels out of system dynamic range.

Concurrence Analysis by Glucose Trend Arrow

Trend Arrow accuracy of the system, as assessed by concurrence analysis, is presented in [Table 7](#). For example, in the clinical study, when the trend arrow indicated that glucose was changing slowly (-1 to 1 mg/dL/min (→)), actual glucose levels in the body were falling quickly (↓) 0.3% of the time, falling (↘) 3.7% of the time, changing slowly (→) 83.0% of the time, rising (↗) 3.9% of the time, and rising quickly (↑) 0.5% of the time.

Table 7. Concurrence Analysis by Glucose Trend Arrow

CGM (mg/dL/min)	YSI (mg/dL/min)							N
	<-2	[-2,-1]	[-1,0]	[0,1]	[1,2]	>2	NA*	
<-2 (↓)	26.3	45.5	10.3	1.3	0.0	0.0	16.7	156
-2 to -1 (↘)	4.3	27.0	54.6	3.8	0.6	0.0	9.7	652
-1 to 1 (→)	0.3	3.7	49.4	33.6	3.9	0.5	8.6	4175
1 to 2 (↗)	0.0	0.6	8.8	38.6	33.3	9.2	9.4	477
>2 (↑)	0.0	0.0	2.8	14.6	34.9	40.6	7.1	212
NA†	2.5	9.1	33.1	20.7	14.0	9.1	11.6	121

* Glucose rate of change not available due to the time difference between glucose readings exceeding 30 minutes.

† Glucose Trend Arrow not available.

Agreement with 'LO' and 'HI' CGM Reading against YSI Reference

The System reports glucose concentrations between 40 and 500 mg/dL. When the System

determines that glucose level is below 40 mg/dL, it will report as 'LO'. When the System determines that glucose level is above 500 mg/dL, it will report as 'HI'. No measurements were obtained above 500 mg/dL in the clinical study. [Table 8](#) displays the concurrence between the CGM and YSI reference glucose when CGM reads 'LO'. For example, in the clinical study, when CGM reading was 'LO' YSI glucose values were less than 40 mg/dL 0.0% of the time, above 40 mg/dL 100.0% of the time, above 50 mg/dL 95.2% of the time, above 60 mg/dL 80.9% of the time, above 70 mg/dL 42.8% of the time, and above 80 mg/dL 19.0% of the time.

Table 8. Concurrence Analysis with 'LO' CGM Reading

	YSI (mg/dL)						Total
	<40	>40	>50	>60	>70	>80	
% of CGM points in YSI range	0.0	100.0	95.2	80.9	42.8	19.0	
Number of CGM points in YSI range	0	21	20	17	9	4	21

Accuracy by Day of Wear

After the 12 hour start-up period, the Sensor can be worn for up to 10 days. To show Sensor performance over time, the absolute relative difference between the System and reference YSI glucose values over the wear is presented in [Table 9](#).

Table 9. Difference Measures by Day (YSI Reference)

Day	Number of CGM-Reference Pairs	Median Absolute Relative Difference (%)	Mean Absolute Relative Difference (%)
1	1497	8.7	10.7
4	1470	7.4	9.6
7	1394	7.4	9.1
10	1411	7.5	9.3

The accuracy of CGM versus YSI reference and BG reference was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for reference values 80 mg/dL and above, and 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL. The results for CGM vs. YSI reference are presented in [Table 10](#).

Table 10. Number and Percent of Results within YSI Reference

Day	Number of CGM-Reference Pairs	Within $\pm 15\%$ / $\pm 15\text{mg/dL}$	Within $\pm 20\%$ / $\pm 20\text{mg/dL}$	Within $\pm 30\%$ / $\pm 30\text{mg/dL}$	Within $\pm 40\%$ / $\pm 40\text{mg/dL}$	Outside $\pm 40\%$ / $\pm 40\text{mg/dL}$
1	1497	76.2	87.4	97.9	99.5	0.5
4	1470	82.3	91.4	97.6	99.3	0.7

7	1394	85.0	93.5	98.6	99.3	0.7
10	1411	85.3	92.3	97.9	99.4	0.6

System Glucose Availability

The System is designed to produce a glucose reading after each user initiated scan that is performed after the start-up time, through the wear period. [Table 11](#) shows the number of available glucose readings reported by all Sensors and the expected number based on the total number of scan attempts. Results are shown for Sensors which produced at least one CGM reading during the clinical study over the total wear period. The percentage of available CGM readings is presented in comparison to the number of expected CGM readings. Overall, 99.5% of CGM readings (9,228 CGM readings out of an expected 9,272) were available.

Table 11. CGM Availability

No. CGM	No. Scan	%
9228	9272	99.5

Detection of Hypoglycemic and Hyperglycemic Events

[Table 12](#) shows the accuracy of the System's Glucose Messages in informing the user of low or high glucose events within 15 minutes before or after the true low or high blood glucose value. Percentages are displayed for three different parameters:

- Detection Rate – amount of time the System displays a Glucose Message correctly.
- Missed Detection Rate – amount of time the System did not display a Glucose Message when it should have.
- False Notification Rate – amount of time the System displays a Glucose Message when it shouldn't have.

For example, in the clinical study, the System was able to detect 85.4% of actual low glucose events (detection rate), but 39.9% of the time a Low Glucose message was displayed in error (false notification rate) and 14.6% of the time a Low Glucose message was not displayed when it should have been (missed detection rate).

Table 12. Detection of Hypoglycemic and Hyperglycemic Events

Type of Notification	Notification Status	15 Minute Interval
Notification of Hypoglycemic Events (Low Glucose message)	Detection Rate (%)	85.4
	Missed Detection Rate (%)	14.6
	False Notification Rate (%)	39.9
	Detection Rate (%)	95.1

Notification of Hyperglycemic Events (High Glucose message)	Missed Detection Rate (%)	4.9
	False Notification Rate (%)	22.1
Impending Notification of Hypoglycemic Events (Glucose Going Low message)	Detection Rate (%)	95.0
	Missed Detection Rate (%)	5.0
	False Notification Rate (%)	46.8
Impending Notification of Hyperglycemic Events (Glucose Going High message)	Detection Rate (%)	97.2
	Missed Detection Rate (%)	2.8
	False Notification Rate (%)	28.4

Precision

Precision of the System was evaluated by comparing the results from two separate Sensors worn on the same subject at the same time. [Table 13](#) provides data from two separate Sensors worn on 47 subjects at the same time, providing 7,319 real-time pairs of CGM measurements, with a mean CV of 6.0%.

Table 13. Overall between Sensor Precision

Mean Glucose (mg/dL)	Median CV	Mean CV	Number of Subjects	Number of Paired Readings
175.3	4.6	6.0	47	7319

Sensor Wear Duration

After the 12 hour start-up period, the Sensor can be worn for up to 10 days. To estimate how long a Sensor will work over the wear duration, 97 Sensors were evaluated in the clinical study to determine how many days of readings each Sensor provided. Of these 97 Sensors, 75 (77.3%) lasted until the final day of use. 84 Sensors (86.6%) lasted more than 5 days. There were 22 (22.7%) Sensors that failed early, of which 11 (11.3%) failed on or before the fifth day of wear.

Adverse Events

No device related serious adverse events occurred during the study. Mild skin irritations, such as erythema, edema, rash, bleeding, itching, induration, and infection were reported around the insertion site and adhesive area by a moderate frequency of subjects (5 out of 48 or 10.4%). Pain was mostly reported as none with only one instance of mild pain.

Performance Characteristics - FreeStyle Libre 14 day Sensor

Clinical Study Overview

Performance of the FreeStyle Libre 14 day Flash Glucose Monitoring System (the System) was evaluated in a clinical study. The study was conducted at 4 centers with a total of 95 subjects with diabetes (84.2% Type 1, 15.8% Type 2). All subjects were aged eighteen and older. Subjects in the study required insulin to manage their diabetes. Each subject wore up to two Sensors on the back of the upper arm. During the study, subjects tested their blood glucose using fingerstick capillary samples at least eight times during each day of the study. Subjects used the blood glucose meter built into the Reader. Additionally, subjects had their venous blood glucose analyzed up to 112 times over four separate visits to the clinical center. Venous blood was analyzed using the Yellow Springs Instrument Life Sciences 2300 STAT Plus™ Glucose & Lactate Analyzer (YSI). YSI is a laboratory glucose and lactate analyzer of whole blood and plasma and is a widely recognized standard in laboratory analysis of blood glucose. Glucose readings obtained from the System were compared to glucose readings obtained from the YSI to evaluate the performance of the System. Three lots of Sensors were evaluated in the study.

Agreement with YSI Levels

Agreement between the System Glucose Measurement (CGM) and venous blood was characterized by using paired CGM and Yellow Springs Instrument measurements (YSI). The accuracy of CGM versus YSI reference was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for reference values 80 mg/dL and above, and 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL when glucose levels are assigned using the YSI values. Overall 90.7% of results were within ± 20 mg/dL / 20% of YSI reference.

Agreement with CGM Glucose Levels

Agreement between CGM and venous blood was characterized by using paired CGM and Yellow Springs Instrument measurements (YSI). The accuracy of CGM versus YSI reference was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for reference values 80 mg/dL and above, and 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL. The results are presented in [Table 14](#) for YSI reference. Overall 90.9% of results were within ± 20 mg/dL / 20% of YSI reference.

Table 14. Number and Percent of Results within YSI Reference

CGM Glucose Level (mg/dL)	Number of CGM-Reference Pairs	Within $\pm 15\%$ / ± 15 mg/dL	Within $\pm 20\%$ / ± 20 mg/dL	Within $\pm 30\%$ / ± 30 mg/dL	Within $\pm 40\%$ / ± 40 mg/dL	Outside $\pm 40\%$ / ± 40 mg/dL
Overall	9725	83.0	90.9	97.3	99.1	0.9
40-50	40	52.5	60.0	80.0	85.0	15.0
51-80	741	58.2	70.9	87.9	94.6	5.4
81-180	6112	82.2	90.9	97.5	99.4	0.6
181-300	2513	91.7	96.9	99.4	99.8	0.2

301-400	291	89.0	95.2	99.7	99.7	0.3
401-500	28	100.0	100.0	100.0	100.0	0.0

Agreement on Day 1 against YSI Reference

The accuracy of CGM versus YSI reference on the first day of Sensor wear was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for reference values 80 mg/dL and above, and within 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL by hourly intervals. The results are presented in [Table 15](#)

Table 15. Number and Percent of Results within YSI Reference

Time Interval (hours)	Number of CGM-Reference Pairs	Within $\pm 15\%$ / $\pm 15\text{mg/dL}$	Within $\pm 20\%$ / $\pm 20\text{mg/dL}$	Within $\pm 30\%$ / $\pm 30\text{mg/dL}$	Within $\pm 40\%$ / $\pm 40\text{mg/dL}$	Outside $\pm 40\%$ / $\pm 40\text{mg/dL}$
(0-2] [*]	182	77.5	87.9	95.6	99.5	0.5
(2-4]	385	71.2	83.1	94.8	98.4	1.6
(4-6]	375	77.6	86.4	91.7	96.3	3.7
(6-8]	373	78.0	87.9	97.3	99.5	0.5
(8-16] [†]	106	74.5	83.0	97.2	100.0	0.0
(16-18]	303	80.2	89.1	96.7	99.7	0.3
(18-20]	344	80.2	86.6	94.8	98.8	1.2
(20-22]	336	81.3	87.8	97.9	98.8	1.2
(22-24]	155	86.5	93.5	96.1	98.7	1.3

* (0-2] interval includes the 1 hour start-up time.

† Time interval is not divided to smaller segments due to small number of data pairs.

Overall Accuracy against YSI reference

Accuracy was measured by comparing the absolute relative difference between the System and reference YSI glucose values. The absolute relative difference measures the level of disagreement between the System and the reference value, but does not tell you whether the System glucose value was, on average, higher or lower than the reference glucose value. The Mean Absolute Relative Difference gives an indication of the average percent disagreement between the CGM and the reference. [Table 16](#) shows the overall absolute difference measure. Overall the Mean Absolute Relative Difference was 9.4% for the comparison with YSI reference. The Median Absolute Relative Difference shows that half of the time the System was within 7.4% of the YSI reference.

Table 16. Difference Measures with YSI Reference

	Median Absolute Relative Difference	
--	-------------------------------------	--

Number of CGM-Reference Pairs	(%)	Mean Absolute Relative Difference (%)
9725	7.4	9.4

Agreement with BG Levels

Agreement between the System and capillary blood glucose values (BG) as measured by the Reader's built-in meter was characterized by using paired System CGM and BG value. The accuracy of CGM versus BG value was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for BG values 80 mg/dL and above, and within 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL. The results are presented in [Table 17](#) for BG values. Overall 87.0% of results were within ± 20 mg/dL / 20% of BG values.

Table 17. Number and Percent of Results within BG Values*

CGM Glucose Level (mg/dL)	Number of CGM-Reference Pairs	Within $\pm 15\%$ / ± 15 mg/dL	Within $\pm 20\%$ / ± 20 mg/dL	Within $\pm 30\%$ / ± 30 mg/dL	Within $\pm 40\%$ / ± 40 mg/dL	Outside $\pm 40\%$ / ± 40 mg/dL
Overall	9234	75.8	87.0	96.1	98.6	1.4
40-50	123	52.0	62.6	78.0	89.4	10.6
51-80	795	60.5	71.9	88.2	95.0	5.0
81-180	4839	73.3	85.9	96.4	98.9	1.1
181-300	2800	83.8	92.8	98.4	99.4	0.6
301-400	563	82.4	93.6	98.0	99.3	0.7
401-500	114	80.7	93.0	97.4	98.2	1.8

* Comparison to BG was performed using the Reader's built-in blood glucose meter. Different performance may be expected when compared to other models of blood glucose meters.

Overall Accuracy against BG values

Accuracy was measured by comparing the absolute relative difference between the System and BG values. The absolute relative difference measures the level of disagreement between the System and the BG value, but does not tell you whether the System glucose value was, on average, higher or lower than the BG value. The Mean Absolute Relative Difference gives an indication of the average percent disagreement between the CGM and the BG value.

[Table 18](#) shows the overall absolute difference measure. Overall the Mean Absolute Relative Difference was 11.4% for the comparison with BG value. The Median Absolute Relative Difference shows that half of the time the System was within 9.1% of the BG value.

Table 18. Difference Measures with BG Value *

Number of CGM-Reference Pairs	Median Absolute Relative Difference (%)	Mean Absolute Relative Difference (%)
9234	9.1	11.4

* Comparison to BG was performed using the Reader's built-in blood glucose meter. Different performance may be expected when compared to other models of blood glucose meters.

[illegible]

CGM	YSI (mg/dL/min)	N
CGM1	10.0	10
CGM2	10.0	10
CGM3	10.0	10
CGM4	10.0	10
CGM5	10.0	10
CGM6	10.0	10
CGM7	10.0	10
CGM8	10.0	10
CGM9	10.0	10
CGM10	10.0	10
CGM11	10.0	10
CGM12	10.0	10
CGM13	10.0	10
CGM14	10.0	10
CGM15	10.0	10
CGM16	10.0	10
CGM17	10.0	10
CGM18	10.0	10
CGM19	10.0	10
CGM20	10.0	10
CGM21	10.0	10
CGM22	10.0	10
CGM23	10.0	10
CGM24	10.0	10
CGM25	10.0	10
CGM26	10.0	10
CGM27	10.0	10
CGM28	10.0	10
CGM29	10.0	10
CGM30	10.0	10
CGM31	10.0	10
CGM32	10.0	10
CGM33	10.0	10
CGM34	10.0	10
CGM35	10.0	10
CGM36	10.0	10
CGM37	10.0	10
CGM38	10.0	10
CGM39	10.0	10
CGM40	10.0	10
CGM41	10.0	10
CGM42	10.0	10
CGM43	10.0	10
CGM44	10.0	10
CGM45	10.0	10
CGM46	10.0	10
CGM47	10.0	10
CGM48	10.0	10
CGM49	10.0	10
CGM50	10.0	10
CGM51	10.0	10
CGM52	10.0	10
CGM53	10.0	10
CGM54	10.0	10
CGM55	10.0	10
CGM56	10.0	10
CGM57	10.0	10
CGM58	10.0	10
CGM59	10.0	10
CGM60	10.0	10
CGM61	10.0	10
CGM62	10.0	10
CGM63	10.0	10
CGM64	10.0	10
CGM65	10.0	10
CGM66	10.0	10
CGM67	10.0	10
CGM68	10.0	10
CGM69	10.0	10
CGM70	10.0	10
CGM71	10.0	10
CGM72	10.0	10
CGM73	10.0	10
CGM74	10.0	10
CGM75	10.0	10
CGM76	10.0	10
CGM77	10.0	10
CGM78	10.0	10
CGM79	10.0	10
CGM80	10.0	10
CGM81	10.0	10
CGM82	10.0	10
CGM83	10.0	10
CGM84	10.0	10
CGM85	10.0	10
CGM86	10.0	10
CGM87	10.0	10
CGM88	10.0	10
CGM89	10.0	10
CGM90	10.0	10
CGM91	10.0	10
CGM92	10.0	10
CGM93	10.0	10
CGM94	10.0	10
CGM95	10.0	10
CGM96	10.0	10
CGM97	10.0	10
CGM98	10.0	10
CGM99	10.0	10
CGM100	10.0	10

(mg/dL/min)	<-2	[-2, -1)	[-1, 0)	[0, 1]	(1, 2]	>2	NA [*]	
<-2 (↓)	39.9	42.7	8.4	0.7	0.0	0.0	8.4	143
-2 to -1 (↘)	4.7	27.2	56.1	5.5	0.3	0.0	6.2	695
-1 to 1 (→)	0.2	2.6	50.2	36.4	4.2	0.6	5.8	7786
1 to 2 (↗)	0.0	0.6	5.6	43.4	36.3	9.8	4.5	717
>2 (↑)	0.0	0.0	1.6	12.6	37.2	43.3	5.3	247
NA [†]	0.6	8.3	27.4	37.5	13.1	4.2	8.9	168

* Glucose rate of change not available due to the time difference between glucose readings exceeding 30 minutes.

† Glucose Trend Arrow not available.

Agreement with 'LO' and 'HI' CGM Reading against YSI Reference

The System reports glucose concentrations between 40 and 500 mg/dL. When the System determines that glucose level is below 40 mg/dL, it will report as 'LO'. When the System determines that glucose level is above 500 mg/dL, it will report as 'HI'. [Table 21](#) displays the concurrence between the CGM and YSI reference glucose when CGM reads 'LO'. For example, in the clinical study, when CGM reading was 'LO', YSI glucose values were less than 40 mg/dL 0.0% of the time, equal to or above 40 mg/dL 100.0% of the time, above 50 mg/dL 83.3% of the time, above 60 mg/dL 55.6% of the time, above 70 mg/dL 33.3% of the time, and above 80 mg/dL 22.2% of the time.

Table 21. Concurrence Analysis with 'LO' CGM Reading

	YSI (mg/dL)						Total
	<40	≥40	>50	>60	>70	>80	
% of CGM points in YSI range	0.0	100.0	83.3	55.6	33.3	22.2	
Number of CGM points in YSI range	0	18	15	10	6	4	18

[Table 22](#) displays the concurrence between the CGM and YSI reference glucose when CGM reads 'HI'. For example, in the clinical study, when CGM reading was 'HI', YSI glucose values were less than or equal to 200 mg/dL 0% of the time, above 200 mg/dL 100% of the time, above 300 mg/dL 100% of the time, above 400 mg/dL 100% of the time, and above 500 mg/dL 7.7% of the time.

Table 22. Concurrence Analysis with 'HI' CGM Reading

	YSI (mg/dL)					N
	≤200	>200	>300	>400	>500	
% of CGM points in YSI	0.0	100.0	100.0	100.0	7.7	

range						
Number of CGM points in YSI range	0	13	13	13	1	13

Accuracy by Day of Wear

The Sensor can be worn for up to 14 days. To show Sensor performance over time, the absolute relative difference between the System and reference YSI glucose values over the wear duration is presented in [Table 23](#).

Table 23. Difference Measures by Day

Day	Number of CGM-Reference Pairs	Median Absolute Relative Difference (%)	Mean Absolute Relative Difference (%)
1	2563	8.3	10.8
6	2545	7.1	8.5
11	2419	7.4	9.3
14	2198	6.6	9.1

The accuracy of CGM versus YSI reference and BG reference was assessed by calculating the percentage of System readings that were within 15%, 20%, 30% and 40% for reference values 80 mg/dL and above, and 15 mg/dL, 20 mg/dL, 30 mg/dL and 40 mg/dL for values below 80 mg/dL. The results for CGM vs. YSI reference are presented in [Table 24](#).

Table 24. Number and Percent of Results within YSI Reference

Day	Number of CGM-Reference Pairs	Within $\pm 15\%$ / $\pm 15\text{mg/dL}$	Within $\pm 20\%$ / $\pm 20\text{mg/dL}$	Within $\pm 30\%$ / $\pm 30\text{mg/dL}$	Within $\pm 40\%$ / $\pm 40\text{mg/dL}$	Outside $\pm 40\%$ / $\pm 40\text{mg/dL}$
1	2563	78.2	87.0	95.6	98.7	1.3
6	2545	87.0	94.5	99.4	99.8	0.2
11	2419	82.3	90.0	98.2	99.9	0.1
14	2198	83.9	91.3	96.2	99.0	1.0

System Glucose Availability

The System is designed to produce a glucose reading after each user initiated scan that is performed throughout the wear period after the start-up time. [Table 25](#) shows the number of available glucose readings reported by all Sensors and the expected number based on the total number of scan attempts. Results are shown for Sensors which produced at least one CGM reading during the clinical study over the total wear period. The percentage of available CGM readings is presented in comparison to the number of expected CGM readings. Overall, 99.6% (18,488 CGM readings out of an expected 18,562) were available.

Table 25. CGM Availability

No. CGM	No. Scan	%
18488	18562	99.6

Detection of Hypoglycemic and Hyperglycemic Events

Table 26 shows the accuracy of the System's Glucose Messages in informing the user of low or high glucose events within 15 minutes before or after the true low or high blood glucose value. Percentages are displayed for three different parameters:

- Detection Rate – amount of time the System displays a Glucose Message correctly.
- Missed Detection Rate – amount of time the System did not display a Glucose Message when it should have.
- False Notification Rate – amount of time the System displays a Glucose Message when it shouldn't have.

For example, in the clinical study, the System was able to detect 77.6% of actual low glucose events (detection rate), but 52.5% of the time a Low Glucose message was displayed in error (false notification rate) and 22.4% of the time a Low Glucose message was not displayed when it should have been (missed detection rate).

Table 26. Detection of Hypoglycemic and Hyperglycemic Events

Type of Notification	Notification Status	15 Minute Interval
Notification of Hypoglycemic Events (Low Glucose message)	Detection Rate (%)	77.6
	Missed Detection Rate (%)	22.4
	False Notification Rate (%)	52.5
Notification of Hyperglycemic Events (High Glucose message)	Detection Rate (%)	84.7
	Missed Detection Rate (%)	15.3
	False Notification Rate (%)	10.6
Impending Notification of Hypoglycemic Events (Glucose Going Low message)	Detection Rate (%)	91.6
	Missed Detection Rate (%)	8.4
	False Notification Rate (%)	58.8
Impending Notification of	Detection Rate (%)	88.8

Hyperglycemic Events (Glucose Going High message)	Missed Detection Rate (%)	11.2
	False Notification Rate (%)	15.5

Precision

Precision of the System was evaluated by comparing the results from two separate Sensors worn on the same subject at the same time. [Table 27](#) provides data from two separate Sensors worn on 95 subjects at the same time, providing 15,283 real-time pairs of CGM measurements, with a mean CV of 5.6%.

Table 27. Overall between Sensor Precision

Mean Glucose (mg/dL)	Median CV	Mean CV	Number of Subjects	Number of Paired Readings
164.2	4.2	5.6	95	15283

Sensor Wear Duration

The Sensor can be worn for up to 14 days. To estimate how long a Sensor will work over the wear duration, 190 Sensors were evaluated in the clinical study to determine how many days of readings each Sensor provided. Of these 190 Sensors, 136 (71.6%) lasted until the final day of use. 168 Sensors (88.4%) lasted at least 7 days. There were 54 (28.4%) Sensors that failed early, of which 22 (11.6%) failed on or before the seventh day of wear.

Adverse Events

No device related serious adverse events occurred during the study. Mild skin irritations, such as erythema, bruising, bleeding, infection and papule were reported around the insertion site and adhesive area by a small number of subjects in the study (8 out of 95 or 8.4%).

Electromagnetic Compatibility

- The Sensor needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual.
- Portable and mobile RF communications equipment can affect the Sensor.
- The use of accessories, transducers and cables other than those specified by Abbott Diabetes Care may result in increased EMISSIONS or decreased IMMUNITY of the Sensor.
- The Sensor should not be used adjacent to or stacked with other equipment and that if

adjacent or stacked use is necessary, the Sensor should be observed to verify normal operation in the configuration in which it will be used.

- The Sensor complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) The Sensor may not cause harmful interference, and (2) the Sensor must accept any interference received, including interference that may cause undesired operation.
- Changes or modifications not approved by Abbott could void the user's authority to operate the equipment.

Guidance and manufacturer's declaration – electromagnetic emissions

The System is intended for use in the electromagnetic environment specified below. The customer or the user of the System should assure that it is used in such an environment.

Emissions test: RF emissions; CISPR 11

Compliance: Group 1

Electromagnetic environment – guidance: The System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Emissions test: RF emissions; CISPR 11

Compliance: Class B

Electromagnetic environment – guidance: The System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity

The Sensor is intended for use in the electromagnetic environment specified below. The customer or the user of the Sensor should assure that it is used in such an environment.

Immunity test: Electrostatic discharge (ESD); IEC 61000-4-2

IEC 60601 test level: ± 8 kV contact; ± 15 kV air

Compliance level: ± 8 kV contact; ± 15 kV air

Electromagnetic environment – guidance: Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Immunity test: Power frequency (50/60 Hz); magnetic field; IEC 61000-4-8

IEC 60601 test level: 30 A/m

Compliance level: 30 A/m

Electromagnetic environment – guidance: Power frequency magnetic fields should be at levels characteristic of a typical location in a typical domestic, commercial, or hospital environment.

Immunity test: Radiated RF; IEC 61000-4-3

IEC 60601 test level: 10 V/m; 80 MHz to 2.7 GHz

Compliance level: 10 V/m

Electromagnetic environment – guidance:

Recommended separation distance

$$d = 1.2 \sqrt{P}$$

80 MHz to 800 MHz

$$d = 2.3 \sqrt{P}$$

800 MHz to 2.5 GHz

P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b

Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Sensor is used

exceeds the applicable RF compliance level above, the Sensor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Sensor.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Sensor

The Sensor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Sensor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Sensor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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©2023 Abbott
The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott. Other trademarks are the property of their respective owners.