**Owner's Booklet  
Manual del Usuario**

# FreeStyle


## Freedom Lite

Blood Glucose Monitoring System  
Sistema de Control de Glucosa en la Sangre

 **Abbott**

# The FreeStyle Freedom Lite® Meter and How It Works

## System Check Screen

This display always appears when the meter is turned on. You should check that your meter matches the example exactly every time your meter turns on. **Do not use the meter if the display check screen does not exactly match the example.** If the screen does

not match, the meter may show an incorrect result. Please contact Customer Care.

## Display Screen

Displays your test results and other important information.

## m (Mode) Button

- Moves to different mode settings
- Scrolls backward
- Silences a reminder alarm

## Sample Areas (dark-colored half-circles on test strip)

Apply blood or control solution to one sample area only. Insert strip into the meter printed side up (see image to the right).

## Data Port

Downloads test results. Requires a compatible data management system. *For more information, please go to our website or contact Customer Care.*

## c (Configure) Button

- Marks a control solution test
- Scrolls forward
- Silences a reminder alarm
- Changes the date, time, alarm settings, sound options, and reviews stored readings

## Test Strip Port

Insert the **Top** end of a new FreeStyle Lite® test strip here. The meter powers on when you insert the test strip.

Top

FreeStyle Lite test strip

Bottom

Sample area

# Table of Contents

|                                                                         |           |
|-------------------------------------------------------------------------|-----------|
| <b>The FreeStyle Freedom Lite® Meter and How It Works</b>               | <b>i</b>  |
| <b>Important Information About Your Blood Glucose Monitoring System</b> | <b>1</b>  |
| Intended Use                                                            | 1         |
| <b>Setting Up the Meter</b>                                             | <b>4</b>  |
| Purpose                                                                 | 4         |
| Set Sound Volume                                                        | 5         |
| Set Time                                                                | 6         |
| Set Date                                                                | 7         |
| <b>Testing Your Blood Glucose</b>                                       | <b>9</b>  |
| Purpose                                                                 | 9         |
| Prepare to Test                                                         | 9         |
| Important Test Strip Information                                        | 11        |
| Performing a Blood Glucose Test                                         | 12        |
| <b>Understanding Blood Glucose Test Results</b>                         | <b>14</b> |
| Low and High Blood Glucose Results                                      | 15        |
| <b>Error Codes</b>                                                      | <b>18</b> |
| <b>Troubleshooting</b>                                                  | <b>21</b> |

|                                                              |           |
|--------------------------------------------------------------|-----------|
| <b>FreeStyle® Control Solution</b> .....                     | <b>23</b> |
| Purpose .....                                                | <b>23</b> |
| Control Solution Information .....                           | <b>23</b> |
| Out of Range Control Solution Results .....                  | <b>24</b> |
| Performing a Control Solution Test .....                     | <b>25</b> |
| <b>Using Your Meter Memory</b> .....                         | <b>28</b> |
| Purpose .....                                                | <b>28</b> |
| Viewing Blood Glucose Averages in Memory .....               | <b>28</b> |
| Viewing Test Results in Memory .....                         | <b>29</b> |
| <b>Setting Reminder Alarms (optional)</b> .....              | <b>30</b> |
| Purpose .....                                                | <b>30</b> |
| Transferring Test Results To a Computer .....                | <b>31</b> |
| <b>Maintaining the Meter</b> .....                           | <b>32</b> |
| Replacing the Battery .....                                  | <b>32</b> |
| Cleaning and Disinfecting the Meter and Lancing Device ..... | <b>34</b> |
| <b>Specifications</b> .....                                  | <b>39</b> |
| <b>Definition of Symbols</b> .....                           | <b>41</b> |

# Important Information About Your Blood Glucose Monitoring System

## Intended Use

The FreeStyle Freedom Lite Blood Glucose Monitoring System is intended for self-testing outside the body (*in vitro* diagnostic use) by people with diabetes at home as an aid in monitoring the effectiveness of a diabetes control program by quantitatively measuring glucose in fresh whole blood from the finger, upper arm and palm. The FreeStyle Freedom Lite Blood Glucose Monitoring System should not be used for the diagnosis of or screening of diabetes or for neonatal use. It is intended to be used by a single person and should not be shared. Alternate site testing should be done only during steady state times (when glucose is not changing rapidly). The FreeStyle Lite Blood Glucose Test Strips are for use with the FreeStyle Freedom Lite Blood Glucose Meter to quantitatively measure glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips, upper arm and palm.

### IMPORTANT:

The device must not be used on more than one person including other family members due to the risk of spreading infection. All parts of the device and its accessories are considered biohazardous and can potentially transmit infectious diseases, even after performing the cleaning procedure.<sup>1,2</sup>

If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be cleaned and disinfected prior to use by the second person.

Use only FreeStyle Lite® test strips and FreeStyle® control solution with this meter. Other products can produce inaccurate results.

<sup>1</sup>"FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010) <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>.

<sup>2</sup>"CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010) <http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

- The expected glucose level for a non-diabetic, non-pregnant fasting adult is under 100 mg/dL.<sup>3</sup>
- Two hours after meals, levels should be less than 140 mg/dL.<sup>3</sup>
- Consult your healthcare professional to determine the range that is appropriate for you.

**Warning:**

- Test on your fingers if you are testing for low blood glucose (hypoglycemia) or if you suffer from hypoglycemia unawareness.
- This system contains small parts that may be dangerous if swallowed.

**Caution:** Please read all the instructions provided in this owner's booklet and practice the testing procedures before using the FreeStyle Freedom Lite® system. Blood glucose monitoring should be done with the guidance of a healthcare professional.

<sup>3</sup>American Diabetes Association. (2015). 2. Classification and Diagnosis of Diabetes. Diabetes Care, 38(Supplement 1), S8-S16.

**Limitations:**

- For *in-vitro* diagnostic use only (external use only).
- Not for diagnosis of or screening for diabetes.
- Not for testing of neonates.
- The test strips are for single use only. Do not reuse test strips.
- The system can give false results greater than 10,000 feet (3,048 meters) above sea level.
- The system can give false results outside the hematocrit range: 15% to 65%.
- Individuals experiencing a hypoxia state may yield falsely low results.
- In situations of decreased peripheral blood flow, capillary blood testing may not be appropriate as it may not represent the true physiological state. Examples would include but are not limited to: severe dehydration caused by diabetic ketoacidosis or the hyperosmolar non-ketotic state, hypotension, shock, or peripheral vascular disease.<sup>4, 5, 6</sup>
- This device should not be used to test critically ill patients.
- Alternative site testing (AST) should not be used to calibrate continuous glucose monitors (CGMs) nor for use in insulin dose calculations.
- AST testing should only be done during steady-state times (when glucose is not changing rapidly).
- For single patient use at home by a lay user.
- Severe dehydration or excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult your physician immediately.
- Cholesterol up to 500 mg/dL or triglycerides up to 3,000 mg/dL do not significantly affect test results. However, glucose values in specimens beyond these levels, should be interpreted with caution.

<sup>4</sup> Atkins SH, Dasmahapatra A, Jaker MA, Chorost MI, Reddy S: Fingerstick glucose determination in shock. Ann Int Med 114: 1020-1024. 1991.

<sup>5</sup> Sandler M, Low-Beer T: Misleading capillary glucose measurements. Practical Diabetes 7: 210. 1990.

<sup>6</sup> Wickham NWR, Achar KN, Cove DH: Unreliability of capillary blood glucose in peripheral vascular disease. Practical Diabetes 3: 100. 1986.

# Setting Up the Meter

## Purpose

The purpose of this section is to set the meter sound volume, time and date.


1. With the meter off, press and **hold** (m) for 3 seconds until SET appears.

2. Press (m) to bypass setting the optional reminder alarms. You may choose to set these later (see Setting Reminder Alarms section).

## Set Sound Volume

With the sound volume on, your meter will beep when:



- The test strip is full of blood
- The test result appears on the display
- An error occurs

You may also choose to silence your meter.

### 3. Set sound volume


| If you prefer: | Press <b>(C)</b> until you see: | Then                       |
|----------------|---------------------------------|----------------------------|
| No sound       |                                 | • Press <b>(m)</b> to save |
| A low beep     |                                 |                            |
| A louder beep  |                                 |                            |

## Set Time

Time Format

**Note:** Your meter can display either the 12-hour (1:24P) or the 24-hour (13:24) time formats. If you prefer the 12-hour format, there is no “A” for AM. If setting a PM time, continue to press **(C)** until you see the “P” displayed for PM.

Hour

### 4. Set hour

- Press **(C)** to display correct hour
- Press **(m)** to save

Minutes

### 5. Set minutes

- Press **(C)** to display correct minutes
- Press **(m)** to save


Time Format

### 6. Set time format

- Press **(m)** to accept time format or
- Press **(C)** to change
- Press **(m)** to save

## Set Date

**Note:** The date format can display either:

The month will flash first

### 7. Set month

- Press **C** to display correct month
- Press **m** to save

### 8. Set day

- Press **C** to display correct day
- Press **m** to save

Date Format

**9.** Set date format

- Press **(M)** to accept date format or
- Press **(C)** to change
- Press **(M)** to save

Year

**10.** Set year

- Press **(M)** to accept year or
- Press **(C)** to change
- Press **(M)** to save

**11.** Press and **hold** **(M)** for **3** seconds until meter turns off. Meter set up is complete.

# Testing Your Blood Glucose

## Purpose

The purpose of this section is to produce an accurate blood glucose test result.

## Prepare to Test

### 1. Select your test site.

**Note:** Avoid moles, veins, bones, and tendons. Bruising may occur at the test site. If bruising occurs, consider selecting another site.

There are differences in testing on fingers versus alternative sites.

### **Warning:**

Test on your fingers if you are testing for low blood glucose (hypoglycemia) or if you suffer from hypoglycemia unawareness.

Check for low blood glucose when:

- You feel symptoms such as weakness, sweating, shakiness, headache, or confusion
- You eat a delayed meal after taking insulin
- You are following your healthcare professional's advice

**IMPORTANT:**

- 2.** Wash your hands and the test site with warm water and soap to ensure accurate results.
- Thoroughly dry your hands and the test site.
  - Do *not* use lotion or cream on the test site.

- 3.** Check test strip expiration date.  
Do *not* use expired test strips; they may produce inaccurate results.

**IMPORTANT:** Use only FreeStyle Lite® test strips. Other test strips can produce inaccurate results. The test strips are for single use only.


## Important Test Strip Information

**Please read the FreeStyle Lite® test strip package insert for more information.**

- Use the test strips within the meter's operating range: 40° – 104° F (4° – 40° C).
- Store test strips in a cool, dry place between 40° – 86° F (4° – 30° C).
- Store test strips away from direct sunlight and heat.
- Store test strips in their original vial only.
- Clean and dry your hands before removing a strip from the vial.
- Close the vial cap tightly immediately after removing a test strip.
- Do *not* transfer test strips to another container.
- Do *not* store test strips outside the vial.
- Do *not* bend, cut, or alter test strips.
- Do *not* apply blood to both sides of the test strip.
- Do *not* press the test strip against the test site.
- Do *not* scrape the blood onto the test strip.
- Do *not* apply blood to the flat side of the test strip.
- Do *not* apply blood to the test strip when the test strip is out of the meter.
- Do *not* put blood or foreign objects into the meter.

**Warning:** Drying agents in the test strip vial may be harmful if inhaled or swallowed and may irritate skin and eyes.


# Performing a Blood Glucose Test


1. Remove test strip from vial.
2. Insert test strip to turn meter on.

**Note:** The meter turns off after 2 minutes of inactivity. Remove and reinsert the unused test strip to restart the meter.

If the meter does not turn on, go to Troubleshooting section.

3. Confirm System Check Screen.
  - This display appears when you turn on the meter. **Do not use the meter if the display check screen does not exactly match the example. Contact Customer Care.**
  - See The FreeStyle Freedom Lite® Meter and How It Works section for more information.
  - If the system check screen does not appear, go to Troubleshooting section.

4. Obtain a blood drop.

When the  and  appears, use the lancing device to obtain a blood drop. (See FreeStyle® lancing device insert for instructions.)

The test strip requires only a 0.3 microliter of blood for accurate results.

**5.** Apply and **hold** blood to 1 sample area of the test strip until:

- You see short lines moving clockwise on the display or
- You hear a beep

This indicates the test strip obtained enough blood and the meter is checking your glucose level. If a short line does not appear after 5 seconds, the sample may be too small. You can add blood to the same side only for up to 60 seconds. If the test does not start after applying blood drop, go to Troubleshooting section.

(Example only)

**6.** View result

The result appears on the display when the test is complete. The time that the meter takes to display a result depends on your blood glucose level. Higher glucose levels require more time.

**Note:** If you see an “**Er1, Er2, Er3, or Er4**” on the display, it is an error code. Go to Error Codes section.

**7.** Remove test strip to turn meter off.

Discard used lancet and test strip properly.

**8.** Wash your hands and test site thoroughly with warm water and soap after handling the meter, lancing device or test strips.

## Understanding Blood Glucose Test Results

The meter displays results in mg/dL. The unit of measurement is preset. You cannot change this setting.

**IMPORTANT:** The meter displays results from 20 – 500 mg/dL. Low or high blood glucose results can indicate a potentially serious medical condition.

## Low and High Blood Glucose Results

| Result                   | What It Means        | What To Do                                                                                                                                                                                                                                                                  |
|--------------------------|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| lower than<br>60 mg/dL   | Low (hypoglycemic)   | <p>Repeat the test with a new test strip, even if you do <i>not</i> feel that you have low blood glucose.</p> <p>If your blood glucose result is still not consistent with your symptoms, contact your healthcare professional and follow his or her treatment advice.</p>  |
| higher than<br>240 mg/dL | High (hyperglycemic) | <p>Repeat the test with a new test strip, even if you do <i>not</i> feel that you have high blood glucose.</p> <p>If your blood glucose result is still not consistent with your symptoms, contact your healthcare professional and follow his or her treatment advice.</p> |

| Display                                                                                                                                | What It Means            | Symptoms                                                   | What To Do                                                                                                                                                                                                                        |
|----------------------------------------------------------------------------------------------------------------------------------------|--------------------------|------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <br>Appears when result is:<br>lower than<br>20 mg/dL | Severe low blood glucose | Weakness<br>Sweating<br>Shakiness<br>Headache<br>Confusion | Follow your healthcare professional's advice to treat low blood glucose.<br>If your blood glucose result does not match how you feel, perform a control solution test to confirm that your meter and strips are working properly. |

**IMPORTANT:**

- Severe dehydration may cause false low results. If you believe you are suffering from severe dehydration, contact your healthcare professional immediately.
- Contact your healthcare professional if you have symptoms that do *not* match your test result, and you have followed the instructions in this owner's booklet.

| Display                                                                                                                                  | What It Means             | Symptoms                                               | What To Do                                                                                                                                                                                                                         |
|------------------------------------------------------------------------------------------------------------------------------------------|---------------------------|--------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <br>Appears when result is:<br>higher than<br>500 mg/dL | Severe high blood glucose | Fatigue<br>Thirst<br>Excess Urination<br>Blurry Vision | Follow your healthcare professional's advice to treat high blood glucose.<br>If your blood glucose result does not match how you feel, perform a control solution test to confirm that your meter and strips are working properly. |

## Error Codes

| Message | What It Means                                                                                                                                                                                             | What To Do                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|---------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Er 1    | <ul style="list-style-type: none"><li>• The sample is too small</li><li>• There may be a problem with the test strip</li><li>• Meter error</li><li>• Very low blood glucose lower than 20 mg/dL</li></ul> | <ol style="list-style-type: none"><li>1. Follow your healthcare professional's advice to treat low blood glucose if you have symptoms such as:<ul style="list-style-type: none"><li>• Weakness</li><li>• Shakiness</li><li>• Sweating</li><li>• Headache</li><li>• Confusion</li></ul></li><li>2. Perform a control solution test using a new test strip. If the test results are within the range printed on the test strip vial, retest using blood and a new test strip.</li><li>3. If the control solution result is out of range or the error reappears, contact Customer Care.</li></ol> |

| Message                                                                           | What It Means                                                                                                                                                                                                                                                                                                 | What To Do                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <ul style="list-style-type: none"><li>• There may be a problem with the test strip</li><li>• Meter error</li><li>• Very high blood glucose higher than 500 mg/dL</li><li>• HIGH control solution was applied when the temperature was too cold. (Applies only to control solution labeled as HIGH).</li></ul> | <ol style="list-style-type: none"><li>1. Follow your healthcare professional's advice to treat high blood glucose if you have symptoms such as:<ul style="list-style-type: none"><li>• Fatigue</li><li>• Thirst</li><li>• Excess urination</li><li>• Blurry vision</li></ul></li><li>2. Perform a control solution test using a new test strip. If the test results are within the range printed on the test strip vial, retest using blood and a new test strip.</li><li>3. If the control solution result is out of range or the error reappears, contact Customer Care.</li></ol> |

| Message | What It Means                                                                                                                                                                                                                                                                                                                                                                                                                                  | What To Do                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|---------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Er3     | <ul style="list-style-type: none"><li>• Incorrect test procedure. For example, putting blood on the test strip before inserting it into the meter or applying blood before  and  appear</li><li>• There may be a problem with the test strip</li><li>• Meter error</li></ul> | <ol style="list-style-type: none"><li>1. Wait for  and  to appear before applying blood or control solution.</li><li>2. Perform a control solution test using a new test strip. If the test results are within the range printed on the test strip vial, retest using blood and a new test strip.</li><li>3. If the control solution result is out of range or the error reappears, contact Customer Care.</li></ol> |
| Er4     | <ul style="list-style-type: none"><li>• There may be a problem with the test strip</li><li>• Meter error</li></ul>                                                                                                                                                                                                                                                                                                                             | <ol style="list-style-type: none"><li>1. Perform a control solution test using a new test strip. If the test results are within the range printed on the test strip vial, retest using blood and a new test strip.</li><li>2. If the control solution result is out of range or the error reappears, contact Customer Care.</li></ol>                                                                                                                                                                                                                                                  |

# Troubleshooting

English

|                                                                            | What It Means                                                                                                                                    | What To Do                                                                                                                          |
|----------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>1. The meter does not enter test mode after inserting a test strip.</b> | <ul style="list-style-type: none"> <li>• Test strip is inserted upside down, or bottom end in, or not fully inserted into the meter</li> </ul>   | 1. Insert test strip print side up, top end in.                                                                                     |
|                                                                            | <ul style="list-style-type: none"> <li>• No battery is installed</li> <li>• Battery is installed incorrectly</li> </ul>                          | 1. Install battery with (+) facing up.                                                                                              |
|                                                                            | <ul style="list-style-type: none"> <li>• Dead battery</li> </ul>                                                                                 | 1. Replace battery. Reset date and time, if necessary.<br>2. If meter still does <i>not</i> enter test mode, contact Customer Care. |
|                                                                            | <ul style="list-style-type: none"> <li>• Defective test strip or meter</li> <li>• Blood or foreign objects put in the test strip port</li> </ul> | 1. Contact Customer Care.                                                                                                           |

|                                                                  | What It Means                                                                                                        | What To Do                                                                                                                                                                                                                                                                                                                                                                           |
|------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>2. The test does not start after applying the blood drop.</b> | <ul style="list-style-type: none"> <li>Blood drop is too small</li> </ul>                                            | <ol style="list-style-type: none"> <li>Add blood to the same sample area within 60 seconds, or repeat the test using a new test strip and a larger blood drop.</li> </ol>                                                                                                                                                                                                            |
|                                                                  | <ul style="list-style-type: none"> <li>Defective test strip</li> <li>Sample applied after meter turns off</li> </ul> | <ol style="list-style-type: none"> <li>Repeat the test using a new test strip. Wait for  and  to appear before applying blood or control solution.</li> <li>If problem persists, contact Customer Care.</li> </ol> |
|                                                                  | <ul style="list-style-type: none"> <li>Defective meter</li> </ul>                                                    | <ol style="list-style-type: none"> <li>Contact Customer Care.</li> </ol>                                                                                                                                                                                                                                                                                                             |

## FreeStyle® Control Solution

### Purpose

FreeStyle control solution is a red liquid that contains a fixed amount of glucose.

The purpose of doing a control solution test is:

- To practice testing *without* using your own blood
- To confirm that the meter and test strips work together properly
- When you are unsure of your blood glucose test results

Contact Customer Care for information on how to obtain control solution.

### Control Solution Information

- Use only FreeStyle control solution with the meter.
- Replace the cap securely on the bottle immediately after use.
- Control solution results are accurate only between 59° – 104° F (15° – 40° C).
- Do *not* add water or other liquid to control solution.
- There are three levels of control solution available (Low, Normal, High).

For information on how to obtain control solution, contact Customer Care.

**IMPORTANT:** Control solution results should fall within the control solution range printed on the test strip vial. Do *not* use this range when testing blood glucose.

## Out of Range Control Solution Results

- Repeat the test if control solution results are outside the range printed on the test strip vial.
- Stop using the meter if control solution results are consistently outside the range printed on the test strip vial. Contact Customer Care.

Causes of out-of-range results include:

- Expired or bad control solution
- Expired or bad test strip
- Testing error
- Watered-down control solution
- Meter malfunction
- Test strip deterioration
- Control solution test done outside 59° – 104° F (15° – 40° C)

## Performing a Control Solution Test

1. Check control solution expiration date.


**IMPORTANT:** Do *not* use control solution past expiration date. Discard control solution 3 months after opening or on the expiration date printed on the bottle whichever comes first. (Example: open April 15, discard July 15; write the discard date on the side of the bottle.)

2. Remove test strip from vial.
3. Insert test strip to turn meter on.

**Note:** The meter turns off after 2 minutes of inactivity. Remove and reinsert the unused test strip to restart the meter.

If the meter does not turn on go to Troubleshooting section.

#### 4. Confirm System Check Screen.

- This display appears when you turn on the meter. **Do not use the meter if the display check screen does not exactly match the example. Contact Customer Care.**
- See The FreeStyle Freedom Lite® Meter and How It Works section for more information.
- If the system check screen does not appear, go to Troubleshooting section.

#### 5. Apply and **hold** control solution to 1 sample area of the test strip until:

- You see short lines moving clockwise on the display or
- You hear a beep

This indicates the test strip obtained enough control solution and the meter is processing. If the test does not start after applying the control solution, go to Troubleshooting section.

**IMPORTANT:** If a short line does *not* appear after 5 seconds, the sample may be too small. Do *not* add control solution to both sides of the strip. The test strips are for single use only. Do *not* reuse test strips. Discard used test strips.

(Example only)

## 6. View result.

The result appears on the display when the test is complete.

## 7. Compare the control solution result to the range printed on the test strip vial. The result should fall within the range.

**Note:** If you see an “Er1, Er2, Er3, or Er4” on the display, it is an error code. Go to Error Codes section.

(Example only)

## 8. Press and **hold** for 2 seconds until appears on the display to mark the result as a control solution result. To unmark a control solution test result, press and **hold** for 2 seconds. A value with the indicates you have marked this reading as a control solution test.

## 9. Remove the test strip to turn meter off.

# Using Your Meter Memory

## Purpose

The purpose of this section is to review your stored test results and 7-, 14-, and 30-day averages.

## Viewing Blood Glucose Averages in Memory

Start with the meter off (no test strip inserted). Press **(m)**. Your 7-day average appears (control solution, LO, and HI readings *not* included).

**Note:** To exit meter memory at any time, press and hold **(m)**.

### Example:

Memory Display

your average

n= number of results included in your average

- Press **(C)** to scroll through 14-, and 30-day averages and all stored test results. The most recent result appears first.
- Press **(m)** to scroll backward.

## Viewing Test Results in Memory

Meter memory stores the last 400 test results with time and date. Press  to scroll through the results. The most recent result appears first.

|                                                                                              |                                                                                                                                                                                                  |
|----------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>mem</b>                                                                                   | on the display indicates that you are viewing stored test results.                                                                                                                               |
| <b>davg</b>                                                                                  | shows the average value of glucose readings from the last 7-, 14-, or 30-days.                                                                                                                   |
|             | on the display indicates that the test result occurred when it was too hot or too cold to use the meter and may <i>not</i> be accurate. This result is included in 7-, 14-, and 30-day averages. |
| <b>LO</b>                                                                                    | a test result of under 20 mg/dL is included in memory as LO.                                                                                                                                     |
| <b>HI</b>                                                                                    | a test result of over 500 mg/dL is included in memory as HI.                                                                                                                                     |
| <b>mem</b>  | indicates a control solution test result.                                                                                                                                                        |

## Setting Reminder Alarms (optional)

### Purpose

The purpose of this section is to set up to four reminder alarms to prompt you to test throughout the day.

Reminder Alarm

1. With the meter off, press and **hold** **(m)** for **3** seconds until **SET** appears and **0** is flashing.
2. Press **(C)**. Reminder alarm **1** flashes.

Hour

**Note:** You are setting time for Reminder Alarm not clock time.

3. Set hour
  - Press **(m)** to set hour
  - Press **(C)** to display correct hour
  - Press **(m)** to save

Minutes

4. Set minutes
  - Press **(C)** to display correct minutes
  - Press **(m)** to save

## 5. Set alarm sound volume

| If you prefer: | Press <b>(C)</b> until you see:                                                   | Then                       |
|----------------|-----------------------------------------------------------------------------------|----------------------------|
| No sound       |  | • Press <b>(m)</b> to save |
| A low beep     |  |                            |
| A louder beep  |  |                            |

6. Press **(C)** and repeat 2, 3, 4 and 5 to set remaining reminder alarms or

7. Press and **hold** **(m)** for 3 seconds until meter turns off.

## Transferring Test Results To a Computer

PC

appears on the display when you insert a data cable in the data port. The computer must have a compatible data management system to complete the transfer. Remove this cable before performing a test.

Contact Customer Care to learn more about the data management system.


## Maintaining the Meter

### Replacing the Battery

The meter comes with a #2032, 3-volt, lithium battery installed. It provides power for about 1000 tests.

- When  appears on the display, the battery is low.
- When  appears, the meter may fail to turn on.
- When you are ready to install a new battery, remove the old battery and install the new battery within 60 seconds to avoid losing time and date settings.
- If time and date settings are lost, both time 12:00 and date 00-00 will flash until reset. Reset time and date to get the correct 7-, 14-, and 30-day averages. If the time and date settings are lost, you may still perform an accurate blood glucose test.
- Memory log and previously stored user settings, such as date and time format, are not affected by removing the battery.

**IMPORTANT:** When  appears, replace the battery immediately.

1. Slide the battery door to the right to open.
2. To remove the battery, pull the plastic tab up to release the battery.

3. Insert a new battery: (+) facing up.

4. To close, slide the battery door in until it snaps into place.

**Note:** The battery in this product should be removed and disposed in accordance with local regulations for separate collection of spent batteries.

## Cleaning and Disinfecting the Meter and Lancing Device

This section describes the cleaning and disinfection tasks for the FreeStyle Freedom Lite Meter and the FreeStyle Lancing Device II.

Cleaning and disinfecting your meter and lancing device is important to prevent the spread of infectious diseases. The meter has a use life of 5 years, which is 522 cleaning and 522 disinfection cycles (2 cycles per week for 5 years). The lancing device has a use life of 2 years, which is 210 cleaning and 210 disinfection cycles (2 cycles per week for 2 years).

Turn off the meter prior to cleaning and disinfection. Users should wash hands thoroughly with soap and water after handling the meter, lancing device, or test strips.

If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be cleaned and disinfected prior to use by the second person.

**Caution:** Do *not* place the meter in water or other liquids. Avoid getting dust, dirt, blood, control solution, water, bleach, or any other substance in the meter's test strip port, data port and battery compartment.

### Assemble cleaning and disinfection supplies.

Clorox Healthcare Bleach Germicidal Wipes, EPA Reg. #67619-12. These disinfectant wipes contain a 0.55% Sodium Hypochlorite (NaOCl) solution and have been shown to be safe for use with the meter and lancing device. They may be purchased at major online retailers, such as Walmart.com, Amazon.com, and OfficeDepot.com.

- Cleaning and disinfection should be performed using Clorox Healthcare Bleach Germicidal Wipes at least once per week.
- The meter has a use life of 5 years and has been validated for 522 cleaning and 522 disinfection cycles (the equivalent of 2 cycles per week for 5 years).
- The lancing device has a use life of 2 years and has been validated for 210 cleaning and disinfection cycles (the equivalent of 2 cycles per week for 2 years).
- Cleaning and disinfecting the meter is important to prevent the spread of infectious diseases.

### **Cleaning and Disinfection**

Cleaning is the physical removal of organic soil from the meter and lancing device surfaces. Keeping the meter and lancing device clean helps ensure that they are working properly and that no dirt gets in the devices. Cleaning allows for successful, subsequent disinfection.

Disinfection is a process that destroys pathogens, such as viruses and other microorganisms, on the meter and lancing device surfaces. Disinfecting the meter and lancing device helps ensure that no infection is passed on when you or others come in contact with the meter and lancing device. The meter and lancing device should be cleaned with a separate bleach wipe prior to each disinfection.

Clean and disinfect your meter at least once a week. Clean and disinfect your lancing device when visibly dirty (at least once a week).

**Note:** Squeeze bleach wipe to remove excess solution before wiping the meter or lancing device.

**IMPORTANT:** If you require assistance or if you notice any signs of deterioration on the meter or lancing device (such as clouding or crazing on the display of the meter, corroding or eroding of the plastic housing, or cracking of plastic housing or display) or if the meter does not turn on or if the lancing device does not work after cleaning and disinfection, discontinue use of the meter or lancing device and contact Customer Care (1-888-522-5226).

## Cleaning

- |               |                                                                                                                                                                                                                                                                                                                                                                                            |
|---------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Step 1</b> | For cleaning, wipe all outside surfaces of the meter and lancing device thoroughly with the bleach wipes until they are visibly clean. Clean the strip port, data port and battery compartment with caution and do not force liquids into sensitive meter electronic openings. Clean your meter at least once a week. Clean your lancing device when visibly dirty (at least once a week). |
| <b>Step 2</b> | Allow the meter and lancing device to air dry completely before performing a blood glucose test.                                                                                                                                                                                                                                                                                           |
| <b>Step 3</b> | When finished, thoroughly wash your hands with soap and water.                                                                                                                                                                                                                                                                                                                             |

**Disinfection**

|               |                                                                                                                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Step 1</b> | For disinfection, wipe all outside surfaces of the meter and lancing device thoroughly with the bleach wipe until they are wet in order to disinfect the meter and lancing device. Clean the strip port, data port and battery compartment with caution and do not force liquids into sensitive meter electronic openings. |
| <b>Step 2</b> | Allow the meter and lancing device surfaces to remain wet for 60 seconds (1 full minute). After 1 full minute, immediately dry with drying towel to remove any residual moisture. Disinfect your meter at least once a week. Disinfect your lancing device when visibly dirty (at least once a week).                      |
| <b>Step 3</b> | Allow the meter and lancing device to air dry completely before performing a blood glucose test.                                                                                                                                                                                                                           |
| <b>Step 4</b> | When finished, thoroughly wash your hands with soap and water.                                                                                                                                                                                                                                                             |

Additional information about the risks for transmitting bloodborne pathogens to persons undergoing fingerstick procedures for blood sampling can be found in the following references:

1. "FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010) [http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025 .htm](http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm)
2. "CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010) <http://www .cdc .gov/injectionsafety/Fingerstick-DevicesBGM.html>

## Specifications

|                                    |                                                                 |
|------------------------------------|-----------------------------------------------------------------|
| <b>Assay method</b>                | Coulometric electrochemical sensor                              |
| <b>Automatic shutoff</b>           | 2 minutes of inactivity                                         |
| <b>Battery life</b>                | 1000 tests                                                      |
| <b>Calibration</b>                 | Plasma equivalent                                               |
| <b>Hematocrit</b>                  | 15% to 65%                                                      |
| <b>Measurement units</b>           | mg/dL                                                           |
| <b>Meter storage temperature</b>   | – 4° to 140° F (–20° to 60° C)                                  |
| <b>Memory</b>                      | 400 blood glucose and control solution tests with date and time |
| <b>Operating relative humidity</b> | 5% to 90% (non-condensing)                                      |

|                                                 |                                                                                                                                  |
|-------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| <b>Operating temperature</b>                    | 40° to 104° F (4° to 40° C)                                                                                                      |
| <b>Power source</b>                             | One CR 2032, 3V lithium battery, replaceable                                                                                     |
| <b>Result range</b>                             | 20 to 500 mg/dL                                                                                                                  |
| <b>Sample</b>                                   | Fresh capillary and venous whole blood samples                                                                                   |
| <b>Sample size</b>                              | 0.3 microliter (300 nanoliters)                                                                                                  |
| <b>Size</b>                                     | 2.0 in (w) x 3.3 in (l) x 0.63 in (d) 5.1 cm (w) x 8.4 cm (l) x 1.6 cm (d)                                                       |
| <b>Weight</b>                                   | 1.4 to 1.6 oz. (39.7 to 45.4 g) including battery                                                                                |
| <b>Meter cleaning and disinfection</b>          | The meter has a use life of 5 years, which is 522 cleaning and 522 disinfection cycles (2 cycles per week for 5 years).          |
| <b>Lancing device cleaning and disinfection</b> | The lancing device has a use life of 2 years, which is 210 cleaning and 210 disinfection cycles (2 cycles per week for 2 years). |

## Definition of Symbols

English

|                                                                                   |                                                                                                                                             |                                                                                   |                                               |
|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------|
|  | Consult instructions for use                                                                                                                |  | For <i>in vitro</i> diagnostic medical device |
|  | Temperature limitation                                                                                                                      |  | Catalog number                                |
|  | Manufacturer                                                                                                                                |  | Recycle                                       |
|  | Batch code                                                                                                                                  |  | Serial number                                 |
|  | Use by                                                                                                                                      |  | Do <i>not</i> reuse                           |
|  | The battery in this product should be removed and disposed in accordance with local regulations for separate collection of spent batteries. |                                                                                   |                                               |

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FreeStyle® and related brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions.
Abbott Diabetes Care Inc.