



FreeStyle
Libre 3

Mobile Device & OS Compatibility

[Click here for English](#)

[Haga clic aquí para español](#)

The circular shape of the sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott. Other trademarks are the property of their respective owners.

©2021 Abbott ART44628-002 Rev. A 09/21



MOBILE DEVICE & OS COMPATIBILITY

Popular mobile devices and operating systems (OS) are regularly tested to evaluate NFC scan performance, Bluetooth connectivity, and app compatibility with Sensors. We recommend checking this guide before installing a new OS version on your phone or before using the App with a new phone.

Recommended app, device, and operating systems

APP / VERSION	DEVICES	OS
FreeStyle Libre 3 (version 3.2)	iPhone X iPhone 12 Pro Max	14.0.1, 14.6

This table will be updated as other devices and operating systems are evaluated. Contact Customer Service for more information on devices and operating systems other than those listed. Beta OS versions are not evaluated or supported. The App cannot be used with jailbroken and rooted devices.



DISPOSITIVO MÓVIL Y COMPATIBILIDAD CON EL SISTEMA OPERATIVO

Los dispositivos móviles y sistemas operativos más populares se prueban regularmente para evaluar el rendimiento del escaneo NFC, la conectividad Bluetooth y la compatibilidad de las aplicaciones con los sensores. Le recomendamos que consulte esta guía antes de instalar una nueva versión del sistema operativo en su teléfono o antes del uso de la aplicación con un nuevo teléfono.

Aplicación, dispositivo y sistemas operativos recomendados

APLICACIÓN/VERSIÓN	DISPOSITIVOS	SISTEMA OPERATIVO
FreeStyle Libre 3 (versión 3.2)	iPhone X iPhone 12 Pro Max	14.0.1, 14.6

Esta tabla se actualizará a medida que se evalúen otros dispositivos y sistemas operativos. Comuníquese con Atención al cliente para obtener más información sobre dispositivos y sistemas operativos distintos de los enumerados. No se evalúan ni se brinda soporte a las versiones beta de los sistemas operativos. La aplicación no se puede usar con dispositivos que tienen hecho el jailbreak o “ruteados”.

