



FreeStyle Libre 2

Compatible Devices, Apps, Software, and Systems

Below is a list of systems, devices, apps, and software that can be used with the FreeStyle Libre 2 Sensor.

Devices:

FreeStyle Libre 2 Reader

Apps:

FreeStyle Libre 2 app

Bigfoot Unity™ App

Software that you can use to upload and store Sensor glucose data from your Reader or App:

FreeStyle Libre Software (Reader only)

LibreView Software (Reader or App)

Systems:

Bigfoot Unity™ Diabetes Management System

This list will be updated as other devices, apps, software, and systems are evaluated by Abbott Diabetes Care. Use of your Sensor with devices, apps, software, and systems that are not listed may cause inaccurate glucose readings.

