

# Mobile Device & OS Compatibility

[Click here for English](#)

[Haga clic aquí para ver la versión en español](#)

The circular shape of the sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott. Other trademarks are the property of their respective owners.

©2021 Abbott ART41556-102 Rev. D 11/21



# MOBILE DEVICE & OS COMPATIBILITY

Popular mobile devices and operating systems (OS) are regularly tested to evaluate NFC scan performance, Bluetooth connectivity, and app compatibility with Sensors. We recommend checking this guide before installing a new OS version on your phone or before using the app with a new phone.

## Recommended app, device, and operating systems

APP / VERSION	DEVICES	OS
FreeStyle Libre 2 (version 2.7.1)	iPhone 7, 7 Plus, 8, 8 Plus, X, XR, XS, XS Max, 11, 11 Pro, 11 Pro Max, 12, 12 Pro, 12 Pro Max, 12 Mini, SE (2020), 13, 13 Pro	iOS: 14.0.1, 14.1, 14.2.1, 14.3, 14.4.1, 14.4.2, 14.5, 14.5.1, 14.6, 14.7, 14.7.1, 14.8, 15.0, 15.1
FreeStyle LibreLink (version 2.5.3)	iPhone 7, 7 Plus, 8, 8 Plus, X, XR, XS, XS Max, 11, 11 Pro, 11 Pro Max, 12, 12 Pro, 12 Pro Max, 12 Mini, SE (2020), 13	iOS: 13.2, 13.3, 13.4, 13.5, 13.6, 13.7, 14.0, 14.1, 14.2, 14.5.1, 14.6, 14.7.1, 14.8, 15.0, 15.1
	Samsung Galaxy Note 8, S7 Edge, S8, S8+, S9, S10, A32 5G, A52 5G	Android: 7.1, 8, 9, 10, 11, 12
	Google Pixel, Pixel 2, Pixel 2 XL, Pixel 3, Pixel 3 XL, Pixel 4, Pixel 4 XL	
	Sony Xperia 1, Xperia 5, Xperia 8, Xperia 10 II	
	LG Nexus 5X	

This table will be updated as other devices and operating systems are evaluated. The contents of this guide are the most up to date and may be different from what is displayed within the FreeStyle Libre 2 app. Contact Customer Service for more information on devices and operating systems other than those listed. Jailbroken devices and beta OS versions are not evaluated or supported.

# COMPATIBILIDAD CON DISPOSITIVOS MÓVILES Y SISTEMAS OPERATIVOS

Los dispositivos móviles y sistemas operativos más populares se prueban regularmente para evaluar el rendimiento del escaneo NFC, la conectividad Bluetooth y la compatibilidad de las aplicaciones con los sensores. Le recomendamos que consulte esta guía antes de instalar una nueva versión del sistema operativo en su teléfono o antes de utilizar la aplicación con un nuevo teléfono.

## Aplicación, dispositivo y sistemas operativos recomendados

APLICACIÓN / VERSIÓN	DISPOSITIVOS	SISTEMA OPERATIVO
FreeStyle Libre 2 (versión 2.7.1)	iPhone 7, 7 Plus, 8, 8 Plus, X, XR, XS, XS Max, 11, 11 Pro, 11 Pro Max, 12, 12 Pro, 12 Pro Max, 12 Mini, SE (2020), 13, 13 Pro	iOS: 14.0.1, 14.1, 14.2.1, 14.3, 14.4.1, 14.4.2, 14.5, 14.5.1, 14.6, 14.7, 14.7.1, 14.8, 15.0, 15.1
FreeStyle LibreLink (versión 2.5.3)	iPhone 7, 7 Plus, 8, 8 Plus, X, XR, XS, XS Max, 11, 11 Pro, 11 Pro Max, 12, 12 Pro, 12 Pro Max, 12 Mini, SE (2020), 13	iOS: 13.2, 13.3, 13.4, 13.5, 13.6, 13.7, 14.0, 14.1, 14.2, 14.5.1, 14.6, 14.7.1, 14.8, 15.0, 15.1
	Samsung Galaxy Note 8, S7 Edge, S8, S8+, S9, S10, A32 5G, A52 5G	Android: 7.1, 8, 9, 10, 11, 12
	Google Pixel, Pixel 2, Pixel 2 XL, Pixel 3, Pixel 3 XL, Pixel 4, Pixel 4 XL	
	Sony Xperia 1, Xperia 5, Xperia 8, Xperia 10 II	
	LG Nexus 5X	

Esta tabla se actualizará a medida que se evalúen otros dispositivos y sistemas operativos. El contenido de esta guía es el más actualizado y puede diferir de lo que se muestra en la aplicación FreeStyle Libre 2. Comuníquese con el Servicio de atención al cliente para obtener más información sobre dispositivos y sistemas operativos diferentes a los enumerados. No se evalúan ni se brinda soporte a los dispositivos que tienen hecho el jailbreak ni a las versiones beta de los sistemas operativos.