

# Sensor Compatibility Guide

[Click here for English](#)

[Haga clic aquí para español](#)

The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott.



Make sure your app and Sensor are compatible.

## Compatible Sensors and Apps

<b>SENSOR</b>	<b>APP</b>
FreeStyle Libre (10 day)	FreeStyle LibreLink
FreeStyle Libre 14 day	FreeStyle LibreLink
FreeStyle Libre 2	FreeStyle Libre 2
FreeStyle Libre 2 Plus	FreeStyle Libre 2
FreeStyle Libre 3	FreeStyle Libre 3
FreeStyle Libre Pro	App not available

Note: These apps may not be compatible if your Sensor is from a different region.

Contact Support [1-855-632-8658](tel:1-855-632-8658)

[www.FreeStyleLibre.com](http://www.FreeStyleLibre.com)

©2023 Abbott ART44626-002 Rev. A 06/23



Asegúrese de que su aplicación y el sensor sean compatibles.

## Sensores y aplicaciones compatibles

<b>SENSOR</b>	<b>APLICACIÓN</b>
FreeStyle Libre (10 días)	FreeStyle LibreLink
FreeStyle Libre 14 días	FreeStyle LibreLink
FreeStyle Libre 2	FreeStyle Libre 2
FreeStyle Libre 2 Plus	FreeStyle Libre 2
FreeStyle Libre 3	FreeStyle Libre 3
FreeStyle Libre Pro	Aplicación no disponible

Nota: es posible que estas aplicaciones no sean compatibles si su sensor es de una región diferente.

Póngase en contacto con la asistencia técnica [1-855-632-8658](tel:1-855-632-8658)

[www.FreeStyleLibre.com](http://www.FreeStyleLibre.com)

©2023 Abbott ART44626-002 Rev. A 06/23

