



## Phone & OS Compatibility

**[Click here](#)**

The sensor housing, Libre, and related brand marks are marks of Abbott. Other trademarks are the property of their respective owners.

©2025 Abbott ART48719-001 Rev. D 03/25




# PHONE & OS COMPATIBILITY

We regularly test popular phones and operating systems (OS) to evaluate App compatibility, NFC scan performance, and Bluetooth connectivity.

Please check this guide before installing a new OS version on your phone or before using the App with a new phone.

**Note:** Refer to the diagram for the NFC scan location on your phone.



## Recommended phones and OS

### Libre app (version 1.0)

OS	PHONES	NFC LOCATION
iOS: 16, 17, 17.0.1, 17.0.2, 17.0.3, 17.1, 17.1.1, 17.1.2, 17.2, 17.3, 17.3.1, 17.4, 17.4.1, 17.5, 17.5.1, 18, 18.0.1	iPhone 8, 8 Plus, X, XR, XS, XS Max, 11, 11 Pro, 11 Pro Max, 12, 12 Pro, 12 Pro Max, 12 Mini, SE (2020, 2nd Generation), SE (2022, 3rd Generation), 13, 13 Pro, 13 Pro Max, 13 Mini, 14, 14 Plus, 14 Pro, 14 Pro Max, 15, 15 Plus, 15 Pro, 15 Pro Max, 16, 16 Plus, 16 Pro, 16 Pro Max	1
Android: 11, 12, 12.1, 13, 14, 15* (*Pixel phones only)	Samsung Galaxy A12, A15 5G, A22 5G, A71 5G	1
	Google Pixel 2, Pixel 2 XL, Pixel 3	
	Dialn G65	
	Samsung Galaxy A20, A42 5G, A52 5G, Note10, Note10+, Note20 5G, Note20 Ultra 5G, S10, S10+, S10e, S20 5G, S20 Ultra 5G, S20+ 5G, S21 5G, S21+ 5G, S22, S22 Ultra, S22+, S23, S23 Ultra, S24, S24 Ultra, S24+, Z Fold2 5G, Z Fold3 5G, ZFold4	2
	Google Pixel 3 XL, Pixel 3a, Pixel 4, Pixel 4a, Pixel 4a 5G, Pixel 4 XL, Pixel 5, Pixel 6, Pixel 6 Pro, Pixel 6a, Pixel 7, Pixel 7 Pro, Pixel 8 Pro	
	Motorola one 5G ace	
	Samsung Galaxy Note10+ 5G, S21 Ultra 5G, Z Flip3 5G, Z Flip4, Z Flip5, Z Fold5	3

We will update this guide as we test other phones and OS. This guide is your most up to date information and may be different from what is shown within the App. Contact Customer Service for more information on phones and OS other than those listed. Beta OS versions are not tested or supported. You cannot use the App with jailbroken and rooted devices.

Using an unsupported phone or OS may affect App security and functionality. You are responsible for the risks of using the App on an unsupported phone or OS.