



## Phone & OS Compatibility

**Click here**

**Haga clic aquí**

The sensor housing, Libre, and related brand marks are marks of Abbott. Other trademarks are the property of their respective owners.


©2025 Abbott ART48719-201 Rev. A 10/25



We regularly test devices and operating systems (OS) to evaluate app compatibility, NFC scan performance, and Bluetooth connectivity.

Please check this guide before installing a new OS version on your device or before using the app with a new device.

**Note:** Refer to the diagram for the NFC scan location on your phone.



Libre app (version 1.2.0)

Recommended phones and OS

| OS                                                                                                                                              | PHONES                                                                                                                                                                                                                                                                                                        | NFC LOCATION |
|-------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| iOS: 16, 17, 17.0.1, 17.0.2, 17.0.3, 17.1, 17.1.1, 17.1.2, 17.2, 17.3, 17.3.1, 17.4, 17.4.1, 17.5, 17.5.1, 18, 18.0.1, 18.4, 18.4.1, 18.5, 18.6 | iPhone 8, 8 Plus, X, XR, XS, XS Max, 11, 11 Pro, 11 Pro Max, 12, 12 Pro, 12 Pro Max, 12 Mini, SE (2020, 2nd generation), SE (2022, 3rd generation), 13, 13 Pro, 13 Pro Max, 13 Mini, 14, 14 Plus, 14 Pro, 14 Pro Max, 15, 15 Plus, 15 Pro, 15 Pro Max, 16, 16 Plus, 16 Pro, 16 Pro Max, 16e                   | 1            |
| Android: 11, 12, 12.1, 13, 14, 15, 16*                                                                                                          | Samsung Galaxy A12, A14 5G, A15 5G, A22 5G, A71 5G                                                                                                                                                                                                                                                            | 1            |
|                                                                                                                                                 | Google Pixel 2, Pixel 2 XL, Pixel 3                                                                                                                                                                                                                                                                           |              |
|                                                                                                                                                 | Dialn G65                                                                                                                                                                                                                                                                                                     |              |
|                                                                                                                                                 | Motorola moto g 5G 2024, moto g stylus 5G 2024                                                                                                                                                                                                                                                                |              |
|                                                                                                                                                 | Samsung Galaxy A20, A42 5G, A52 5G, Note10, Note10+, Note20 5G, Note20 Ultra 5G, S10, S10+, S10e, S20 5G, S20 Ultra 5G, S20+ 5G, S21 5G, S21+ 5G, S21 FE 5G, S22, S22 Ultra, S22+, S23, S23 Ultra, S23+, S23 FE, S24, S24 Ultra, S24+, S24 FE, S25, S25 Ultra, S25+, Z Fold2 5G, Z Fold3 5G, Z Fold4, Z Fold6 | 2            |
|                                                                                                                                                 | Google Pixel 3 XL, Pixel 3a, Pixel 4, Pixel 4a, Pixel 4a 5G, Pixel 4 XL, Pixel 5, Pixel 6, Pixel 6 Pro, Pixel 6a, Pixel 7, Pixel 7 Pro, Pixel 8 Pro, Pixel 8a, Pixel 9 Pro, Pixel 9 Pro XL, Pixel 9 Pro Fold                                                                                                  |              |
|                                                                                                                                                 | Motorola one 5G ace                                                                                                                                                                                                                                                                                           |              |
|                                                                                                                                                 | Samsung Galaxy Note10+ 5G, S21 Ultra 5G, Z Flip3 5G, Z Flip4, Z Flip5, Z Fold5, Z Flip6                                                                                                                                                                                                                       | 3            |

Recommended watches and OS

| OS                                    | WATCHES                                                                                                                                  |
|---------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| watchOS: 10.6.1, 11.0.1, 11.5, 11.6.1 | Apple Watch: Series 4, Series 5, SE 1st generation, Series 6, Series 7, SE 2nd generation, Series 8, Ultra, Series 9, Ultra 2, Series 10 |

\*Pixel phones only.


We will update this guide as we test other devices and OS. This guide is your most up to date information and may be different from what is shown within the app. Contact Customer Service for more information on devices and OS other than those listed. Beta OS versions are not tested or supported. You cannot use the app with jailbroken and rooted devices.

Using an unsupported device or OS may affect app security and functionality. You are responsible for the risks of using the app on an unsupported device or OS.

Probamos periódicamente dispositivos y sistemas operativos para evaluar la compatibilidad de las aplicaciones, el rendimiento del escaneo NFC y la conectividad Bluetooth.

Consulte esta guía antes de instalar una nueva versión del sistema operativo en su dispositivo o antes del uso de la aplicación con un nuevo dispositivo.

**Nota:** Consulte el diagrama para conocer la ubicación del escáner NFC en su teléfono.



Aplicación Libre (versión 1.2.0)

Teléfonos y sistema operativo recomendados

| SISTEMA OPERATIVO                                                                                                                               | TELÉFONOS                                                                                                                                                                                                                                                                                                     | UBICACIÓN NFC |
|-------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|
| iOS: 16, 17, 17.0.1, 17.0.2, 17.0.3, 17.1, 17.1.1, 17.1.2, 17.2, 17.3, 17.3.1, 17.4, 17.4.1, 17.5, 17.5.1, 18, 18.0.1, 18.4, 18.4.1, 18.5, 18.6 | iPhone 8, 8 Plus, X, XR, XS, XS Max, 11, 11 Pro, 11 Pro Max, 12, 12 Pro, 12 Pro Max, 12 Mini, SE (2020, 2nd generation), SE (2022, 3rd generation), 13, 13 Pro, 13 Pro Max, 13 Mini, 14, 14 Plus, 14 Pro, 14 Pro Max, 15, 15 Plus, 15 Pro, 15 Pro Max, 16, 16 Plus, 16 Pro, 16 Pro Max, 16e                   | 1             |
| Android: 11, 12, 12.1, 13, 14, 15, 16*                                                                                                          | Samsung Galaxy A12, A14 5G, A15 5G, A22 5G, A71 5G                                                                                                                                                                                                                                                            | 1             |
|                                                                                                                                                 | Google Pixel 2, Pixel 2 XL, Pixel 3                                                                                                                                                                                                                                                                           |               |
|                                                                                                                                                 | Dialn G65                                                                                                                                                                                                                                                                                                     |               |
|                                                                                                                                                 | Motorola moto g 5G 2024, moto g stylus 5G 2024                                                                                                                                                                                                                                                                |               |
|                                                                                                                                                 | Samsung Galaxy A20, A42 5G, A52 5G, Note10, Note10+, Note20 5G, Note20 Ultra 5G, S10, S10+, S10e, S20 5G, S20 Ultra 5G, S20+ 5G, S21 5G, S21+ 5G, S21 FE 5G, S22, S22 Ultra, S22+, S23, S23 Ultra, S23+, S23 FE, S24, S24 Ultra, S24+, S24 FE, S25, S25 Ultra, S25+, Z Fold2 5G, Z Fold3 5G, Z Fold4, Z Fold6 | 2             |
|                                                                                                                                                 | Google Pixel 3 XL, Pixel 3a, Pixel 4, Pixel 4a, Pixel 4a 5G, Pixel 4 XL, Pixel 5, Pixel 6, Pixel 6 Pro, Pixel 6a, Pixel 7, Pixel 7 Pro, Pixel 8 Pro, Pixel 8a, Pixel 9 Pro, Pixel 9 Pro XL, Pixel 9 Pro Fold                                                                                                  |               |
|                                                                                                                                                 | Motorola one 5G ace                                                                                                                                                                                                                                                                                           |               |
|                                                                                                                                                 | Samsung Galaxy Note10+ 5G, S21 Ultra 5G, Z Flip3 5G, Z Flip4, Z Flip5, Z Fold5, Z Flip6                                                                                                                                                                                                                       | 3             |

Relojes y sistema operativo recomendados

| SISTEMA OPERATIVO                     | RELOJES                                                                                                                                  |
|---------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| watchOS: 10.6.1, 11.0.1, 11.5, 11.6.1 | Apple Watch: Series 4, Series 5, SE 1st generation, Series 6, Series 7, SE 2nd generation, Series 8, Ultra, Series 9, Ultra 2, Series 10 |

\*Solo teléfonos Pixel.

Actualizaremos esta guía a medida que probemos otros dispositivos y sistemas operativos. Esta guía es la información más actualizada y puede diferir de lo que se muestra en la aplicación. Póngase en contacto con Atención al cliente para obtener más información sobre dispositivos y sistemas operativos distintos a los de la lista. No se hacen pruebas ni se brinda soporte a las versiones beta de los sistemas operativos. No puede usar la aplicación con dispositivos que tienen hecho el jailbreak o “ruteados”.

Usar un dispositivo o un sistema operativo no compatibles puede afectar la seguridad y la funcionalidad de la aplicación. Usted es responsable de los riesgos de utilizar la aplicación en un dispositivo o sistema operativo no compatible.