

Sensor Compatibility Guide

[Click here for English](#)

[Haga clic aquí para español](#)

The circular shape of the sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott.



Make sure your app and Sensor are compatible.

Compatible Sensors and Apps

SENSOR	APP
FreeStyle Libre (10 day)	FreeStyle LibreLink
FreeStyle Libre 14 day	FreeStyle LibreLink
FreeStyle Libre 2	FreeStyle Libre 2
FreeStyle Libre 3	FreeStyle Libre 3
FreeStyle Libre Pro	App not available

Note: These apps may not be compatible if your Sensor is from a different region.

Contact Support [1-855-632-8658](tel:1-855-632-8658)

www.FreeStyleLibre.com

©2021 Abbott ART44626-001 Rev. A 09/21



Asegúrese de que su aplicación y el sensor sean compatibles.

Sensores y aplicaciones compatibles

SENSOR	APLICACIÓN
FreeStyle Libre (10 días)	FreeStyle LibreLink
FreeStyle Libre 14 días	FreeStyle LibreLink
FreeStyle Libre 2	FreeStyle Libre 2
FreeStyle Libre 3	FreeStyle Libre 3
FreeStyle Libre Pro	Aplicación no disponible

Nota: es posible que estas aplicaciones no sean compatibles si su sensor es de una región diferente.

Póngase en contacto con la asistencia técnica [1-855-632-8658](tel:1-855-632-8658)

www.FreeStyleLibre.com

©2021 Abbott ART44626-001 Rev. A 09/21

